

Item #	Seat	UOM
COBRAGT420CS	20" Captain's Seat	1/cs
COBRAGT422CS	22" Captain's Seat	1/cs



Cobra GT4 Heavy Duty Scooter

- New, innovative design for a distinctive, stylish appearance.
- One-touch, infinite tiller angle adjustment allows the user to easily position the tiller into the ideal driving position.
- Ergonomic delta control handles makes it easier for users with limited strength and dexterity to operate the scooter.
- Large, air-filled, low-profile tires and full suspension for smoothest outdoor ride quality.
- “Blade” style alloy wheel rims.
- Efficient, extra bright LED lighting.
- Choice of Full size, Sport Style Captain’s Seats in 20” or 22” widths.
- Seats are height adjustable, swivel and feature a fold-down, reclining backrest. Seats can be adjusted forward or backward on auto-style slides.
- Padded armrests feature easy turn cylinders for angle adjustment.
- Positioning belt included.
- “High/Low” switch for instant speed control.
- Anti-tip wheels for greater safety.
- Easy, freewheel mode

specifications

HCPCS	N/A	Batteries (pair)	12V x 75AH
Weight Capacity	450 lbs.	Charger	5A Offboard
Top Speed	10 mph	Brakes	Electromagnetic
Maximum Range	35 miles	Front Wheels	14.3" x 4" Pneumatic
Turning Radius	82"	Rear Wheels	14.3" x 4" Pneumatic
Climbing Angle	8°	Base Weight	191 lbs.
Ground Clearance	7"	Battery Weight (pair)	94 lbs.
Dimensions	63" (L) x 28.8" (W)	Seat Weight	55 lbs.
Floor-to-Seat Height	27" - 30"	Total Weight	340 lbs.
Seat Dimensions	20" (W) x 20" (D) 22" (W) x 20" (D)	Warranty of Frame	Lifetime
Controller	Dynamic Rhino II 160A	Warranty on Electronic Controller/Drive Train Components	14 Months
Motor	24V x 950W x 5200 rpm	Warranty on Batteries	6 months