



Item #	Seat	Battery	UOM
SPITFIRE132016FS12	16" Folding Seat	12AH	1/cs
SPITFIRE132016FS21	16" Folding Seat	21AH	1/cs



HCPCS
K0800



Spitfire EX 1320 Compact Travel Scooter

- Includes interchangeable color panels in Red, Blue and Silver.
- Quick and easy disassembly. (Figure A)
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height adjustable swivel seat with two-tone upholstery.
- Flat-free, non-marking tires.
- Large, plastic carry basket with removable liner bag.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Available with 12AH or 21AH batteries.
- 21AH batteries increase scooter's maximum range.
- Swivel seat with fold-down backrest.
- Delta tiller.
- Armrest-mounted reflectors.
- Easy freewheel operation.

specifications

HCPCS	K0800	Charger	2A Offboard
Weight Capacity	300 lbs.	Brakes	Electromagnetic
Top Speed	5 mph	Front Wheel	7.5" x 2" Flat-free
Maximum Range	9 miles (12AH) 15 miles (21AH)	Rear Wheels	8.5" x 3" Flat-free
Turning Radius	46"	Base Weight	60 lbs.
Climbing Angle	6°	Battery Weight (pair)	22 lbs. (12AH) 28 lbs. (21AH)
Ground Clearance	3"	Seat Weight	17 lbs.
Dimensions	39" (L) x 19.5" (W)	Total Weight	99 lbs. (12AH) 105 lbs. (21AH)
Floor-to-Seat Height	19" - 21"	Heaviest Piece Weight	34 lbs. (Front Section)
Seat Type	Folding	Warranty on Frame	Lifetime
Seat Dimensions	17" (W) x 16" (D)	Warranty on Electronic Controller/Drive Train Components	14 Months
Batteries (pair)	12V x 12AH 12V x 21AH	Warranty on Batteries	6 Months
Motor	24V x 270W x 4700 rpm		