

Med-Aire 5" Alternating Pressure and Low Air Loss Overlay System

HCPCS: E0372



FIGURE A



FIGURE B

- Active mattress overlay system provides both alternating pressure and low air loss to optimize pressure redistribution, shear/friction reduction, and microclimate control
- Designed to prevent, treat and heal pressure ulcers in the home or long term care setting
- Control unit offers 10 minute alternating cycles
- Adjustable patient weight settings allow for optimal immersion and patient comfort & compliance
- Pillow feature maintains the 3 air cells at the head of the bed in static mode for patient comfort
- 5" deep air cells are easily removed for convenient servicing
- Removable Nylon cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture
- CPR valve allows for rapid deflation
- Tubing quickly and easily disconnects with EZ Lock/Quick Release design

PRODUCTS

Description

14025N	Med Aire Low Air Loss, 5" 1/cs, RTL
--------	-------------------------------------

ACCESSORIES

Description

14025M	5"Mattress only for 14025,1/ea, EA
14025-CELL	Bladder Only, 14025,1/ea", EA
14025-CELL 1	Individual Bladder with Lazer Holes,1/ea, EA
14025NP	Pump for 14025N, 1/ea, EA
14025C	Cover for 14025,1/ea, EA