

# Long Term Care Empowering Better Care



# Beautiful, Functional Designs to Attract New Residents

You know the challenges. We know the solutions. Our full suite of Long Term Care (LTC) solutions was crafted specifically to meet the needs of your residents and make it easy for families to choose your facility. The functional, stylish designs exceed the rigorous challenges and expectations that many LTC facilities face while offering unique benefits to the industry:

- Clinical representatives with the research and data to find the best solutions for your facility
- · Local representatives ready to provide training, education, and support
- Unsurpassed value across our full line of products for Long Term Care facilities:
  - LTC Beds
  - · Therapeutic Support Surfaces
  - · Resident Room Furniture
  - Resident Lifts
  - Mobility

- Durable Medical Equipment
- Respiratory
- Bariatric
- Fall Prevention
- Showering and Bathing Systems · Recliners
  - Dining & Common Room Furniture

Equip your facility with products that elevate the daily experience for your residents and staff and enjoy expert support every step of the way.



#### Meet the Clinical Care Team

Our dedicated clinical care team is committed to providing exceptional support and expertise. Our team consists of highly skilled professionals who specialize in selecting, managing, and optimizing medical equipment to enhance the quality of life for your residents.

#### **Clinical Care Leadership Team**



Jason Moskowitz
VP, Sales - Clinical Care

jmoskowitz@drivemedical.com

Jason has been with Drive Medical for fifteen years. He started his career in Homecare and has spent the past four years leading the LTC strategy

for Drive Medical.



Lance Lerch
VP, Chain Accounts

llerch@drivemedical.com

Lance has been in healthcare sales for over 35 years with a focus on equipment for the LTC industry. He joined Drive Medical 12 years ago to build the LTC business, working closely with LTC chains and key distributor partners.



Karen Lerner VP, Clinical

klerner@drivemedical.com

Karen is a clinician with wide-ranging experience including DME, wound care, hospice, LTC, and rehab. She develops programs to help improve and expand aptitude in clinically-based equipment selection. She began her career as a critical care clinical specialist in major teaching hospitals and has been with Drive Medical since 2013.



**Daniel Bozeman** RVP, Corporate Accounts

dbozeman@drivemedical.com

Daniel joined Drive Medical over 8 years ago. He helps develop and implement the LTC strategy.



Wil Fraser
VP, Corporate Accounts
wfraser@drivemedical.com

Wil has 40 years of industry experience and over two decades dedicated to Drive Medical. He brings unparalleled experience and leadership to the team.

#### **Clinical Care Regional Sales Managers**



Aaron Sivley
RSM - Clinical Care
asivley@drivemedical.com

Covered Territories: AL, AR, LA, MS, OK, TX



Fernando DeCordova RSM - Clinical Care fdecordova@drivemedical.com

Covered Territories: CO, IA, ID, ND, KS, NE, SD, UT, WY



**Jerry Lyda** *RSM - Clinical Care*jlyda@drivemedical.com

Covered Territories: IL, MN, MO, WI



John Diefenbach

RSM - Clinical Care

jdiefenbach@drivemedical.com

Covered Territories: WA, OR, MT



**Justin Bono** *RSM - Clinical Care*jbono@drivemedical.com

Covered Territories: FL, GA, NJ, NY, PR, SC



Matt Herbert
RSM - Clinical Care
mherbert@drivemedical.com

Covered Territories: CT, DC, DE, MA, MD, ME, NC, NH, PA, RI, VA, VT, WV



Rachel Sorensen
RSM - Clinical Care
rsorensen@drivemedical.com

Covered Territories: IN, KY, MI, OH, TN



Sophia Mateo RSM - Clinical Care smateo@drivemedical.com

Covered Territories: AK, AZ, CA, HI, NM, NV

# PRIMECARE Model #P750

#### **FEATURES**

- · 750 lb. weight capacity
- Exceptional travel range, 11"-30"
- Integrated push-button width expansion 39"-42"-48"
- · Length 80"–84"–88" expansion
- Handheld backlit pendant control with Auto Contour and cardiac comfort chair
- · Under-mount LED lighting
- Color-coded, dual-pedal floor lock
- · Synchronized Hi/Lo motors
- · Rolls easily at any height
- · Limited Warranty: 15 year frame, lifetime welds, 5 year electronics

#### **OPTIONS**

- · Staff control with Trendelenburg/Reverse Trendelenburg
- · Deluxe assist rails, assist rails, assist bars
- Bed ends available in 16 woodgrain laminate finishes and 2 different styles\* Recommend using Quick Ship OSLO 42"
- · Trapeze, battery backup, bed scale and bed transporter

\*Quick Ship Bed ends available in 3 finishes: OSLO Oak, OSLO Cherry, OSLO Figured Mahogany and available in 2 bed end styles: Brownlee Woods, Mill Creek

#### The **NEW** Drive Medical PrimeCare<sup>™</sup> P750

A full feature, 750 lb. weight capacity, bariatric bed specific to Long Term Care facilities

Now with improved height range, higher net patient weight capacity, one touch central braking system, backlit hand pendant, under bed LED light, cardiac comfort chair and an industry-best FIVE (5) year electronic warranty.



(Shown with optional Staff Control)



Cardiac Comfort Chair



Push Button Width Integrated Length Expansion



Extension



Optional Trapeze



Optional Bed Scale



#### **Expand the Possibilities**

With the unparalleled flexibility of width and length expansion, the P903 Expandable Low Bed can accommodate residents of all shapes and sizes. The P903 simplifies carefully transferring a resident while keeping in mind the safety of your staff. The integrated width expansion adjusts from 36", 39", 42" while the length adjusts to 80", 84" or 88". In its low position, the bed deck rests only 7" inches off the ground and can be raised to a high position of 30" off the ground.





Integrated Width Expansion



Integrated Length Extension



Exceptional Travel Range



Cardiac Comfort Chair



Optional Trapeze

# PRIMECAF Model #P703

#### **FEATURES**

- · 550 lb. weight capacity
- Height adjustable from 7"–30"
- Sleep surface length adjusts to 76", 80" or 84
- · Two-way red/green color-coded floor lock system
- · Eight-function handheld pendant control with auto contour
- 4 DC motors synchronized hi/lo control system keeps bed level
- Sealed bearings secured in "G" track
- IV pole holders
- Extended heel lift
- Eight 3" heavy duty casters roll at any height, caster guides
- Fixed adjustable wall bumper
- · Limited Warranty: 15 year frame, lifetime welds, 5 year electronics

#### **OPTIONS**

- Width expansion kit easily expands from 39"-42"
- Length extension kit easily extend deck length to 88".
- Staff control with Trendelenburg/Reverse Trendelenburg and Cardiac Chair
- Deluxe assist rail, assist rails and assist bars are interchangeable with the P903
- Bed Ends available in 16 woodgrain laminate finishes and 2 different styles
- Trapeze, battery backup and bed transporter

\*Quick Ship Bed Ends available in 3 finishes: OSLO Oak, OSLO Cherry, OSLO Figured Mahogany and available in 2 bed end styles: Brownlee Woods, Mill Creek

#### **Bigger is Not Always Better**

With most of the bells and whistles of the P903, the P703 is reduced in size to accommodate a smaller LTC resident room. The weight capacity is 550 lb. and the length adjusts to 76", 80", 84". In its low position, the bed deck rests only 7" off the ground and can be raised to a high position of 30" off the ground.





550 lb. Weight Capacity



Optional Width Expansion Kit 39"-42"



76". 80". 84" Standard



Length Extension Red/Green Two-way Floor Lock



Optional Trapeze

# PRIMECARE

## Model #P503

#### **FEATURES**

- · 450 lb. weight capacity
- Best-in-class range of motion, 7.8"–30"
- · Integrated 76" and 80" length extension
- Handheld, 8-function hand pendant with Auto Contour
- Eight, 360-degree swivel casters that roll at any height, four casters lock
- · Sturdy, quick-release bed ends
- Sealed bearings provide smooth, reliable and True Vertical Lift operation
- Limited Warranty: 10 year frame, 15 year welds, 2 year electronics

#### **OPTIONS**

- Assist rails and assist bars also fit the P703 and P903
- Bed Ends available in 16 woodgrain laminate finishes and 2 different styles\*

\*Quick Ship Bed Ends available in 3 finishes: OSLO Oak, OSLO Cherry, OSLO Figured Mahogan, and available in 2 bed end styles: Brownlee Woods, Mill Creek

#### Comfort. Safety. Flexibility.

The P503 Low Bed from Drive Medical is an affordable option for facilities wanting to provide residents and staff with a safe and comfortable bed while addressing tight capital budgets.

The contemporary, new design will accommodate your residents with unmatched features and quality not seen in this category of beds.





450 lb. Weight Capacity



Bed Height Range 7.8"-30"



Hand Pendant



Floor Lock and Wall Bumper

#### PrimePLUS™ Bariatric Bed

Model #P2002

This is the ultimate bed! Industry-leading, high-capacity bed has a single-step lock-steer which is ideal for LTC facility staff.



1,000 lb. Weight Capacity



Optional Scale



39", 42", 48", 54" Built In Width Expansion



\*Quick Ship Bed Ends (Sculptured) available in 3 finishes: OSLO Oak, OSLO Cherry, OSLO Figured Mahogany

#### **Features**

- 1000 lb. weight capacity
- 19" low height and 32" high working height
- Width expands from 39", 42", 48" and 54"
- $\cdot$  Sleep surface length adjusts to 80". 84" or 88
- Standard one-button, drop down padded assist rails
- Quiet and easy-to-use hand pendant
- Cardiac comfort chair
- Edema foot ratchet elevates the foot of the bed
- Four heavy-duty locking casters roll easy at any height and lock to secure
- Fixed adjustable wall bumper
- Limited Warranty: 15 year frame, lifetime welds 3 year electronics

#### **Options**

- Embedded staff contro with Trendelenburg/ Reverse Trendelenburg
- Sculptured Bed End available in three woodgrain laminate finishes\*
- Bed accommodates set of padded assist rails
- · Trapeze
- Integral weigh scale (P2002-WS)
- Battery Backup

#### **Accessories**

#### PrimeCare Standard Beds & PrimePlus Bariatric Beds

#### For Long Term Care PrimeCare Standard Beds Only



Rotating Assist Bar Item# PLTCAB, 1/pr

Compatible with the P903, P703, P750 and P503



Rotating Assist Rail Item# PLTCAR, 1/pr

Compatible with the P903, P703, P750 and P503



Deluxe Rotating
Assist Rail
Item# PLTCDR,
1/pr

Compatible with the P903 P703 and P750 Bed



Optional Width Extension Kit Item# P703WE

Compatible with the P703 Bed

## For Long Term Care PrimePlus Bariatric Beds Only Fits P2002 and discontinued P1752 only



Folding Padded Rails

Item# PBASR



### **Create Your Own Style**

#### Headboards and Footboards - Custom

#### **Custom Bed End Options Brownlee Woods**



35BWBE-Custom

35", Brownlee Woods 42BWBE-Custom 42". Brownlee Woods

35" & 42" Custom **Bed Ends** 13 finishes



Kensington Maple



Solar Oak



Wild Cherry



Montana Walnut



Windsor Mahogany



Figured Mahogany



Hemlock



**Pinnacle** Walnut



Veranda Teak



Cordovan



Millwork Cherry



Spiced Cherry



Warm Walnut

#### Mill Creek



35MCBE-Custom 35", Mill Creek 42MCBE-Custom 42", Mill Creek\*

\*42" Mill Creek Bed Ends - note the top arch (t-mold) color will vary from the 13 standard finishes. Vendor will use standard brown.

#### **Never Compromise on Style**

#### Headboards and Footboards - Quick Ship

#### Quick Ship OSLO Beds ONLY Quick Ship (QS) Bed Ends **Brownlee Woods**



35", OSLO Cherry 35BWBE-A 35BWBE-B 35". OSLO Oak

35BWBE-FM 35", OSLO Fig. Mahogany

42BWBE-A 42", OSLO Cherry 42BWBE-B 42", OSLO Oak

42BWBE-FM 42". OSLO Fig. Mahogany







**OSLO** Cherry

OSLO Oak

**OSLO Figured** Mahogany

#### Staff Control Panel (optional)

All collections offer the choice of adding a Staff Control. The staff control allows the caregiver to operate the bed from the foot panel and to individually lock out functions on the handset and/or

staff control as desired.

#### Mill Creek



35MCBE-A 35", OSLO Cherry 35MCBE-B 35", OSLO Oak

35MCBE-FM 35", OSLO Fig. Mahogany

42MCBE-A 42", OSLO Cherry 42MCBE-B 42". OSLO Oak

42MCBE-FM 42", OSLO Fig. Mahogany

#### Bariatric - Available in 38" | Sculptured



PBABE-A Cherry. Use with 1750/52, 2000/02 PBABE-B Oak, Use with 1750/52, 2000/02

PBABE-EM Fig. Mahogany, Use with 1750/52, 2000/02

## GRAVITY™ 9,8&7

## Pressure Redistribution Mattresses

Optional raised rails and multiple sizes available.

#### GRAVITY 9 | Model #15970

- · 450 lb. weight capacity
- 4 Layers of therapeutic foam including full length visco memory foam top layer and heel slope
- Limited Warranty: 15 year mattress, 3 year cover

#### **GRAVITY 8 | Model #15870**


- · 400 lb. weight capacity
- 3 layers of therapeutic foam including full width die cut top layer and heel slope
- Limited Warranty: 10 year mattress, 3 year cover

#### GRAVITY 7 | Model #15770

- · 350 lb. weight capacity
- 3 layers of therapeutic foam including full width die cut top layer
- Limited Warranty: 7 year mattress, 3 year cover

#### Alleviate the Pressure

Designed specifically for pressure redistribution, shear reduction and skin protection, the Gravity series therapeutic foam mattresses include unique combinations of multi-layered, die cut, channel cut, and hybrid visco memory foam.




#### Dynamic Elite™ Pressure Redistribution Long Term Care Mattress

Model #6500-DE-2-FB

#### **Premium Pressure Redistribution**

Optimal pressure redistribution coupled with heat dissipation and increased breathability make the Dynamic Elite the premium option for residents at risk of developing pressure injuries. Gel-infused foam top layer keeps skin cool and dry and enhances patient comfort.



- 400 lb. weight capacity
- 4 Layers of therapeutic foam with full
   width gel-infused memory foam top layer and heel slope
- Fluid- and tear-resistant stretch cover; low shear, easy-to-clean
- $\cdot$  Meets flammability standards (CFR 1632 and CFR 1633)
- Limited Warranty: 15 year mattress, 3 year cover

#### Multi-Ply ShearCare<sup>™</sup> 1500 Pressure Redistribution Foam Mattress

Model #1500SC-42-FB | Model #1500SC-48-FB

#### **Bariatric Pressure Redistribution**

High-density, bariatric, therapeutic foam mattress provides pressure redistribution and maximum support.




- 650 lb (42") 750 lb (48"
- Foam heel section offloads pressure from vulnerable heels
- Removable premium, multi-stretch cover is fluid resistant, low shear and vapor permeable helping protect the skin from friction and moisture.
- Meets flammability standards (CFR 1632 and CFR 1633)
- Limited Warrantv: 5 year mattress, 3 year cover

## Balanced Aire™ Self Adjusting Convertible Mattress

Model #BA9600-NP | BA9600-NP-84 | BA9600-NP-42

Finding the Right Balance for Optimal Pressure Redistribution

Made with revolutionary hybrid visco elastic memory foam, it offers self-adjusting pressure redistribution along with durability and versatility. Combine with a control unit for a dynamic powered air system.



- $\cdot$  500 lb. weight capacity for 35" W | 700 lb. weight capacity for 42" W
- Premium foam filled air cells self adjust to the resident's needs
- Vapor-permeable fluid and tear resistant low shear stretch cover.
- Latex free cover is treated with an antimicrobial to guard against degradation from microorganisms.
- Meets flammability standards (CFR 1632 and CFR 1633)
- Available in various sizes including bariatric
- Limited vvarranty: 5 year mattress, 5 year cover

#### Harmony<sup>™</sup> True Low Air Loss Tri-Therapy Mattress Replacement System

Model #14200N | 14200N-42

#### **Pure Harmony**

Feature-filled True Low Air Loss System for optimal micro-climate control and resident comfort, the intuitive 1300 LPM blower system has alternating pressure and pulsation for the prevention and treatment of any stage pressure injury, including unstageable and deep tissue injuries.



- $\cdot$  550 lb. weight capacity for 36" W
- 650 lb weight capacity for 42" W
- 8" H air cells fully inflate in under two minutes
- PreserveTech™ removable quilted four way stretch cover is manufactured with an anti-microbial\* agent and is fluid resistant, low shear and vapor permeable helping to protect the skin from friction and moisture.
- Limited Warranty: 2 year mattress, 2 year cover, 2 year control unit

## Med-Aire Edge™ Alternating Pressure & Low Air Loss Mattress Replacement System

Model #14360-B | 14360XP

#### **Therapeutic Power**

Experience the therapeutic power of air along with the safety and transfer ease of foam with Drive Medical's new Med-Aire Edge!

Designed to prevent, treat and heal pressure injuries in the home or long term care setting.



- 450 lb weight capacity
- CPR Valve allows for rapid deflation
- 4-way stretch quilted mattress cover is fluid-resistant, low shear and vapor permeable, protecting the skin from friction and moisture. Cover is treated with an antimicrobial to guard against degradation from microorganisms.
- 20 individual air cells to maximize immersion and envelopment, including 9 ventilated air cells for micro-climate control
- 5" deep air cells are easily removed with quick disconnects for convenient servicing
- $\cdot$  3" firm foam perimeter for safe ingress, egress and patient transfers
- $\cdot$  3" firm foam base for protection during power outages
- Limited Warranty: 18 month mattress, 18 month cover, 18 month control unit

## Med-Aire™ Plus 8" Alternating Pressure and Low Air Loss Mattress System

Model #14029 | 14029DP | 14029-84

#### **Drive Medical's Most Popular Powered Air Support Surface**

Combination therapy system provides both alternating pressure and low air loss to optimize pressure redistribution and manage skin maceration for the prevention and treatment of any stage pressure injury. Cell-on-cell design of 8" deep air cells prevents "bottoming out" and provides at least 24 hours of power outage protection.



- 450 lb. weight capacity
- Multi-stretch, fluid-resistant, vapor-permeable, low shear, quilted cover
- $\cdot$  Options available with 10" deep defined perimeters or 84" length
- Advanced clinical and caregiver features: Auto-firm, seat inflate, microprocessing feedback, lock-out and more!
- Limited Warranty: 18 months mattress, 18 months cover, 18 months control unit

## Med-Aire 8" Alternating Pressure and Low Air Loss Mattress System

Model #14027

#### **Just the Right Pressure**

Combination therapy 8" mattress replacement system provides both alternating pressure and low air loss to prevent, treat and manage any stage pressure injury.



- 350 lb weight capacity
- Multi-stretch, fluid-resistant, vapor-permeable, low shear, quilted cover
- Adjustable patient comfort settings allow for optimal immersion and patient satisfaction
- Limited Warranty: 1 year mattress, 1 year cover, 1 year control unit



#### **INTRODUCING THE BRAND NEW**

# COMMON ROOM & DINING ROOM COLLECTIONS Basic & Deluxe Options Available

Elevate your facility with complete furniture solutions for senior living

#### **KEY FEATURES & BENEFITS**

- Highly Customizable: Tailor your selection to your facility's needs and budget with seven fabric grades, more than 2,000 unique patterns, and customizable finish options
- Durable Construction: Commercial-grade materials hold up under heavy use, reducing replacements and maintenance costs
- Ergonomic Comfort: Designed to support resident mobility and seating needs
- Safety-Focused Design: Rounded edges and stable bases help reduce the risk of accidents
- Low-Maintenance Materials: Smooth surfaces and resilient upholstery simplify cleaning and care
- Made in America Quality: Local manufacturing supports consistent quality and supply reliability

#### Create a Space That Feels Like Home

Residents notice the difference. So will the families who visit your facility.

The new, premium deluxe and performance-focused basic collections transform common and dining areas into spaces that feel warm and inviting, just like home. With seven fabric grades, thousands of patterns to choose from, and customizable wood finishes, each piece can be tailored to complement any decor and balance durability and budget for your facility.

Contact us today to discuss your custom furniture request or to learn about quick-ship options!



#### **Dining Room Furniture Collection**

#### Bringing residents together, one meal at a time.

Transform daily meals into dining experiences. Our chairs, tables, and accessories are built to be accessible, withstand daily use, and simplify cleaning and storage. Custom fabrics and finishes allow you to create dining spaces that encourage conversation and connection while matching your facility's style and budget.



#### **Basic Dining Collection**

Durable, practical furniture built for busy dining areas.

#### **Basic Dining Chairs**

DC001	Basic Dining Chair			
Basic Dining Table Tops Available in Square or Round				
DT001-36S	Basic Dining Table: 36" Countertop			
DT001-42S	Basic Dining Table: 42" Countertop			
DT001-48S	Basic Dining Table: 48" Countertop			
Basic Dining Table Bases				
TB001	Basic Table Base			
TB001-W	Basic Table Base, Wheelchair Height			
TB001-36	Basic 36" Table Base			
TB001-36-W	Basic 36" Table Base - Wheelchair Height			

Chairs & Basic Dining Table Finish Options



Warm Walnut



Millwork Cherry

#### **Deluxe Dining Collection**

Elegantly designed with premium fabrics, high-performance materials, and convenient features for movement and storage.

#### **Deluxe Dining Chairs**

5				
DC002	Deluxe Dining Chair			
Deluxe Dining Table Tops + Bases				
DT002-S-1-BK	Deluxe Dining Table - Square, 1 Castor, Black			
DT002-S-4-BK	Deluxe Dining Table - Square, 4 Castor, Black			
DT002-S-1-BN	Deluxe Dining Table - Square, 1 Castor, Brown			
DT002-S-4-BN	Deluxe Dining Table - Square, 4 Castor, Brown			
DT002-R-1-BK	Deluxe Dining Table - Round, 1 Castor, Black			
DT002-R-4-BK	Deluxe Dining Table - Round, 4 Castor, Black			
DT002-R-1-BN	Deluxe Dining Table - Round, 1 Castor, Brown			
DT002-R-4-BN	Deluxe Dining Table - Round, 4 Castor, Brown			

#### Deluxe Dining Table Finish Options



Spice Walnut



Millwork Cherry

#### **Common Room Furniture Collection**

#### Spaces that invite residents to gather and socialize

Create welcoming areas where residents can relax, interact, and form connections. Our common room furniture combines comfort, durability, and easy maintenance to support both residents and staff.



#### **Basic Dining Collection**

High-performance, accessible furniture that invites every resident to sit down and stay awhile

#### Basic Sofas & Love Seats

BSF001	Basic Sofa			
BLS001	Basic Loveseat			
Basic Lounge Chairs				
BLC001	Basic Lounge Chair			
Basic Coffee & End Tables				
ET001-T	Basic End Table - Tall			
ETOO1-W	Basic Coffee Table - Wide			

#### **Deluxe Dining Collection**

Sophisticated seating with premium fabric options and refined finishes, bringing style and comfort to shared spaces that foster socialization

#### Deluxe Sofas & Love Seats

Delake Solas & Love Seats				
DSF001	Deluxe Sofa			
DLS001	Deluxe Loveseat			
Deluxe Lounge Chairs				
DLC001	Deluxe Lounge Chair			
Deluxe Coffee & End Tables				
ET002-S	Deluxe End Table - Square			
ET002-R	Deluxe Coffee Table - Rectangle			

#### Most Popular Upholstery Options



Burch Acadia Aloe



Burch Mimic Cinder



Burch Acadia Silt

#### 1,000s of Options Available

# OSLO OLLECTION

BEDSIDE CABINETS, CHESTS & WARDROBES



- Precision-engineered full extension drawer glides and rubber stops \ reduce jarring to provide quiet closure
- European-style hinges are hidden from view while providing ease of operation, 110° accessible opening and self-closing feature
- Door and drawer handles are easy to grasp and smooth to the touch, have an attractive contemporary design and quality brushed nickel finish
- Ventilated back allows air flow throughout the cabinet interior
- Removable drawer fronts on cabinets and chests are easy to replace if damaged, helping to ensure years of service
- All exposed interior surfaces are finished with matching color for aesthetic appeal and ease of cleaning. Finished edges and sealed bottom help keep moisture out
- Durable plastic floor glides facilitate ease of movement and resist marring

#### **Best in Class & Value**

## The Oslo Collection is a stylish, high-quality line of LTC furniture that will help your residents feel right at home.

Offering a full suite of bedside cabinets, chests and wardrobes, the Oslo Collection is equipped with high-pressure laminate top and drawer face, matching melamine sides, drawer interior and ventilated back. Also included are precision full extension drawer glides, self closing Euro hinges.



## Available Finish Options



OSLO Cherry



OSLO Oak



OSLO Figured Mahogany

#### Chests

NOC3D-C	3 Drawer, Cherry		
NOC3D-O	3 Drawer, Oak		
NOC3D-FM	3 Drawer, Figured Mahogany		
NOC4D-C 4 Drawer, Cherry			
NOC4D-O 4 Drawer, Oak			
NOC4D-FM 4 Drawer, Figured, Mahogany			

Bedside Cabinets				

# FURNITURE SUITES RIVERSIDE, MILL CREEK,

RIVERSIDE, MILL CREEK & BROWNLEE WOODS,

#### **FEATURES**

- · Drawer slides are epoxy coated and have a 100 lb. rating
- Four-sided, white vinyl wrapped drawer boxes make it easier for resident to identify personal property
- Ventilated backs facilitate air tirculation to help prevent bacteria growth and reduce odor

#### **OPTIONS**

- Choice of 8 popular woodgrain laminate finishes complement most any décor. Seven custom finish options available with specific quantities and pricing.
- Option of 9 fixed drawer/door handles are beautifully crafted and easy to grasp
- Full extension drawer glides
- · Drawer locks

#### First Impressions Leave a Lasting Impression

An LTC facility is home to your residents. Welcoming them with an inviting bedroom, styled with a safe and functional bed, contemporary furniture and a common area, furnished in a way that evokes the feelings of home, will create a warm environment for long term care.

Contact your local sales representative about our new and improved wardrobe collection



#### **Standard Finish Options**



#### Additional Custom Finish Options\*

\*Custom colors available with specific quantities & pricing. Contact your Sales Representative for more information



#### **Riverside Collection**

This premium furniture line has solid hardwood face frames that reflect the mission style.



#### Bedside Cabinets and Chests

#### Wardrobes

Model #	Description	Model #	Description
RS-BDSC-11	1 Door 1 Drawer Bedside	RS-WDRB-24-10	1 Door Wardrobe
RS-BDSC-03	3 Drawer Bedside	RS-WDRB-24-12	1 Door, 2 Drawer Wardrobe
RS-CHST-03	3 Drawer Chest		
RS-CHST-03W	3 Drawer Chest Wide	RS-WDRB-30-20	2 Door Wardrobe
R3-CH31-U3VV	3 Diawei Criest Wide	DC WDDD 7C 20	2 Door Wardrobe Wide
RS-CHST-04	4 Drawer Chest	R3-VVDRB-30-20	2 Door Wardrobe Wide
	- Drawer chest		2 Door, 2 Drawer Wardrobe
RS-CHST-06	6 Drawer Chest	RS-WDRB-30-22	
RS-HUTC-1	Bedside Hutch	RS-WDRB-36-22	2 Door, 2 Drawer Wardrobe Wide

- Tops feature hardwood edge and high-pressure laminate for maximum durability
- Sides are high-pressure thermofused
   laminate with PVC Edgeband and hardwood face frame for enhanced appeal
- Drawers and doors are high-pressure thermofused laminate

# Mill Creek Collection

This premium furniture line features solid wood edge drawer fronts and tops that will beautifully complement any Long Term Care Facility.



# Bedside Cabinets and Chests

# Wardrobes

Model #	Description	Model #	Description
MC-BDSC-11	1 Door, 1 Drawer Bedside	MC-WDRB-10	1 Door, No Drawer
MC-BDSC-03	3 Drawer Bedside	-	Wardrobe
MC-CHST-03	3 Drawer Chest	MC-WDRB-12	1 Door, 2 Drawer Wardrobe
MC-CHST-03W	3 Drawer Chest X Wide	MC M/DDD 22	2 Door, 2 Drawer
MC-CHST-04	4 Drawer Chest	MC-WDRB-22	Wardrobe
MC-CHST-06	6 Drawer Chest (3 x 3)	MC-WDRB-20W	2 Door, No Drawer Wardrobe Wide
MC-HUTC-1	Hutch - Bedside	MC-WDRB-22W	2 Door, 2 Drawer Wardrobe Wide

- Tops feature high-pressure laminate for maximum durability with solid wood edges for enhanced appeal; cabinet finish is high-pressure thermofused laminate
- Drawers and doors feature solid wood edges

# **Brownlee Woods Collection**

Featuring a high-pressure laminate top with a T-mold leading edge, the Brownlee Woods Collection has a luxurious look at an affordable price.



#### **Bedside Cabinets and Chests**

#### Wardrobes

Model #	Description	Model #	Description
BW-BDSC-11	1 Door, 1 Drawer Bedside	BW-WDRB-10	1 Door, No Drawer Wardrobe
BW-BDSC-03	3 Drawer Bedside	BW-WDRB-12	1 Door, 2 Drawer Wardrobe
BW-CHST-03	3 Drawer Chest	BW-WDRB-22	2 Door, 2 Drawer Wardrobe
BW-CHST-03W	3 Drawer Chest X Wide	BW-WDRB-20W	2 Door, No Drawer
BW-CHST-04	4 Drawer Chest		Wardrobe Wide
BW-CHST-06	6 Drawer Chest (3 x 3)	BW-WDRB-22W	2 Door, 2 Drawer Wardrobe Wide
BW-HUTC-1	Hutch 21" Bedside	RS-WDRB-36-22	2 Door, 2 Drawer Wardrobe Wide

- Tops and sides feature high-pressure thermofused laminate
- Drawers and doors with PVC protective edge

# **Furniture Options**

# Standard Features, Options and Hardware

# Standard Features





Standard Drawer Slide



Euro Door Hinge

# **Options**



Upgraded Drawer Slide



Drawer Lock

Interior Finish

# Hardware



Bronze Pull, Oil Rubbed



Revere Pull. Brushed Nickel



Arc Pull



Aspire Pull, Brushed Nickel



Bar Pull. Stainless Steel



Modern Pull. Brushed Nickel



Hutton Pull. Nickel



X Pull. Satin Nickel



ECP Pull. Antique Brass



Straight Toe Kick



Crescent Toe Kick



Scalloped Toe Kick

# **Overbed Tables**

# Choices made simple.

Colors, contours, and cutouts... Select the features and finishes that satisfy your particular requirements throughout the continuum of care. A wide array of custom finishes and designs, colors and styles, lets you complement your room décor for a coordinated aesthetic look.

Contact us today with special requests and ideas.



Model #13003 Non-Tilt Over Bed Table with H base

Model #13067 Non-Tilt Overbed Table with H base

Model #13069 Overbed Table with Plastic Top and H base

# Model #13003

- Walnut wood grain tabletop locks securely when height adjustment handle is released
- It can be raised with the slightest upward pressure
- Chrome base
- 2" swivel casters for easy maneuverability

# Model #13067

- Walnut wood grain low-pressure laminate top
- Table top can be raised or lowered between 28"-45"
- Silver vein "H" easy-tomaintain and less institutional looking

#### Model #13069

- Tabletop can be raised or lowered between 29"–42" with light upward pressure
- Tabletop locks securely when height adjustment handle is released
- Swivel casters for easy maneuverability
- Quick and easy assembly

# **Custom Order Requests**

Custom order requests can include changes in:

Color of base, shape of top, spill containment options, special order laminates & special order thermoform.



# **Available Stock Pressure Laminates**



# Model #16012

- Spring loaded lift mechanism positions top in 1" increments from 28"-45"
- 3/4" x 15" x 30" high-pressure laminate top with T-molding
- Available with and without vanity
- Available with "H" or "U" base in chrome or opal powder coat

# PrimeMat™ 2.0 Impact Reduction Fall Mats

Model # PM20BA (brown) | PM20GA (gray)

# **Fall Mitigation**

The PrimeMat<sup>™</sup> 2.0 is the only choice when developing a comprehensive fall prevention program. The molded polyurethane foam inner core construction provides maximum absorption and reduces the impact in the event of a fall.



# √ Egg Drop Tested

In which an egg is dropped from 5 ft above the mat without damage

- Full 72" length provides maximum coverage alongside bed
- Passed 5' "Egg Drop Test"—an egg was dropped from 5' above the mat without damage
- Low-profile with tapered edges allows equipment such as overbed tables or patient lifts to roll over it easily
- Anti-slip bottom netting keeps the mat secure to help prevent accidental tripping
- Available in brown and gray
- · Limited Warranty: 1 year

# **Resident Lifts**

Available with Manual or Power Base (shown)

#### Introducing a suite of resident lifts

Contemporary soft curve, stylish and equipped with innovative features, such as Smart & Safety™ Technology, these floor lifts provide an efficient, cost-effective solution for resident care.



Model #	Description
FLNP500	$Levantar^{\text{TM}} \ Floor \ Lift, \ Battery-Powered, \ Manual \ Base, 500 \ lb., 6 \ point \ cradle$
FLP500	Levantar Floor Lift, Battery-Powered, Power Base, 500 lb., 6 point cradle
FLNP600	Gravis <sup>™</sup> Floor Lift, Battery-Powered, Manual Base, 600 lb., 6 point cradle
FLP600	Gravis Floor Lift, Battery-Powered, Power Base, 600 lb., 6 point cradle

- Lift control module and battery are smaller and complement design of the lift
- Smart & Safety Technology
  - Redundant lift controls and visual battery fuel indicator
  - Audible alert battery fuel indicator
  - O.O.S. (Object Obstruction Sensor) stops lift if resident or caregiver accidentally gets in the way of the boom arm.
  - Audible alert if the lift detects weight overload
- Maintenance monitoring:
  - Lift Control Station displays lift utilization and service data. Will also show usage history and has a service reminder to keep you on track for compliance.
  - Utilizing data is key in facilities who look for solutions in meeting compliance with their Safe Handling Programs
- Limited Warranty: Lifetime frame, 5 year electrical, 1 year battery

# Sit to Stand Lifts

Available with Manual or Power Base (shown)

#### Introducing new Sit to Stand lifts

An elegant, contemporary design with soft lines. Packed full of smart and safety features typically only seen on premium Sit to Stand lifts, BUT at an entry-level price point.



Model #	Description
STSP450	Sit to Stand Lift with Power Base
STSM450	Sit to Stand Lift with Manual Base

- 450 lb. Weight Capacity
- Best in Class Limited Warranty
- Dual controls never have a down lift because of inoperable hand pendant
- Audible and visual alert if exceeds weight capacity of 450 lb
- Object Obstruction Sensor lift stops if detects patient or caregiver is in the way
- Adjustable knee pad
- Removable foot plate
- Lowest floor clearance in class
- Vertical and horizontal position ergonomic handles
- Battery can be charged on and off the lift
- Audible/visual alert when battery is running low
- STSP450 power base includes LCD screen for visual read out and maintenance monitoring
- Limited Warranty: Lifetime frame, 5 year electrical, 1 year battery





Full Body Sling with Commode Opening



Divided Leg Sling



Sit-to-Stand Sling

# **Available Slings**

Model #	Description	Model #	Description
13220L	U-Sling 4 Point Padded, Large	13222M	Full Body, Solid, Medium
13220M	U-Sling 4 Point Padded, Medium	13262L	Divided Leg Sling with Headrest, Large
13220S	U-Sling 4 Point Padded, Small	13262M	Divided Leg Sling with Headrest, Medium
13220XL	U-Sling 4 Point Padded, XL	13263A	Patient Stand-Up Sling, Large
13221L	Commode Opening, Mesh, Large	13263B	Patient Stand-Up Sling, XL
13221M	Commode Opening, Mesh, Medium	13263E	Patient Sit-to-Stand Sling, Medium
13221XL	Commode Opening, Mesh, XL	13263F	Patient Sit-to-Stand Sling, Large
13221XXL	Full Body, Mesh, XXL Large	13263H	Patient Sit-to-Stand-Sling, Strap
13222L	Full Body, Solid, Large	13264A	Patient Transport Sling

# **Resident Room Recliners**

# Recovery Made Easy

Support residents as they recover with clinically designed, comfortable recliners.

# **3-Position Recliners**

Model #	Description	
D574-CHOC	Chocolate	
D574-CHAR	Charcoal	
D574-TAN	Tan	

# **3-Position Bariatric Recliners**

Model #	Description	
D574EW-CHOC	Chocolate	
D574EW-CHAR	Charcoal	
D574EW-TAN	Tan	

#### **4-Position Recliners**

Model #	Description
D577-CHOC	Chocolate
D577-CHAR	Charcoal
D577-TAN	Tan







Chocolate Charcoal

Tan







Recline



Full Recline



Trendelenburg



# Silver Sport Wheelchair

# Make Every Second of Every Day Count for Your Residents

Designed for easy transport and user comfort. It features a sturdy steel frame with a powder-coated finish, offering long-lasting strength and resistance to wear. With its smooth-rolling wheels and rear wheel locks, it ensures both maneuverability and safety, making it an ideal choice for individuals seeking reliable daily mobility support.



# **Features**

- 250 lb. weight capacity
- Padded fixed armrests
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Carry pocket on backrest provides additional convenience
- · Comes with swing-away footrests

# **Cushion Features**

- Cover is water-resistant, lowshear stretch to reduce heat and moisture build-up
- Convenient 12" tie straps secure the cushion to the wheelchair
- Meets flammability standard TR117
- · Additional sizes available
- Limited Warranty: 18 Month

# Viper Plus GT Wheelchair with Universal Armrests

# **Superior Customization**

Allows the user to easily convert the armrests into desk arm, half arm, and full arm with the push of a button. The armrests also adjust to eight individual heights in ½" increments for complete customization.



#### **Features**

- 450 lb. weight capacity
- Adjustable height and length, flip-back, removable, padded arms
- Built-in seat rail extensions and extendable upholstery easily adjusts seat depth from 16" to 18"
- Adjustable angle back easily adjusts from 5° – 20°
- Adjustable height back adjusts from 17" 19" in 1" increments, adaptable for individual patient needs

# **Cushion Features**

- Available in 4" and 2" adjustable air cell heights and various sizes
- Machine washable, low-shear stretch cover with adjustable locking straps and non-skid bottom
- Hand pump included to adjust the cushion to individual patient comfort
- Meets flammability standard TB117
- Limited Warrantv: 2 Year

# iGo®2 Portable Oxygen Concentrator with Optional Bluetooth Connectivity

Model# 125D | # 125D-BT

Oxygen Optimized, It's the peace of mind you have been waiting for!

With patented SmartDose® Technology, the new iGo®2 automatically adjusts oxygen production to match a patient's oxygen needs based on breath rate. Now your patients can enjoy life more freely without worrying about changing oxygen settings with activity, and you'll have peace of mind knowing your patients are receiving clinically-effective oxygen therapy no matter the location, time of day or level of activity.



- Revolutionary SmartDose Technology automatically adjusts oxygen production & delivery in real time to meet patient needs
- Unmatched conserver trigger sensitivity (0.05 cm H2O)
- $\cdot$  Lightweight and easy to carry in its convenient over-the-shoulder carry bag
- Durable, high-quality design
- Bluetooth connectivity (125D-BT only) helps reduce costs and increase operational efficiency with Remote Fleet Monitoring

# **5L and 10L Oxygen Concentrator**

Model #525DS | Model #1025DS

#### **Breathe More Easily & Effectively**

Designed with enhanced quality, durability and simplicity in mind, the 5 Liter Oxygen Concentrator is built upon a legacy of reliable, field-proven concentrator models.

The 10 Liter Oxygen Concentrator is designed to provide optimal oxygen delivery for residents with varying flow requirements. With a weight of 42 lb. and an identical footprint to the 5 Liter Concentrator, the 10 Liter Concentrator is substantially smaller and lighter than competitive 10 liter concentrators.





#### **Features**

#### 5 Liter

- Readily accessible controls, protected cannula fitting and recessed humidifier nook to prevent damage
- Alarms audible and visual high/ low pressure, low flow, low oxygen, power fail, OSD® (Oxygen Sensing Device)
- Transfill Caddy sold separately

#### 10 Liter

- Capable of delivering 87% to 96% of O2 purity from 2 to 10 LPM
- Equipped with accessible controls, bright LEDs that can be seen from a distance, protected cannula fitting and recessed humidifier nook to prevent damage
- Transfill Caddy sold separately

# Vacu-Aide® Suction Machine

Model #18600

# **Powerful Suctioning**

This powerful and durable suction pump combines the performance expected for an institution with the ease of operation needed for use in long term care.



- Includes 800 cc suction canister, 6' suction tube, 10" suction canister tubing, hydrophobic filter, plastic elbow connector and manual
- Vacuum levels up to 560 mmHg

# Vacu-Aide® Suction Machine

Model #7325P-D

# **Lightweight & Lithium-Ion Battery Power Option**

Quiet, powerful and easy to use, the Vacu-Aide® 7325P-D features a long-lasting lithium-ion battery and weighs less than 4 pounds.



- Compact footprint enables portability while using standard canisters
- Ultra-quiet for peaceful treatment at ~59 dBA
- Lightweight, rechargeable lithium-ion battery provides lasting power and includes an industry-leading one-year warranty
- 2-Year Limited Warranty on Suction Unit

# BATHING SYSTEMS

COLORADO, SARATOGA & NIAGARA TUBS

# Colorado Spring™ Bathing System Model #DB-COLORADO-1100

The Drive Medical Colorado is a free-standing tub that allows the caregiver 360° unobstructed access when bathing a resident. The low threshold and wide side-entry door allow easy access for residents with various degrees of mobility, and the built-in contoured seat allows the resident to sit comfortably and safely in a slightly reclined position.

# Saratoga Spring™ Bathing System Model #DB-SARATOGA-1100

The Drive Medical Saratoga is designed for healthcare facilities looking for a more time efficient and labor-saving bathing solution. The Saratoga has a unique tilting function which allows residents with challenging conditions to recline in a contoured seat and quickly become immersed in warm water.



- · All tubs are 100% made in the USA
- · Best in class limited warranty
- · Thermostatic mixing valve with anti-scald protection
- · Built-in Disinfection System and freestanding tub design
- Limited Warranty: Lifetime door seal, 5 year external finishes and composite shell, 3 year moving and stationary parts

# Niagara Spring™ Bathing System Model #DB-Niagara-1500

The Drive Medical Niagara Tub is designed to bathe highly dependent residents safely and comfortably, and is compatible with a wide range of transport and lifting devices. The key-hole shape enhances patient comfort with extra room for the shoulders and arms, plus allows the caregiver closer and greater access.





# www.drivemedical.com

©2025 Medical Depot, Inc. dba Drive Medical, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.