

**COMPASSHEALTH**

*Improving Lives*



# AEROSOL THERAPY/ SUCTION CATALOG



800.376.7263 | [www.compasshealthbrands.com](http://www.compasshealthbrands.com)



## Nebulizers

We take great pride in offering you a wide selection of products to meet all your patient's unique needs. Our high-quality, economical aerosol compressors feature sleek designs, durable shells and built-in handles that are easy to pick up and go.

### AEROSOL THERAPY

Nebulizer Compressor Systems	3
Portable Nebulizer	4
Pediatric Nebulizer Kits	5-9
Nebulizer Kits	10-11
Aerosol Chambers & Masks	11
Aerosol Supplies & Accessories	12-15

### SUCTION THERAPY

Aspirators	16
Portable Suction	17
Aerosol Compressors	18-19
Suction Canisters & Kits	20
Suction Filters & Accessories	21-23



## Rite-Neb 5 NEBULIZER COMPRESSOR SYSTEM

The **Roscoe RITE-NEB5 Nebulizer System** offers reliability, high performance and convenient size with a convenient carry handle. The simple one-button operation and lightweight design makes it extremely portable and patient friendly. Weighing under 3 lbs, the Rite-Neb 4 includes the compressor, nebulizer, tee, mouthpiece, corrugated reservoir tube, tubing, five extra filters, instruction manual and a five year warranty. The Roscoe RITE-NEB5 is also available with a reusable neb kit option.



### FEATURES:

- One-button operation
- Delivers better output for faster treatments and better patient compliance
- Compact design makes it portable
- Quiet at less than 65 dBA
- Compact and lightweight at 2.7 lbs
- Metal piston for durability
- Attractive retail packaging
- Comes complete with compressor, nebulizer, tee, mouthpiece, reservoir tube, tubing and five extra filters
- TruNeb Reusable Neb Kits available – purchase with the compressor or alone
- HCPCS Code: E0570\*
- HCPCS Code: A7003\* for nebulizer kits
- HCPCS Code: A7005\* for reusable nebulizer kits
- HCPCS Code: A7013\* for felt filters

### PRODUCT SPECIFICATIONS:

Dimensions:	6.1" L x 5.5" W x 3.5" H
Weight:	2.7 lbs
Method of Nebulization:	Piston compressor
Maximum Pressure:	30 to 50 psi
Operating Pressure Range:	9.5 PSI - 16 PSI
Maximum Flow:	6 - 10 LPM
Sound Level:	<65 dBA
Medication Capacity:	5 mL
Particle Size:	< 0.5 µm
Electrical Requirements:	120 VAC/60Hz
Operation Cycle:	20 min. on, 30 min. off
Warranty:	Limited 5 year

### OPTIONAL ACCESSORIES:



Item #	Description	Qty
RITE-NEB5	RITE-NEB5 Nebulizer Compressor System with Disposable Neb Kit	each
RITE-NEB5-REUS	RITE-NEB5 Nebulizer Compressor System with Reusable Neb Kit	each
RITE-NEB5-REUSPM-PB	RITE-NEB5 Nebulizer Compressor System with one disposable neb kit, one reusable neb kit, one pediatric mask with piggyback serial labels	each
Optional Accessories:		
CIF-1128PK6	Foam Filter, 5/8" diameter x 1/2" thick	6/pk
50016	TruNeb Reusable Nebulizer Kit	10/cs
BAG-NEBR	Nylon Carry Bag	each
BAG-DRW	Nylon Drawstring Carry Bag	each

\* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

# nebGO Portable Handheld Nebulizer

## ULTRASONIC MESH NEBULIZER



The Roscoe nebGO will be your portable nebulizer solution for on-the-go aerosol treatments. With a small, lightweight design that comfortably fits in the palm of your hand, it delivers powerful, quiet, and discreet neb treatments anywhere using vibrating mesh technology. Weighing in at just 3.7 ounces (with batteries), the nebGO provides quick, effective treatments without the bulk of a standard tabletop compressor. The system includes the main unit, medication cup, mouthpiece, adult mask, pediatric mask, lanyard, USB cord, and power plug, stylish hard shell carrying case, and user manual.



### PRODUCT SPECIFICATIONS:

Dimensions:	3" L x 1.6" W x 2.9" H (7.65cm x 4.1cm x 7.3cm)
Weight:	2.6 oz/74g (without batteries) 3.7 oz/104.9g (with batteries)
Particle Size:	MMAD < 5µm
Average Nebulization Rate:	≥0.25mL/min
Operating Cycle:	20 minutes on, 20 minutes off
Battery Life:	Up to 90 minutes continuous use; 4 days if used 20 min/day
Operating Temperature/ Humidity Range:	50°F to 104°F (10°C to 40°C) / 30-75%RH
Storage Temperature/ Humidity Range:	-4°F to 158°F (-20°C to 70°C) / 20-80%RH
Power:	Two AA Alkaline Batteries; USB Power Adapter
Power Input:	AC100-240V/50-60Hz
Power Consumption:	<1.5W
Medication Capacity:	6mL
Vibrating Frequency:	117kHz (+/- 15%)
Warranty:	Two Year Limited Warranty
HCPCS:	E0574

### ADDITIONAL ACCESSORIES:



Item #	Description	UPC	UOM
NEB-GO	nebGO Ultrasonic Handheld Nebulizer	092237621738	each
<b>Accessories:</b>			
X1GOSH01	nebGO Medication Cup Replacement	092237621745	each
X1GOSH02	nebGO Mouthpiece Replacement	092237621752	each
X1GOSH03	nebGO Adult Mask Replacement	092237621769	each
X1GOSH04	nebGO Pediatric Mask Replacement	092237621776	each
X1GOSH05	nebGO Carry Case Replacement	092237621783	each
X1GOSH06	nebGO Power Cord Replacement	092237621790	each
X1GOSH07	nebGO Wall Outlet Replacement	092237621806	each

\* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

# Respiratory – Aerosol Therapy

## DINOSAUR PEDIATRIC NEBULIZER SYSTEM



Designed to reduce anxiety during treatments and encourage compliance, the child-friendly **Roscoe Dinosaur Compressor Nebulizer** provides effective aerosol treatments with a state-of-the-art piston compressor. Delivers 42 PSI at maximum pressure and 13 PSI operating pressure (liter flow range 5 - 8 LPM). Equipped with a convenient nebulizer parking post. Comes complete with compressor, disposable nebulizer, child mask, angled mouthpiece, 7-foot kink-resistant tubing, five extra filters and a five year warranty. Also available with a reusable neb kit option (Model No. DINO-TRU).



### FEATURES

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- Nylon carry bag available with BAG-DINO
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone
- DINO-TRU includes TruNeb Reusable Neb Kit
- HCPCS\* Code: E0570\* for nebulizer with compressor  
 A7003\* for disposable nebulizer kits  
 A7005\* for reusable nebulizer kits  
 A7013\* for foam filters  
 A7015\* for aerosol mask

### SPECIFICATIONS

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	0.5 µm to 10 µm
Unit Dimensions:	6" (L) x 6" (W) x 9" (H)
Product Weight:	3.4 lbs.
Maximum Pressure:	42 PSI
Operating Pressure:	13 PSI
Liter Flow:	8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

Item #	Description	UOM
NEB-DINO	Dinosaur Pediatric Nebulizer System with one Disposable Neb Kkit	each
DINO-TRU	Dinosaur Pediatric Nebulizer System with one Disposable Neb Kit, one TruNeb Reusable Neb Kit	each
DINO-TRUWB-PB	Dinosaur Pediatric Nebulizer System with one Disposable Neb Kit, one Reusable Neb Kit and Black Nylon Carry Bag with Piggyback Serial Labels	each
<b>Accessories</b>		
BAG-DINO	Dinosaur Nylon Carry Bag	each
50016	TruNeb Reusable Neb kit with Angle Mouthpiece and Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 in. Diameter x 1/2 in. Thick	6/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

\*\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

# Respiratory - Aerosol Therapy

## BUNNY PEDIATRIC NEBULIZER SYSTEM



Designed to reduce anxiety during treatments, the child-friendly **Roscoe Bunny Nebulizer** provides effective treatments with a state-of-the-art piston compressor. Delivers 42 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty. The Roscoe Bunny Nebulizer Compressor is also available with the TruNeb Reusable Neb Kit option.



### SPECIFICATIONS

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	Less than 0.5 µm
Unit Dimensions:	8 4/5" (L) x 5 2/5" (W) x 9" (H)
Product Weight:	3.4 lbs.
Maximum Pressure:	42 PSI
Operating Pressure:	16 PSI
Liter Flow:	5 - 8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

### FEATURES

- Child-friendly design helps to reduce anxiety during treatments
  - Provides smooth, quiet treatments
  - Equipped with a convenient nebulizer parking post
  - Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
  - Nylon carry bag available BAG-BUNNY
  - TruNeb Reusable Neb Kit available for purchase with the compressor or alone
  - BUNNY-TRU includes TruNeb Reusable Neb Kit
- HCPCS\* Code: E0570\* for nebulizer with compressor  
 A7003\* for disposable nebulizer kits  
 A7005\* for reusable nebulizer kits  
 A7013\* for foam filters  
 A7015\* for aerosol mask

Item #	Description	UOM
NEB-BUNNY	Bunny Pediatric Nebulizer System with one Disposable Neb Kit	each
BUNNY-TRU	Bunny Pediatric Nebulizer System with one Disposable Neb Kit, one Truneb Reusable Neb Kit	each
BUNNY-TRUWB-PB	Bunny Pediatric Nebulizer System with one Disposable Neb Kit, one Reusable Neb Kit and Black Nylon Carry Bag with Piggyback Serial Labels	each
<b>Accessories</b>		
BAG-BUNNY	Bunny Nylon Carry Bag	each
50016	TruNeb reusable neb kit with Angle Mouthpiece and Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 in. Diameter x 1/2 in. Thick	6/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.  
 \*\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

# Respiratory - Aerosol Therapy

## FROG PEDIATRIC NEBULIZER SYSTEM



Designed to reduce anxiety during treatments, the child-friendly **Roscoe Pediatric Frog Nebulizer** provides effective treatments with a state-of-the-art piston compressor. Delivers 42 psi at maximum pressure and 16 psi operating pressure (liter flow range 5 - 8 LPM) and weighs only 3.4 pounds. Equipped with a built-in carrying handle and a convenient nebulizer parking post. Comes complete with compressor, nebulizer, child mask, tubing, angled mouthpiece, five extra filters and a five year warranty. The Roscoe Frog Nebulizer Compressor is also available with the TruNeb Reusable Neb Kit option.



### SPECIFICATIONS

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	Less than 0.5 µm
Unit Dimensions:	7 5/8" (L) x 6 1/4" (W) x 5 1/3" (H)
Product Weight:	3.4 lbs.
Maximum Pressure:	42 PSI
Operating Pressure:	16 PSI
Liter Flow:	5 - 8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

### FEATURES

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- Nylon carry bag available 50002
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone
- FROG-TRU includes TruNeb Reusable Neb Kit
- HCPCS\* Code: E0570\* for nebulizer with compressor  
 A7003\* for disposable nebulizer kits  
 A7005\* for reusable nebulizer kits  
 A7013\* for foam filters  
 A7015\* for aerosol mask

Item #	Description	UOM
50000	Frog Pediatric Nebulizer System with one Disposable Neb Kit	each
FROG-TRU	Frog Pediatric Nebulizer System with one Disposable Neb Kit, One Truneb Reusable Neb Kit	each
FROG-TRUWB-PB	Frog Pediatric Nebulizer System with one Disposable Neb Kit, one Reusable Neb Kit and Black Nylon Carry Bag with Piggyback Serial Labels	each
<b>Accessories</b>		
50002	Frog Nylon Carry Bag	each
50016	Truneb Reusable Neb Kit with Angle Mouthpiece and Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 In. Diameter X 1/2 In. Thick	6/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

\*\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

# Respiratory - Aerosol Therapy

## DOG PEDIATRIC NEBULIZER SYSTEM



Designed to reduce anxiety during treatments and encourage compliance, the child-friendly **Roscoe Dog Compressor Nebulizer** provides effective aerosol treatments with a state-of-the-art piston compressor. Delivers 43 PSI at maximum pressure and 16 PSI operating pressure (liter flow range 5 - 8 LPM). Equipped with a convenient nebulizer parking post. Comes complete with compressor, disposable nebulizer, child mask, angled mouthpiece, 7-foot kink-resistant tubing, five extra filters and a five year warranty. Also available with a reusable neb kit option (Model No. DOG-TRU).



BAG-NBDG  
Nylon Carry Bag

### FEATURES

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- Nylon carry bag available BAG-NBDG
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone
- DOG-TRU includes TruNeb Reusable Neb Kit
- HCPCS\* Code: E0570\* for nebulizer with compressor  
A7003\* for disposable nebulizer kits  
A7005\* for reusable nebulizer kits  
A7013\* for foam filters  
A7015\* for aerosol mask

### SPECIFICATIONS

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	Less than 0.5 µm
Unit Dimensions:	8 4/5" (L) x 5 2/5" (W) x 9" (H)
Product Weight:	3.4 lbs.
Maximum Pressure:	43 PSI
Operating Pressure:	16 PSI
Liter Flow:	5 - 8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

Item #	Description	UOM
NEB-DOG	Dog Pediatric Nebulizer System with one Disposable Neb Kit	each
DOG-TRU	Dog Pediatric Nebulizer System with one Disposable Neb Kit, one TruNeb Reusable Neb Kit	each
DOG-TRUWB-PB	Dog Pediatric Nebulizer System with one Disposable Neb Kit, one Reusable Neb Kit and Black Nylon Carry Bag with Piggyback Serial Labels	each
<b>Accessories</b>		
BAG-NBDG	Dog Nylon Carry Bag	each
50016	TruNeb Reusable Neb Kit with Angle Mouthpiece and Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 In. Diameter X 1/2 In. Thick	6/pk

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

\*\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

# Respiratory – Aerosol Therapy



## SOCCKER BALL NEBULIZER COMPRESSOR SYSTEM

Designed to reduce anxiety during treatments and encourage compliance, the child-friendly Roscoe Soccer Compressor Nebulizer provides effective aerosol treatments with a state-of-the-art piston compressor. Delivers 29 PSI at maximum pressure and 8 to 13 PSI operating pressure (liter flow range 5 - 8 LPM). Equipped with a convenient nebulizer parking post. Comes complete with compressor, disposable nebulizer, mask, mouthpiece, 7-foot kink-resistant tubing, five extra filters and a five year warranty. Also available with a reusable neb kit option (Model No. SOCCER-TRU & SOCCER-TRUWB-PB).



### FEATURES:

- Child-friendly design helps to reduce anxiety during treatments
- Provides smooth, quiet treatments
- Equipped with a convenient nebulizer parking post
- Complete with compressor, nebulizer, child mask, angled mouthpiece, tubing, five extra filters and instruction manual
- TruNeb Reusable Neb Kit available for purchase with the compressor or alone
- SOCCER-TRU & SOCCER-TRUWB-PB includes Reusable Neb Kit

### OPTIONAL ACCESSORIES:



Item #	Description	UOM
NEB-SOCCER	Soccer Ball Nebulizer System with Disposable Neb Kit	each
SOCCER-TRU	Soccer Ball Nebulizer System with Reusable TruNeb Kit	each
SOCCER-TRUWB-PB	Soccer Ball Nebulizer System with Disposable Neb Kit, Reusable TruNeb Kit, Carry Bag, and Piggyback Serial Labels	each
Accessories		
50016	TruNeb Reusable Neb Kit With Angle Mouthpiece And Tubing	10/cs
CIF-1128PK6	Foam Filter, 5/8 in. Diameter x 1/2 in. Thick	6/pk

### PRODUCT SPECIFICATIONS:

Method of Nebulization:	Piston Compressor
Medication Capacity:	5 ml
Particle Size:	0.5 µm to 10 µm
Unit Dimensions:	5.25" L x 6" W x 7.5" H)
Product Weight:	3.3 lbs.
Maximum Pressure:	29 PSI
Operating Pressure:	8 – 13 PSI
Liter Flow:	8 LPM
Electrical Requirement:	120 VAC/60Hz
Warranty:	Limited 5 Year

### HPCPS CODES:

E0570*	Nebulizer with Compressor
A7003*	Disposable Nebulizer Kits
A7005*	Reusable Nebulizer Kits
A7013*	Foam Filters
A7015*	Aerosol Mask

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

\*\* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

## Nebulizer Kits

Did you know? By setting up patient resupply schedules for nebulizer kits, you could increase your annual reimbursement revenue per patient by 30%. Additionally, by switching disposable to reusable neb kits, you can save on costs and increase your profit per patient.



## Roscoe Medical® Nebulizer Kits

Our Roscoe Medical® Nebulizer Kits are available in either disposable or reusable models to meet the needs of a wide range of patients. They deliver fast and effective treatment efficiently while the easy-seal, threaded cap and anti-spill jar prevent leakage of costly medication. Provides 80% of particles in respirable range.



#NEB-ADLTM



#NEB-ROSKT



#50016

Item #	Description	Qty	HCPCS Codes
NEB-ROSKT	Roscoe® Disposable Nebulizer Kit with angled mouthpiece, 7 ft. tube, "T" adapter, 6 in. reservoir tube and disposable nebulizer	50/cs	A7003: Disposable Nebulizer Kits A7005: Reusable Nebulizer Kits
50016	Roscoe® TruNeb Reusable Neb Kit with nebulizer, angled mouthpiece, tubing and 10 cc dishwasher-safe medication cup included	10/cs	
NEB-PEDM	Roscoe® Nebulizer Kit with pediatric mask, neb cup, mask and tubing included	50/cs	
NEB-ADLTM	Roscoe® Nebulizer Kit with adult mask, neb cup, mask and tubing included	50/cs	

## Salter® Nebulizer Kits

Salter's® T-piece Jet Nebulizers deliver fast and effective treatment and feature a one-piece jet design that delivers high output rates, resulting in faster treatment times and consistent medication delivery. A full 10 cc capacity provides medication dosing flexibility. The anti-spill design helps prevent medication waste and allows the unit to be used at an angle to accommodate various patient positions.



#8900

Item #	Description	Qty
8900	Salter® Disposable Nebulizer with 6" reservoir tube and 7 ft. tubing	50/cs
8924	Salter® Disposable Nebulizer with adult mask and 7 ft. tubing	50/cs
8906	Salter® Disposable Nebulizer with pediatric mask, and 7 ft. tubing	50/cs

### HCPCS CODES

A7003: Disposable Nebulizer Kits

## Westmed® VixOne® Nebulizer Kits

The VixOne® Small Volume Jet Nebulizer is one of the enduring standards in respiratory therapy. A full 10 cc capacity provides medication dosing flexibility. The anti-spill design helps prevent medication waste and allows the unit to be used up to a 45 degree angle to accommodate various patient positions.



#0382

Item #	Description	Qty
0382	PediNeb Nebulizer Kit with pacifier, 45° elbow and 7 ft. tubing	10/cs

### HCPCS CODES

A7003: Disposable Nebulizer Kits

## BreatheRite™ Rigid Valved Holding Chamber and Masks

The BreatheRite™ Valved Holding Chamber (VHC) combines design advantages with portability to provide both effective medication delivery and encourages increased patient use.



M10/M20

Item #	Description	Qty
BR300	BreatheRite™ Valved Holding Chamber	20/cs
M10	Aerosol Chamber Mask – Adult	10/cs
M20	Aerosol Chamber Mask – Pediatric	10/cs



BR300

## Aerosol Masks



Item #	Description	Qty
<b>Adult Masks</b>		
ADLT-MASK	Roscoe® Adult Aerosol Mask without tubing	50/cs
8100	Salter® Adult Elongated Aerosol Mask without tubing, elastic strap	50/cs
<b>Pediatric Masks</b>		
PED-MASK	Roscoe Pediatric Aerosol Mask with tubing	50/cs
1120	Salter® Pediatric Aerosol Mask without tubing	50/cs

**HCPCS CODES**

A7015: Aerosol Mask

## Large Volume Nebulizers

The Large Volume Non-Prefilled Nebulizer is a disposable nebulizer designed for use with clinician-provided water that delivers humidified gas in concentrations ranging from 28% to 98%.



CF002002

Item #	Description	Qty
CF002002	Large Volume Nebulizer, empty, 350mL, 35% - 100%	25/cs

**HCPCS CODES**

A7007: Large Volume Nebulizer, Unfilled

## Disposable Corrugated Tubing

Corrugated Tubing for aerosol and IPPB applications offers flexibility and convenience to meet all therapeutic needs. For heated and non-heated use, tubing comes complete with easy-fit tubing ends. Segmented every six inches with cuffs on all tube ends that lock to 7/8 inches (22 mm) connections.



001405

Item #	Description	Qty
001405	Corrugated Aerosol Tubing w/ 22 mm Cuffs, 100' length segmented every 6"	1/ea

**HCPCS CODES**

A7010: Disposable Corrugated Tubing

## PEP Therapy/Vibratory Clearance Devices

- Home care
- Loosens mucus
- Uses vibration to open airways
- Small, reusable, non-invasive, drug-free



21-1530-EA



21-1015-EA

Item #	Description	Qty
21-1015-EA	Acapella DM Low Flow (Blue) Vibratory PEP Therapy System with Mouthpiece	10/cs
21-1530-EA	Acapella DM High Flow (Green) Vibratory PEP Therapy System with Mouthpiece	10/cs

### HCPCS CODES

E0484: Non-elec oscillating pep dvc.

## Trach Mist Aerosol Drainage Bag with Y



Item #	Description	Qty
001561	Trach Mist Aerosol Drainage Bag with Y	50/cs

## Trach Mask



Item #	Description	Qty
001225	Adult Trach Mask	50/cs

## Oxygen Enrichment Connector






Item #	Description	Qty
AG1642	Oxygen Enrichment Connector	1/ea

## Incentive Spirometers





Item #	Description	Qty
BAX001902A	Incentive Spirometers	1/ea

## Replacement Filters for Nebulizer Compressors

Item #	Unit	Qty	Dimensions	OEM	AG	HCPCS	
PB8NEB5FLT		Roscoe® RITE-NEB5 Nebulizer Compressor Filter	25/pk	1/2" x 1/8" Foam Filter	PB8NEB5FLT	N/A	A7013
CIF-1128PK6		Works with Roscoe Pediatric Nebulizer Systems	6/pk	5/8" x 1/2" Foam Filter	CIF-1128PK6	N/A	A7013
90325		For NEB50-PSI	each	1 13/16" x 3/4" Felt Filter	90325	N/A	A7014

## Miscellaneous Nebulizer Filters

CIF-1128PK6		Drive / MedQuip Nebulizers	6/pk	5/8" x 1/2" Foam Filter	18090	AF1036898	A7013
CIF-1128PK6		Salter AIRE Plus	6/pk	5/8" x 1/2" Foam Filter	8501	F153	A7013



## High Vacuum / Low Flow Aspirator



The **Roscoe Aspirator** is the perfect solution for all your suction needs. It is simple to use and has a durable design to withstand the rigors of homecare use. It produces a maximum vacuum of 675 mmHg and features a high-quality vacuum regulator with anti-vibration vacuum gauge setting, an anti-bacteria filter design helps prevent contamination, oil-less, permanent lubrication, and a collection bottle with overflow protection. Designed for easy portability at just 7.7 pounds, its easy to carry and easy to clean.

### FEATURES:

- Filter design helps prevent contamination
- High quality vacuum regulator and anti-vibration vacuum gauge setting
- Oil-free, permanent lubrication
- Collection bottle with overflow protection
- Superior performance and high reliability for continuous operation
- Protective gauge guard prevents damage
- Easy carry and cleaning
- HCPCS Code: E0600\*
- Respiratory suction pump, home model, portable or stationary, electric

### PRODUCT SPECIFICATIONS:

Dimensions:	13" L x 6" W x 8.25" H
Weight:	7.7 lbs
Maximum Vacuum:	-675 +/- 3% mm Hg
Max Air Flow:	-15-17 LPM
Power Cord Length:	6'
Electrical Requirements:	110-120V/220-240V (50-60Hz)
Pressure Display:	Analog Gauge
Filter Pore Size:	0.45 microns
Collection Capacity:	800 cc
Operating Temperature:	32° - 104°F
Operating Humidity:	95% RH
Warranty:	Two Year Limited

### OPTIONAL ACCESSORIES:



Item #	Description	Qty
ASP-ROS	Roscoe High Vacuum / Low Flow Aspirator	each
90386	Suction Kit (Includes elbow, filter, collection jar, 6' tubing and two 7" tubes)	each
<b>Replacement Parts:</b>		
1160B-1	800 cc Disposable Collection Jar	each
ILB-ROS	Filter, 3/8" to 1/2" barbed	each
ROS-610E	Elbow	20/cs

\* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



## truAIRE-PASP Portable Aspirator



The truAIRE-PASP Portable Aspirator is engineered for rapid response and reliable performance, the truAIRE-PASP delivers powerful suction in a compact, lightweight design. Ideal for emergency medical services and home healthcare, this unit offers exceptional portability without compromising suction strength. With a simple flow control adjustment knob and long-lasting rechargeable battery, this device is ready to go when it matters the most.



Easy to read vacuum gauge and variable suction settings



True portability - Unit can operate while inside carry bag



Unit can be charged while in use

### PRODUCT FEATURES:

- Compact and portable design, perfect for oral, nasal, and tracheal suctioning
- Can be operated while in carry bag for on-the-go use
- Variable vacuum regulation from 50 to 550 mmHg
- Complete with carry bag, dual elbow tubing assembly, 800cc collection jar, bacteria filter, 6' suction tubing, and user manual
- Replaceable, lithium battery has a 90-minute duration with a working cycle of 30 minutes on and 30 minutes off
- Comes complete with A/C D/C power supply with 90° connector
- 2 Year warranty on the device, 1 Year warranty on the battery
- HCPCS Code: E0600\*

### TECHNICAL SPECIFICATIONS:

Dimensions:	8.25" L x 7.5" W x 7.75" H
Weight:	3.6 lbs
Maximum Vacuum:	600mmHg (Recommend Range 50 – 550mmHg)
Flow Range:	18 LPM +/- 10%
Sound Level:	≤ 60 dBA
Electrical Requirements:	100-240V AC, 50-60Hz
Operating Temperature:	32° - 104°F (0°-70°C)
Operating Humidity:	0 - 95%

Item #	Description	UPC	Qty
ASP-ROSP	truAIRE-PASP Portable Aspirator	092237626672	2/CA
ASP-KIT	Disposable Suction Kit for ASP-ROSP (Patient Kit contains Canister, 2-Elbow Tubes, 6' Tube & Bacteria Filter)		1/ea
1160B-1	Bemis 800cc Suction Canister & LID		70
55N56A	6' Suction Tubing		50
X1PASP01	Elbow Tubing Assembly for ASP-ROSP		1/ea
X1PASP03	Bacteria Filter for ASP-ROSP		1/ea
PORT-CARADP	DC Car Charger		1/ea
X1PASP04	Replacement Battery - ASP-ROSP		1/ea

## Cliq 50 PSI Compressor



The CX-500 is a 50 PSI compressor designed for efficient high pressure performance for extended periods of time. It is intended to administer humidified air or medication into the airways to treat respiratory disorders.

Designed for use with heated nebulizers, compatible hand-held nebulizer systems, dual action humidifiers/nebulizers, mist tents, high output pneumatic nebulizers and deluxe nebulizers with oxygen diluter. Includes U-shaped DISS connection for single neb bottle and 2 extra air filters.



### FEATURES

- Powerful, efficient, dependable and quiet
- Durable and easy to clean with a portable design for easy transport
- Designed for extended periods of time aerosol and humidification applications
- For use with heated nebulizers, compatible hand-held nebulizer systems, dual action humidifiers/nebulizers, mist tents, high output pneumatic nebulizers and deluxe nebulizers with oxygen diluter
- Excellent airflow characteristics mean lower operating temperatures and longer life
- HCPCS Code: E0565\* Compressor, Air Poer Source

### SPECIFICATIONS

Power:	110 – 120 V/50 – 60 Hz
Operating Temperature:	41° - 104°F
Operating Pressure:	0 – 95 PSI
Operating Flow Rate:	0 psi: 37 LPM/50 psi: 17 LPM
Unit Dimensions:	14" (L) x 7.5" (W) x 12.5" (H)
Weight:	14.9 lbs.
Warranty:	1 Year

Item #	Description	Qty
CX-500	Cliq 50 PSI Compressor	each

## EasyAir Compressor

### WITH FLOW OR PRESSURE GAUGE


EasyAir Compressor is a high performance portable medical air compressor designed to reliably supply compressed air 24 hours a day.

EasyAir is ideal for continuous use in any clinical or homecare environment. Its sound dampened case is easy to clean, allows quiet operation and is durable enough for years of heavy use. The oil-less compressor runs cooler and provides clean, breathable air. A precise dial and lock pressure regulator ensures air flow and pressure remain constant.

The unit is available with either a pressure gauge or flow gauge.

HCPCS: E0565

Compressor, air power source for equipment which is not self-contained or cylinder driven.



Item #	Description	Qty
PM15F	Easy Air Compressor with Flow Gauge	each
PM15P	Easy Air Compressor with Pressure Gauge	each

\* This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

## Precision EasyVac Heavy Duty Aspirator

EasyVac Aspirator provides reliable and safe suction of fluids, using up to 635 mmHg of vacuum at continuous operation. EasyVac provides additional safety that guards against the suction canister from overflowing. The lightweight device is easy to transport and store when necessary and won't overheat due to a resetting thermal limit switch.



Item #	Description	Qty
PM60	Precision EasyVac Heavy Duty Aspirator	each

#### HCPCS CODES

E0600: Respiratory Suction Pump, Home Model, Portable Or Stationary, Electric

## Precision EasyGo Portable Suction Machine

The EasyGoVac battery-powered aspirator combines ease of use with powerful suction to create a portable unit that is ideal for hospitals, nursing homes, ambulances and homecare. EasyGoVac produces vacuum up to 483 mmHg and provides dependable continuous suction power to remove respiratory fluids. The high-performance battery delivers 2 hours of run time, more than any aspirator of comparable size.



Item #	Description	Qty
PM65	Precision EasyGo Portable Suction Machine	each

#### HCPCS CODES

E0600: Respiratory Suction Pump, Home Model, Portable Or Stationary, Electric

Replacement Parts		
Item #	Description	Qty
503106	Precision RPL Power Cord for PM65	each
503974	Precision AC Battery Charger	each

## Bemis Suction Canister with float top, 800cc

These Hi-Flow Suction Canisters are a high-quality way to meet fluid collection and disposal needs. They provide a high flow rate, so suctioning is quick and efficient during procedures.



1160B-1

Item #	Description	Qty
1160B-1	Bemis Suction Canister with float top, 800cc	1/ea or 70/case

## Suction Kits

Convenient suction kits are available that take the guess work out of suction supplies. See popular kits available below or call us for more information on building your own custom kit!

*\* To keep your heavy-duty aspirator performing optimally, replace its parts regularly and easily with the Roscoe Suction Kit.*



Item #	Description	Qty
90386	Suction Kit for ASP-ROS & ASP-HD; Kit includes: 50038, ILB-ROS, ROS-610E, 7" Silicone Tubing (2), 6' Suction Tubing	each
<b>Kits</b>		
KIT-7305D	Suction Kit for 7305d & 7305p; Kit includes: 50038, ILB-ROS, 7305D-023, 4" Silicone Tubing, 6' Suction Tubing	each
KIT-50004	Suction Kit for 50004 & 50006; Kit includes: 50038, ILB-38B, ROS-610E, 6" Silicone Tubing (2), 6' Suction Tubing	each
KIT-018D	Suction Set Up Kit- Devilbiss; Kit includes: 50038, ILB-38B, ROS-610E, 16" Vinyl Tubing, 6' Suction Tubing, Bulb-Tip Yankauer w/ Control Vent	each
KIT-019	Suction Setup Kit; Kit includes: 50038, ILB-38, ROS-610E, 14" Vinyl Tubing, 6' Suction Tubing, Bulb-Tip Yankauer w/ Control Vent	each
KIT-020	Suction Setup Kit; Kit includes: 50038 (2), 14" Vinyl Tubing, 6' Suction Tubing (2), 20" Suction Tubing, Bulb-Tip Yankauer w/ Control Vent	each
KIT-021	Suction Kit for PM60; Kit includes: 14" Vinyl Tubing, ROS-610E	each
KIT-022	Suction Kit for Devilbiss; Kit includes: 12" Silicone Tubing, ROS-610E	each

## Suction Filters

The replacement Bacteria Filter prevents fluids and other contaminants from entering and damaging your suction pump.



ILB-ROS



ILB-38



ILB-38B

Item #	Description	Qty
ILB-ROS	Suction Bacteria Filter with 3/8" Stepped Barb, for 7305D	each
ILB-38	Suction Bacteria Filter with 1/4" Stepped Barb x 1/8" MPT, for Roscoe Suction Kits 50004, 50006 and other suction machines	each
ILB-38B	Suction Bacteria Filter, 1/4"-3/8" Stepped Barb x Barb, for Roscoe Suction Kits 50004 and 50006, SUC-COMP, and 601 and 605 DeVilbiss 7305D, and also other suction machines	each

## Accessories

### Elbow for DeVilbiss 7305D Suction Canister

The Replacement Elbow allows users to easily and securely connect tubing to their suction canister and prevents tubing from kinking. It is designed for use with the DeVilbiss 7305D Suction Canister only.



Item #	Description	Qty
7305D-023	Elbow for DeVilbiss 7305D Suction Canister	each

### Elbow for ASP-ROS

The Replacement Elbow allows users to easily and securely connect tubing to their suction canister and prevents tubing from kinking. It is designed for use with the 601, 605, and ASP-ROS suction machines only.



Item #	Description	Qty
ROS-610E	Elbow for ASP-ROS	each

## Silicone Tubing

High quality, medical grade silicone tubing cut to the perfect length for use between the aspirator, suction filter and collection canister.

Item #	Description	Qty
TK-MA	Silicone Tubing for ASP-ROS, 7"	10/pk
TK-SUC4	Silicone Tubing, 4"	12/pk



## Suction Tubing

This Suction Tubing is durable and safe, with a "no memory" design that lies flat for exceptional performance. It features scallop female connectors on both ends to fit a wide range of equipment and is latex-free.

Item #	Description	Qty
55N52A	Suction Tubing, 3/16" x 18"	50/ca
55N56A	Suction Tubing, 3/16" x 6'	50/ca



## Bulb-Tip Yankauer

The Bulb-Tip Yankauer allows for precise suctioning without damaging the surrounding tissue. It features a slip-resistant handle and smooth inner surface for high-quality performance.

Item #	Description	Qty
55K86	Bulb-Tip Yankauer without Control Vent	50/ca
55K86V	Bulb-Tip Yankauer with Control Vent	50/ca



55N52A

## Trach Care Kit

The Trach Care Kit features everything needed for convenient and efficient infection control during a tracheostomy procedure. It includes two vinyl gloves, a moisture-proof drape, four 4" x 4" 4-ply gauze squares, two cotton tip applicators, a pliable cleaning brush, and 4" x 4" 6-ply gauze trach dressing in a 3-compartment tray with a peel-down Tyvek lid.



Item #	Description	Qty
DYND40582	Trach Care Kit	20/ca

## Equipment Bags and Covers

Meets accreditation requirements! Our equipment bags/covers are ideal for home and long-term care environments.



#BAG-1015

#BAG-2430R

Item #	Description	Qty
BAG-1015	Clear Plastic Nebulizer Bag, 1 ml (10 in. x 15 in.)	1000/cs
BAG-15924	Clear Plastic Nebulizer Bag, 1 ml (15 in. x 9 in. x 24 in.)	500/roll
BAG-2430R	Equipment Cover, 1.5 ml (24 in. x 30 in.)	500/roll



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

[www.compasshealthbrands.com](http://www.compasshealthbrands.com)

rev\_03\_251013