

# 2026 Bath Safety Catalog

Enhancing quality  
of life is our  
driving force

**drive**<sup>®</sup>   
M E D I C A L



# Bathroom Safety

## Bathroom Safety Solutions

We at Drive Medical know that the bathroom is one of the most hazardous rooms in the house. To that end, building a safer bathroom is an important part of maintaining independence and a higher quality of life. Our bath safety products were clinically designed to meet the needs of individuals facing mobility challenges in both the long and short term. From shower chairs, to bath lifts, we provide a full spectrum of bath safety products to provide customers with customized solutions that help enhance their safety and comfort in daily life. With a focus on durability, portability, ease of assembly, and easy to clean surfaces, we've thoughtfully designed our products to suit the busy lifestyles of customers and caregivers.



# Bathroom Safety

## M+D X-Fold Shower Chair



Item #	Color
MD50100WH	Pearl
MD50100TP	Stone
MD50100GR	Storm

- Innovative X-Fold™ design with 3 height settings, no buttons or tools
- Extra-roomy seat (~32% larger than standard) for comfort & stability
- Removable armrests for easier transfers & personalized support
- Angled anti-slip feet with water-dispersing treads for wet-surface traction
- Advanced drainage system reduces pooling & speeds drying
- Modern color options to match contemporary bathrooms

Coming Soon



## M+D Folding Shower Chair



Item # MD770525WH

- Easily set down and access your shower spray with built in cut-outs in the seat
- Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15" to 22"
- Push button, snap together design requires no tools and assembles in just minutes

Coming Soon



# Bathroom Safety

## Deluxe Aluminum Bath Chair

Item # RTL12202KDR

- Easy, tool-free assembly of back seat and legs
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Made using recycled materials

[Shop Now](#)



## Deluxe Aluminum Shower Chair

Item # 12202KD-4

- Easy, tool-free assembly of back, seat and legs
- Aluminum frame is lightweight, durable and corrosion proof
- New pinch-free cover allows for push pins to be depressed without pinching finger
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Support collar prevents leg movement and reduces rattle
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Shower Chair with Back and Removable Padded Arms

Item # 12445KD-1

- Removable padded armrests provide comfort and added security
- Blow-molded bench and back provide comfort and strength
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Leg height adjusts in 1" increments
- Made using recycled materials



[Shop Now](#)

## PreserveTech™ 360° Swivel Bath Chair

Item # RTL12A001-GR

- Compact bath chair featuring 360 degree swivel action and locking at every 90 degrees
- Each leg includes a non-skid, non-slip tip
- Padded armrests provide superior support
- Pivoting armrests can be raised out of the way
- Seat height adjustments on this range from 15" to 22"
- Drainage holes help to prevent water from pooling
- Push button, snap together design requires no tools and assembles in just minutes



[Shop Now](#)

# Bathroom Safety

## PreserveTech™ Aquachair Bathing System with Bidet

Item # RTL12A004-WH

- Unique shower chair with bidet and shower spray enables full-body cleansing without standing
- Padded non-slip seat and backrest offer support and comfort
- Adjust the height with the push of a button to achieve the ideal seated position
- Legs include a non-skid, non-slip tip for additional security on wet surfaces
- Lightweight and durable-rust-resistant aluminum frame supports up to 350 lbs
- Anti-microbial treated to help prevent stains and odors from buildup of bacteria, mold and mildew
- No electrical connection or batteries necessary

[Shop Now](#)



## PreserveTech™ Aquachair Transfer Accessory with Curtain Guard Protection

Item # RTL12A004-TR

- Easily installs to the PreserveTech™ Aquachair Bathing System with Bidet in minutes
- Bench design allows user to slide over side of tub from a safer seated position
- Patented "U" shape cut-out design keeps shower curtain liner in the tub and prevents water spillage
- Adjust leg height with the push of a button
- Nonslip feet offer added security on bathroom floors

[Shop Now](#)



# Bathroom Safety

## Shower Chair with Folding Back

Item # RTL12606

- Adds safety and support in the bath tub and shower
- Easy, tool-free assembly and height adjustable legs
- Folding back and removable legs make this shower chair great for storage and transport
- Lightweight design with stable legs and non-slip feet
- Armrests and soap tray add comfort and convenience
- Made using recycled materials

[Shop Now](#)



## Deluxe Bariatric Shower Chair with Cross-Frame Brace

Item # 12021KD-1

- Frame has width/depth cross braces attached with aircraft-type rivets
- Durable, blow-molded plastic bench
- Aluminum frame is lightweight, durable and corrosion proof Drainage holes reduce slipping
- Euro-style clip for leg height adjustment Height adjusts in 1" increments
- Made using recycled materials

[Shop Now](#)



## Deluxe Bariatric Shower Chair with Cross-Frame Brace

Item # 12029-2

- Frame has width/depth cross braces attached with aircraft-type rivets
- Large seat with 16.5" depth
- Seat height adjusts in 1" increments
- Durable, blow-molded plastic bench
- Tool-free removable back
- Aluminum frame is lightweight, durable and corrosion proof
- Drainage holes reduce slipping
- Made using recycled materials

[Shop Now](#)



## Deluxe Aluminum Shower Bench without Back

Item # RTL12203KDR

- Tool-free assembly of seat and legs
- Drainage holes in bench reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Deluxe Aluminum Shower Bench without Back



Item # 12203KD-4

- Tool-free assembly of seat and legs
- Drainage holes in bench reduce slipping
- Frame is lightweight, durable and corrosion proof
- New pinch-free cover allows for push pins to be depressed without pinching finger
- Legs are height adjustable in 1" increments, and angled with suction style tips to provide additional stability
- Support collar prevents leg movement and reduces rattle
- Made using recycled materials

[Shop Now](#)



## Swivel Seat Shower Stool



Item # RTL12061M

- The rotating padded seat reduces twisting and reaching while in the shower
- Seat rotates 360° and locks every 90°
- The rotating top of the seat can be removed for easy cleaning
- Easy tool-free assembly
- Ideal for individuals who have balance or mobility issues
- Seat height can be adjusted from 16" to 22", in 1.18" increments
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Folding Shower Chair

Item # 12486

- Product available with back and without
- Conveniently folds flat
- Blow-molded bench provides comfort and strength
- Drainage holes in seat reduce slipping
- Aluminum frame is lightweight, durable and corrosion proof
- Angled legs provide additional stability
- Made using recycled materials



[Shop Now](#)

## Shower Stool

Item # RTL12004KD

- Aluminum frame is lightweight, durable and corrosion proof
- Impact-resistant, composite seat
- Legs are height adjustable in 1" increments, and are crack proof and tarnish resistant
- Easy, tool-free assembly
- Made using recycled materials



[Shop Now](#)

# Bathroom Safety

## M+D Shower Stool

Item # MD770535WH

- Easily set down and access your shower spray with built in cut-outs in the seat.
- Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 15" to 22".
- Push button, snap together design requires no tools and assembles in just minutes.



Coming Soon

## M+D Transfer Bench

Item #	Color
MD50200WH	Pearl
MD50200GR	Storm

- 360° rotating seat with 90° lock points and side-to-side sliding transfer allows users to sit securely outside the tub and glide into position.
- 21 height adjustments accommodate tubs with walls from 14.5" to 20.5", creating a customized, level fit across a wide range of bathroom setups.
- Angled anti-slip feet with water-dispersing treads help channel water away and maintain stronger traction on shower and bathroom floors.
- Durable frame supports up to 350 lb and features a dual locking system with four secure seat positions and audible click engagement.



Coming Soon

# Bathroom Safety

## Transfer Tub Bench

Item # 12011KD-2



- Tool-free assembly of back, legs and arm
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

[Shop Now](#)



## Three Piece Transfer Tub Bench

Item # RTL12031KDR



- Tool-free assembly of back, legs, and arm
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Splash Defense™ Transfer Bench with Curtain Guard Protection

Item # RTL12032KDR

- New "U" shape cut-out design enables the shower curtain liner to stay within the bath tub and avoiding water spillage outside onto the bathroom floor
- Reversible to accommodate any bathroom
- Tool-free assembly of back, legs and arm
- Height adjusts in ½" increments
- Durable blow-molded plastic bench and backrest
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

[Shop Now](#)



## Folding Universal Sliding Transfer Bench

Item # RTL12075

- Transfer Bench slides left and right along frame for easy entry into and out of tub
- Seat height adjusts in ½" increments
- Folds flat for transport and storage
- Comes standard with removable soap dish
- Comes with safety net to catch dropped personal items and prevent bending
- Comes with seat belt for added safety
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Padded Transfer Bench

Item # 12005KD-1

- Comfortable cushioned seat and backrest
- Tool-free assembly of back, legs and arm
- New pinch-free cover allows for push pins to be depressed without pinching fingers
- A-frame construction provides stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs



[Shop Now](#)

## Combination Padded Transfer Bench/Commode

Item # 12005KDC-1

- Combines transfer bench and commode into one
- Tool-free assembly of back, legs and arm
- Comfortable cushioned seat and backrest
- Aluminum A-frame construction is lightweight and provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- No exposed hardware to injure patient
- Extra-large suction cups provide added safety
- New pinch-free cover allows for push pins to be depressed without pinching finger



[Shop Now](#)

## Bariatric Transfer Bench

Item # 12025KD-1



- Tool-free assembly of back, legs and arm
- A-frame provides additional stability
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- All transfer benches are reversible to accommodate any bathroom
- Euro-style clip ensures a secure fit
- No exposed hardware to injure patient
- Large, locking suction cups for added safety
- Weight Capacity 500 lbs
- Made using recycled materials



[Shop Now](#)

## 2-in-1 Locking Raised Toilet Seat with Tool-free Removable Arms

Item # RTL12027RA

Retail

Item # 12027RA-4BULK



- Adds 5" to height of toilet seat
- Tool-free arms can be removed or added
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Easily attaches to toilet bowl, tool-free
- Easy-to-clean and fits most toilets
- Made using recycled materials



[Shop Now](#)

# Bathroom Safety

## Premium Plastic, Raised, Elongated Toilet Seat with Lock

Item # 12013

- Adds 5" to height of toilet seat
- Tool-free arms are 21" apart and can be removed or added as needed
- Wide opening in front and back
- Heavy-duty molded plastic construction provides additional strength and durability
- Standard locking mechanism ensures safety
- Lightweight and portable
- No tools required for installation, fits most elongated toilets
- Made using recycled materials

[Shop Now](#)



## PreserveTech™ Universal Raised Toilet Seat

Item # RTL12C002-WH

- Accommodates three different seat heights (2", 4" and 6") by maneuvering just two components
- Pivoting armrests can be raised out of the way or removed depending on your needs or bathroom environment
- Fits comfortably on most standard and elongated toilets
- Includes a lid with recessed grip for better hygiene and aesthetic and effortless opening
- Contoured to provide a safe and comfortable elevation

[Shop Now](#)



# Bathroom Safety

## PreserveTech™ Secure Lock Raised Toilet Seat



Item # RTL12C003-WH

- Three locking points of security located on the front, left and right side providing the highest level of stability and security
- This toilet seat adds an extra 5" of height
- Padded armrests provide support and create a comfortable, non-slip surface for a secure grip
- Easily adjust armrest height to your needs
- Fits comfortably on most standard and elongated toilets
- No Tools Required! Simply turn the adjustable knobs to tighten the locking plate
- Made using recycled materials



[Shop Now](#)

## Raised Toilet Seat with/without Lid



Item # 12065

With Lid

Item # RTL12064

Without Lid

- Designed for individuals who have difficulty sitting down or standing up from the toilet
- Heavy-duty molded plastic construction
- Locking device with larger, heavy-duty "worm screw" and locking plate for a safe, secure fit
- Lightweight and portable
- Fits most toilets
- No tools required for installation
- Easy-to-clean
- Made using recycled materials



[Shop Now](#)

# Bathroom Safety

## PreserveTech™ Raised Toilet Seat with Bidet With Warm Water

Item # RTL12C006-WH

- Control panel activates and adjusts water pressure and cleansing modes
- Bidet provides front or female and back cleansing depending on needs
- Retractable bidet nozzle and guard for a more sanitary experience
- Contoured for safe and comfortable elevation
- Antimicrobial treated raised toilet seat helps minimize the growth of bacterial odors, mold, and fungi on treated surfaces

[Shop Now](#)



## PreserveTech™ Raised Toilet Seat with Bidet With Ambient Water

Item # RTL12C005-WH

- Control panel activates and adjusts water pressure and cleansing modes
- Bidet provides front or female and back cleansing depending on needs
- Retractable bidet nozzle and guard for a more sanitary experience
- Contoured for safe and comfortable elevation
- Antimicrobial treated raised toilet seat helps minimize the growth of bacterial odors, mold, and fungi on treated surfaces

[Shop Now](#)



# Bathroom Safety

## Premium Raised Toilet Seat with Removable Arms



Item # 12403

- Allows individual to use existing toilet seat and lid
- Tool-free arms are 19.75" apart and can be removed or added as needed
- Arm sleeves and receivers are made of metal
- Heavy-duty molded plastic construction
- Lightweight and portable
- Easily attaches to toilet bowl
- No tools required for installation
- Easy to clean
- Accommodates most elongated toilets
- Made using recycled materials

[Shop Now](#)



## Hinged Toilet Seat Riser

Item # RTL12607

Regular

Item # RTL12608

Elongated



- Raises the height of the toilet 3.5" for individuals that require assistance sitting on and standing from the toilet
- Hinge design allows riser to be moved out of the way when not in use
- Heavy-duty molded plastic construction provides strength and durability
- Lightweight, portable, and easily attaches to toilet bowl
- Two models to accommodate either standard or elongated toilet bowls
- Made using recycled materials

[Shop Now](#)



# Bathroom Safety

## Toilet Safety Frame with Padded Arms

Item # RTL12000

- Anodized aluminum is sturdy and lightweight
- Padded armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions)
- Powder-coated aluminum bracket easily attaches frame to toilet
- Bracket secures safety frame in place

[Shop Now](#)



## Toilet Safety Frame

Item # 12001-4

- Anodized aluminum is sturdy and lightweight
- Plastic armrests provide additional comfort and support
- Comes with height-adjustable (1" increments) and width-adjustable arms (2 positions)
- Powder-coated aluminum bracket easily attaches frame to toilet
- Bracket secures safety frame in place

[Shop Now](#)



# Bathroom Safety

## Free-standing Toilet Safety Rail

Item # RTL12079

- Free standing around toilet
- Easy, tool-free assembly using wing nuts
- Fits standard and elongated toilets
- Comes with padded armrests for added comfort
- Comes with a magazine rack for added convenience



[Shop Now](#)

## Chrome Knurled Grab Bar

Item # 12112-3	12 inch
Item # 12116-3	16 inch
Item # RTL12118	18 inch
Item # 12124-3	24 inch
Item # 12132-3	32 inch



- Knurled grab bar diameter provides a no-slip grip
- Mounts vertically or horizontally
- 1 ½" distance from wall



[Shop Now](#)

# Bathroom Safety

## White Powder-Coated Grab Bar

Item # RTL12018      18 inch  
Item # 12024-3      24 inch

- Durable, white powder-coated steel is attractive and easy to clean
- Mounts vertically or horizontally
- 1 ½" distance from wall
- Excess water should be wiped off grab bar after each use to prevent rusting



[Shop Now](#)

## Portable Bath Step

Item # RTL12068

- Makes getting in and out of the bath tub easier and safer
- Unique modular design allows the multiple steps to be stacked
- Step height is 4"
- When stacking, each additional step adds 2"
- Bath step comes with four removable covers on the surface
- When stacking, remove covers to allow steps to lock together
- Attractive, textured surface is slip resistant and comfortable



[Shop Now](#)

# Bathroom Safety

## Suction Cup Grab Bar

Item # RTL13083

- Provides additional stability when entering and exiting the tub or shower
- Easy to install or remove with locking latches that secure the grab bar into place
- Color indicators ensure safe and secure installation. Green display means the grab bar is secure and ready for use, red will be displayed if the seal is not airtight
- Featuring a rubber non-slip grip and molded handle finger grips on the underside for additional safety and security

[Shop Now](#)



## Suction Cup Grab Bar

Item # RTL13095

- Provides additional stability when entering and exiting the tub or shower
- Easy to install or remove with locking latches that secure the grab bar into place
- Color indicators ensure safe and secure installation. Green display means the grab bar is secure and ready for use, red will be displayed if the seal is not airtight
- Featuring a rubber non-slip grip and molded handle finger grips on the underside for additional safety and security

[To purchase, reach out to your sales rep](#)



# Bathroom Safety

## Clamp-On Tub Rail

Item # RTL12036-ADJ

- Durable, white powder-coated steel construction
- 14.5" - 17" height adjustment
- Width adjusts from 3" - 7"
- Stainless steel locking mechanism and non-wearing rubber pads provide scratch-proof security
- Fits most tubs (Tighten carefully when using with Fiberglass tubs)
- Tool-free installation

[Shop Now](#)



## Plastic Tub Rail

Item # 12044

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip
- Handle is perpendicular to tub
- 16.5" in height
- Tool-free, easy to install
- Width adjusts from 3" - 6.25"
- Made using recycled materials

[Shop Now](#)



## M+D Multi-Use Tub Rail

Item # MD785340WH



- Two independent gripping surfaces allow this safety rail to adapt to a wide variety of tubs.
- Rust-resistant steel construction offers solid support up to 300 lbs.
- Rotating handle for multi-directional support
- Non-marring rubber pads on all contact points for a solid yet bathtub-friendly grip.
- Fits bathtub walls between 3" and 7" wide
- Installs securely without the use of tools.

Coming Soon



## M+D Massage Bath Mat

Item # MD785500WH

- Innovative design combines safety for the whole family with invigorating massage zones
- Gentle, low profile, textured surface that allows you to be comfortable standing or sitting
- Innovative design accommodates specialized bath seating without compromising stability
- Contoured ends allow for better tub coverage without blocking the drain
- A full-length bath mat with over 400 suction cups to provide an active, secure grip that conforms to your tub, making it safer for everyone

Coming Soon



# Bathroom Safety

## Deluxe Handheld Shower Spray with Diverter Valve

Item # 12037

- On/Off switch built into handle for easy access
- Extra-long 80" white reinforced nylon hose provides extra convenience
- Comes with diverter valve that allows either the regular shower head to operate or be used as a handheld shower spray
- Useful when shower chair, bench or stool is being used

[Shop Now](#)



## Deluxe Handheld Shower Massager with Three Massaging Options

Item # RTL12045

- Adjustable massaging head
- Chrome shower-arm mount and extra-long 84" chrome/mylar reinforced tangle-free hose
- Can be used as handheld shower spray or used in handle mount as traditional showerhead
- Useful when shower, chair, bench or stool is being used

[Shop Now](#)



# Bathroom Safety

## Bellavita Dive Bath Lift

Item # 477400252

- New removable backrest with 50° recline
- Larger seat with hygiene cutout
- Improved dual column design provides optimal stability for up to 300 pounds
- Seat height adjusts from 2.6" to 18.9"
- Easy-to-use floating remote control
- Installs easy in almost any bathtub—no plumbing or electrical work required



[Shop Now](#)

## Whisper Bath Lift Chair

Item # 477150312

- Designed to fit just about any style of bathtub, the Whisper shower chair for bathroom is easy to install and can support up to 300 pounds
- Durable plastic and is comfortable thanks to an ergonomically designed padded seat and backrest
- Designed to fit your unique needs, our bathtub patient lift can be adjusted to raise, lower, and recline the lift up to 50 degrees
- The bath chair's unique controls feature large, color-coded buttons for easy use; control floats in water for easy accessibility



[To purchase, reach out to your sales rep](#)

# Bathroom Safety

## Bellavita Comfort Cover

Item #	Color
460900252	White
460900312	Blue

- Cover set for Bellavita bath lift features a backrest and seat cover that come in an attractive color to match your bathroom decor
- The seat and backrest of the bath chair lift are padded for added comfort during while bathing
- Covers can be removed and washed separately to maintain a hygienic environment
- Easy to Install: Easily enhance your bath lift with a perfect fit

[Shop Now](#)



## Bellavita Comfort Cover, Blue

Item # 460900312

- Cover set for Bellavita bath lift features a backrest and seat cover that come in an attractive color to match your bathroom decor
- The seat and backrest of the bath chair lift are padded for added comfort during while bathing
- Covers can be removed and washed separately to maintain a hygienic environment
- Easy to Install: Easily enhance your bath lift with a perfect fit

[Shop Now](#)



# Bathroom Safety

## Vitaturm Swivel and Transfer Aid

Item # 412100200

- Simplifies transfer, positioning and rotating when using a bath lift
- Soft, mobile seat surface glides smoothly
- Slip-resistant surface improves safety
- The skin-friendly coating is easy-to-clean

[Shop Now](#)



## Universal Bather Pouch

Item # 10267-1

- Universal pouch for bath & shower chairs and benches
- Mesh pocket to securely hold personal items soap, shampoo, etc
- Easy-to-clean nylon

[Shop Now](#)





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  [www.drivemedical.com](http://www.drivemedical.com)

CANADA OFFICE  256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7  866-752-1363

©2026 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.