



# CPAP CATALOG





## CPAP PRODUCTS

Roscoe® has consistently offered tools that help sleep professionals provide better care and achieve better results. We are committed to help you break through emerging clinical and industry barriers by offering more options in CPAP Therapy than ever before. Find everything you need at Roscoe to manage your patients and your business.

# CPAP Products

Patient Interfaces	3-6
CPAP Systems	7
CPAP Private Label	8
CPAP Filters	9-10
CPAP Tubing	11-12
CPAP Replacement Headgear	13-14
CPAP Accessories	15-18

# DreamEasy™ 2 Full Face CPAP Mask and Fit Pack

The DreamEasy 2 Full Face CPAP Mask provides a superior seal with its removable Comfort Cushion that conforms to the facial contours around the nose and mouth. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. The unique keyhole design forehead bar provides stability while eliminating pressure. Quick-release clips allow for easy removal of the mask without re-adjusting the headgear. Dishwasher safe for easy cleaning (top rack only).



## FEATURES & BENEFITS

### 360° Rotational Elbow

Allows for more freedom of movement when sleeping with standard 22 mm tubing.

### Eliminate Pressure

The unique keyhole design forehead bar provides stability while eliminating pressure.

### Quick-Release

Clips allow for easy removal of the mask without re-adjusting the headgear.

### Dishwasher Safe

for easy cleaning (top rack only).

### Full Face Design

Perfect for patients whose therapy is compromised by mouth breathing.

### Lightweight Design

Durable and lightweight designed face mask with hook-and-loop fastened headgear allows for maximum movement.

### Comfort Cushion

Made with medical liquid silicone, the cushion is removable, odorless and provides an airtight seal.

### CO<sub>2</sub> Washout

Elbow vent holes provide effective CO<sub>2</sub> washout.



Compass Health offers a 30-day, no questions asked, money-back guarantee on all Roscoe CPAP Masks, effective from the original date of purchase.

Item #	Description	Size	HCPCS Code	Qty
CPM-DEF2S	DreamEasy 2 Full Face CPAP Mask with Headgear	Small	A7030	Each
CPM-DEF2M	DreamEasy 2 Full Face CPAP Mask with Headgear	Medium	A7030	Each
CPM-DEF2L	DreamEasy 2 Full Face CPAP Mask with Headgear	Large	A7030	Each
CPM-DEF2SK	DreamEasy 2 Full Face CPAP Mask with Headgear	All Sizes Fit Pack	A7030	Each
X1DEF2HS01	DreamEasy 2 Full Face Replacement Seal	Small	A7031	Each
X1DEF2HS02	DreamEasy 2 Full Face Replacement Seal	Medium	A7031	Each
X1DEF2HS03	DreamEasy 2 Full Face Replacement Seal	Large	A7031	Each
X1DEF2HS04	DreamEasy 2 Full Face Replacement Headgear, Gray	One Size	A7035	Each
X1DEF2HS05	DreamEasy 2 Full Face Replacement Headgear Clips	N/A	N/A	Pair

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

## Roscoe Medical® DreamEasy™ Nasal Mask

The DreamEasy™ Nasal CPAP Mask provides a superior seal with its removable Comfort Cushion™ that conforms around the nose and facial contours for a comfortable night's rest. Equipped with a 360° rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping. Comfortable forehead silicone pad diffuses any pressure. Quick-release clips allow for easy removal of the mask without readjusting the headgear. Dishwasher safe for easy cleaning. Replacement Comfort Cushion seals available.



### PRODUCT FEATURES

- Flexible and lightweight design allows for a comfortable, secure fit
- Headgear with Velcro® fasteners provide maximum versatility and movement
- 360° rotational elbow with standard 22 mm tubing connector
- The Comfort Cushion™ made with medical liquid silicone is completely odorless
- Comfort Cushion can be removed for easy cleaning
- Ergonomic design is airtight
- Vent hole provides effective CO<sub>2</sub> washout
- Latex free
- Mask is dishwasher safe (top rack) for easy cleaning
- 90-day limited warranty

ITEM #	DESCRIPTION	QTY.
CPM-DENS	DreamEasy™ Nasal Mask with Headgear, Small	each
CPM-DENM	DreamEasy™ Nasal Mask with Headgear, Medium	each
CPM-DENL	DreamEasy™ Nasal Mask with Headgear, Large	each
CPM-DENSK	DreamEasy™ Nasal Mask with Headgear, Starter kit	each
CPAP-SZGD	Nasal/Full Face Mask Sizing Guide	each
90091	DreamEasy™ Replacement Nasal Seal, Small	each
90093	DreamEasy™ Replacement Nasal Seal, Medium	each
90095	DreamEasy™ Replacement Nasal Seal, Large	each
90395	Replacement Headgear Clips	2/pk
90377	Replacement Headgear	each

### BILLING TIPS (see page 19)

#### HCPCS CODES

- A7034: Mask Interface, Nasal or Cannula Style
- A7032: Replacement Nasal Cushion
- A7035: Headgear

## Roscoe Medical® Zzz - Mask Full Face CPAP Mask

The Roscoe® Full Face CPAP Mask's one piece design is flexible, lightweight, and very easy to clean. The soft medical grade silicone mask allows for a comfortable fit and excellent seal without the need of a forehead support. Equipped with a 360 degree flexible rotational elbow and standard 22 mm tubing connector allowing for freedom of movement when sleeping.



### PRODUCT FEATURES

- Made of high-quality, lightweight, clear silicone with a flexible frame
- Dual 360° swivel at the tubing attachment point
- Headgear with Velcro® fasteners provide maximum versatility
- Latex free
- Dishwasher safe (top rack) for easy cleaning
- Full face mask comes complete with headgear (black)
- 90-day limited warranty

ITEM #	DESCRIPTION	QTY.
PB7800S	ZZZ-Mask Full Face Mask with Headgear, Small	each
PB7800M	ZZZ-Mask Full Face Mask with Headgear, Medium	each
PB7800L	ZZZ-Mask Full Face Mask with Headgear, Large	each
CPAP-SZGD	Nasal/Full Face Mask Sizing Guide	each
PB780S	ZZZ-Mask Full Face Seal, Small	each
PB780M	ZZZ-Mask Full Face Seal, Medium	each
PB780L	ZZZ-Mask Full Face Seal, Large	each

### BILLING TIPS (see page 19)

#### HCPCS CODES



- A7030: Mask Interface, Full Face Style
- A7035: Headgear
- A7031: Replacement Full Face Cushion

# Roscoe Medical® Shadow™ Nasal Pillows CPAP Mask

The Roscoe® Shadow™ Nasal Pillows CPAP Mask features a sleek design that conforms naturally to the contours of the face providing a comfortable fit with minimal facial contact. This permits for almost any sleep position without mask displacement. This innovative mask allows the user to anchor the tube at the top of the headgear or to leave it free to track. The anatomically shaped soft silicon pillows range in size from extra-small to large for a perfect fit. A complete CPAP solution right out of the box!  
**Replacement seals available.**

## PRODUCT FEATURES

- Lightweight sleek design conforms to facial contours with minimal contact for maximum comfort and compliance
- Anatomically shaped soft silicone pillows provide maximum comfort with an excellent seal
- 360° rotational elbow design allows for tubing can be anchored at the top of the headgear or left free to track
- A built-in air vent disperses air quietly to prevent sleep disturbance
- Fit pack comes complete with 4 sets of pillows: extra-small, small, medium and large
- Replacement pillows available
- Latex free
- 90-day limited warranty



ITEM #	DESCRIPTION	QTY.
70116	Shadow™ Nasal Pillows Mask Fit Pack with 4 Sets of Pillows and Headgear	each
90385	Shadow™ Replacement Pillows, extra small	pair
90387	Shadow™ Replacement Pillows, small	pair
90389	Shadow™ Replacement Pillows, medium	pair
90391	Shadow™ Replacement Pillows, large	pair

BILLING TIPS (see page 19)
HCPCS CODES
A7034: Mask Interface, Nasal or Cannula Style
A7033: Replacement Pillows
A7035: Headgear

## Sleepnet® Ascend® CPAP Mask Series

The Sleepnet® Ascend® CPAP Mask Series features advanced cushion technology with the replaceable AIR°gel® cushion. For added convenience, patients can replace the cushion without having to replace the whole mask. The patented AIR°gel provides maximum comfort and minimum weight helping to reduce pressure points and skin irritation for better patient compliance. The flexible spacebar has a unique bendable design that provides the patient with superior stability without contacting the skin.



Ascend® Nasal Mask

Ascend® Full Face Mask

### PRODUCT FEATURES

- Detachable cushion allows you to replace your cushion without replacing your entire mask
- Flexible spacebar is designed to provide superior stability without contacting the skin
- Active headgear connector moves when the patient moves to help maintain an effective seal
- Available in 3 sizes (small, medium and large) with simple and breathable quick-release EZ Fit headgear

- 360° rotating swivel for optimal tube positioning and easier mobility
- Patented AIR°gel® with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal
- Available in your choice of full face or nasal mask

ITEM #	DESCRIPTION	QTY.
50825	Sleepnet® Ascend® Full Face Mask System with Headgear, small, medium and large cushion included	each
50207	AIR°gel® Full Face Replacement Cushion, Small	each
50208	AIR°gel® Full Face Replacement Cushion, Medium	each
50209	AIR°gel® Full Face Replacement Cushion, Large	each
50174	Sleepnet® Ascend® Nasal Mask System with Headgear, small, medium and large cushion included	each
50201	AIR°gel® Nasal Mask Replacement Cushion, Small	each
50202	AIR°gel® Nasal Mask Replacement Cushion, Medium	each
50203	AIR°gel® Nasal Mask Replacement Cushion, Large	each

### BILLING TIPS (see page 19)

#### HCPCS CODES

A7034: Mask Interface, Nasal or Cannula Style

A7030: Mask Interface, Full Face Style

A7031: Replacement Full Face Cushion

A7032: Replacement Nasal Cushion

A7035: Headgear

## Sleepnet® MiniMe® 2 Pediatric Nasal Mask

Designed specifically for children, the MiniMe® 2 is available in two sizes for a greater range of comfort for children between ages two and twelve. Enhanced Custom Fit Technology™ enables the mask to be easily molded to a child's face maximizing comfort. The MiniMe 2 also features Advanced Cushion Technology™ (A.C.T.), Sleepnet's new ergonomic cushion technology delivers an effective seal yet minimizes pressure on small faces.



### PRODUCT FEATURES

- Specifically engineered with children in mind, the Sleepnet® MiniMe® 2 provides unmatched comfort, a non-invasive and customized design, easy adjustment, and a pediatric size frame for the extra small faces we care so much about
- Touchless Spacebar™ eliminates any contact with the patient's forehead
- Custom Fit Technology™ means the shell of the mask can be gently molded to match your needs

- Designed specifically for children, the MiniMe® 2 is available in 2 sizes for a greater range of comfort for children ages 2-12
- Lightweight, breathable headgear
- Designed and built in the USA
- Does not contain any natural rubber latex

ITEM #	DESCRIPTION	QTY.
60217	MiniMe® 2 Pediatric Mask with Headgear, Small/Medium	each
60218	MiniMe® 2 Pediatric Mask with Headgear, Medium/Large	each

### BILLING TIPS (see page 19)

#### HCPCS CODES

A7034: Mask Interface, Nasal or Cannula Style

A7035: Headgear

## Resvent iBreeze Auto CPAP & Bilevel Systems



Introducing the Resvent iBreeze Auto CPAP and Bilevel Systems with Heated Humidification, feature-rich PAP devices that were granted EUA (Emergency Use Authorization) by the FDA in 2020.\* Some of the features include IPR (Intelligence Pressure Release), an algorithm during inspiration and expiration offering patients a more comfortable breathing experience; an auto humidifier which adjusts the heat based on room temperature and humidity to prevent rainout in the hose or mask; and finally, a 2-in-1 humidifier design which monitors the water level and when low, will signal a visible alert and automatically switch off the heat to prevent damage to the device. The iBreeze Systems come complete with SD Card, Carry Bag, 6' CPAP Tube, Tube Connector, Power Cord, and User Manual. *Wi-Fi equipped devices optional.*

### SPECIFICATIONS

H2O Therapy Pressure:	4-20 cm (APAP) 4-30 cm (BPAP)
Ramp Time:	0-60 minutes (5-minute increments)
Dimensions:	9.25" w x 5.8" d x 4.9" h (device); 14.25" x 9.75" x 7" (shipping)
Weight:	3.4 lbs (device); 6.2 lbs (shipping)
Water Capacity:	290 mL
Humidifier Settings:	Auto / 0-8 (Default: 3)
Humidifier Output:	≥10 mg H2O/L
Warranty:	Two Years Limited (device)



### FEATURES & BENEFITS

<p><b>Large 3.5" LCD Display</b> Easily adjust settings and view last session to encourage user compliance.</p>	<p><b>Quiet Operation</b> with real time audible alert for tubing/mask leaks and automatic leakage/altitude compensation.</p>
<p><b>Intelligence Pressure Release</b> with 3 comfort settings and inspiration trigger for automatic startup.</p>	<p><b>Compliance Data</b> access compliance reporting via SD card or optional Wi-Fi.</p>
<p><b>Embedded Memory</b> stores previous night and 365 night user raw data.</p>	<p><b>iMatrix/ResAssist</b> Register your patients to maintain records and generate sleep reports on iMatrix desktop software or ResAssist cloud platform.</p>

ITEM	DESCRIPTION	HCPCS	UOM	QTY
IBREEZE20A	iBreeze Auto CPAP System with Heated Humidification	E0601	Each	1
IBREEZE20A-W	iBreeze Auto CPAP System with WiFi and Heated Humidification	E0601	Each	1
IBREEZE25STA	BiLevel Auto S/T Stem (4-25) with WiFi and Heated Humidification	E0471	Each	1
IBREEZE30STA	BiLevel Auto S/T System (4-30) with WiFi and Heated Humidification	E0471	Each	1
<b>Replacement Parts</b>				
05205-000001-00	iBreeze Replacement 2-Piece Power Adapter with Brick	N/A	Each	1
10001-010067-01	iBreeze Replacement Humidifier Chamber	A7046	Each	1
CPF-IBREEZE	iBreeze Replacement Disposable Filter	A7038	Pack	2
69999-000032-00	iBreeze APAP Replacement Knob	N/A	Each	1
69999-000033-00	iBreeze BPAP Replacement Knob	N/A	Each	1
10001-010043-00	iBreeze Replacement Tube Adapter	N/A	Each	1

# Roscoe Medical® Private Label Program

With Roscoe Medical's CPAP Private Label Program, you have the ability to build your brand awareness with CPAP replenishment supplies and improve your repeat patient business. Create a professional image and improve patient loyalty by private labeling your CPAP tubing, filters and customized tubing/filter kits. Include an insert with helpful replacement tips and related cash items, all branded to your company with your information. You will enjoy repeat supply orders and additional cash sales directly from your patients.

## PROGRAM BENEFITS

- Low minimum order quantity of 50 each/pack with quick turnaround time
- Grow your repeat business for replenishment supplies
- Boost cash sales for CPAP-related accessories (Volt™ Battery Pack, pillows, wipes, REMZzzz's, hose wraps, etc...)
- Helps maintain profitability in the wake of Round 2 reimbursement cuts
- Provides the ability to brand replenishment supplies such as:
  - › CPAP Tubing
  - › Filters
  - › Tubing/Filter Kits
- Create a professional image for your business
- Include helpful replacement tips informing patients how often tubing and filters should be replaced



Sample Private Label




DreamEasy2 Full Face Mask & Easy-Flex6 Tube







Ascend Full Face Mask, Easy-Flex6 Tube, O2 Bleed-in Adapter & 25' Oxygen Tube

SAMPLE PRIVATE LABEL KIT		QTY
CPM-DEF2SK	DreamEasy2 Full Face Mask	1
EASY-FLEX6	Easy-Flex 6' CPAP Tubing	1
50825	Sleepnet Ascend Full Face Mask w/ EZFit Headgear	1
AG1642	O2 Bleed-in Adapter	1
TUB-ROS25G	25' Oxygen Tube Green	1




## iBreeze Filters For Resvent Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
CPF-IBREEZE	 Resvent	2/pk	6.1" x 3.94" x 0.47"	N/A	N/A	N/A	A7038

## Disposable Filters For Philips Respironics Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
CPF-DSD2	 DreamStation	2/pk	2 1/4" x 1"	1122518	N/A	N/A	A7038
CPF-PR1	 PR System One	2/pk	1 3/4" x 7/8"	1063096	AG1063096MED	CF1009NT-2	A7038
CPF-F30PK2	 Remstar Pro M Series, SleepEasy	2/pk	1 3/4" x 1 1/4" (includes tab)	1029331	AG1029331M	CF1007-2	A7038
CPF-F19	 Remstar Plus, Pro, BiPAP Pro 2 with Bi-Flex, Remstar Autosmart	2/pk	3 5/8" x 2" (includes tab)	1005945	AG5945	CF1006	A7038


## Reusable Filters For Philips Respironics Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
CPF-DSR1	 Dream Station	1/pk	2 1/4" x 7/8" x 5/8"	1122446	A7039	N/A	N/A
CPF-F31PK1	 PR System One, Remstar M Series, SleepEasy	1/pk	1 3/4" x 7/8"	1029330	AG1029330M	CF1007F-1	A7039
CPF-F28PK1	 Remstar Plus, Pro, BiPAP Pro 2 w/Bi-Flex, Remstar Autosmart, Legacy C-Flex	1/pk	3 5/8" x 1 5/8"	1005964	AG5964M	CF1006F	A7039



## Filters For ResMed® Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
CPF-A11	 Airsense 11	2/pk	4.65" x 3.78" x 1.5"	36851	AG36851	CF2107-2	A7038
CPF-S9PK2	 S9 Series, AirSense 10, AirCurve 10	2/pk	2 1/8" x 1 7/16"	36851	AG36851	CF2107-2	A7038
CPF-S9HYPK2	 S9 Series Hypoallergenic, AirSense 10, AirCurve 10	2/pk	2 1/8" x 1 7/16"	36856	AG36856	CF2107HA	A7038
CPF-S8PK2	 S8 Series	2/pk	1 5/16" x 1 3/8"	33916	AF33916MED	CF2105-2	A7038


## Filters for DeVilbiss® Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
ULTRA FINE FILTER							
CPF-F32PK2	 IntelliPAP Standard, Standard Plus, Auto Adjust and Bi-level S	2/pk	1 5/8" x 1 3/16"	DV51D-603	F603MED	CF5004	A7038

## Filters for 3B Medical Units

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
LG1020	 Luna/Luna II	6/pk	2" x 1 1/4"	Same	N/A	CF4003-2	A7038
CPF-3BPK2	 RESmart	2/pk	3 1/4" x 1 5/8"	N/A	N/A	CF8004	A7038

## Bacteria Filters

ITEM #	UNIT	QTY	DIMENSIONS	OEM	AG	SUNSET	HCPCS CODE
AG7178	 Multiple	each	22 mm	Same	Same	BF2200	N/A



## CPAP TUBING & ACCESSORIES

Our Roscoe Medical® high quality CPAP Tubing is offered in a variety of styles. The tubing features a smooth interior which allows for excellent air flow and is compatible with virtually all CPAP, APAP and Bi-Level Systems. We also offer a wide range of headgear and accessories to meet all your patient's needs.

*\*Please note: All Roscoe Medical® CPAP tubing is intended for single-patient use.*

## Easy-Flex™ Lightweight CPAP Tubing

The Roscoe Medical® Easy-Flex™ Lightweight CPAP Tubing has a soft, flexible feel providing patients with the high quality they expect. Equipped with 22 mm cuffs, our improved tubing is lighter and more flexible than comparable tubing. This tubing is dark gray in color and available in a variety of lengths.



#EASY-FLEX6

ITEM #	DESCRIPTION	QTY.
EASY-FLEX6*	Easy-Flex™ Lightweight Tubing, 6 ft. Dark Gray with 22 mm Cuffs	each
EASY-FLEX8*	Easy-Flex™ Lightweight Tubing, 8 ft. Dark Gray with 22 mm Cuffs	each
EASY-FLEX10*	Easy-Flex™ Lightweight Tubing, 10 ft. Dark Gray with 22 mm Cuffs	each
EASY-FLEX18*	Easy-Flex™ Lightweight Tubing, 18 in. Dark Gray with 22 mm Cuffs	each

BILLING TIPS (see page 19)

HCPCS CODES

A7037: CPAP Tubing

## FlexiSlim™ CPAP Tubing

The Roscoe Medical® FlexiSlim™ CPAP Tubing is a small diameter air tube that is exceptionally light and flexible. FlexiSlim tubing works well with smaller lightweight masks as it is less likely to pull on the mask and break the seal.



### PRODUCT FEATURES

- 6 ft. dark gray, small diameter air tube
- Ergonomic 22 mm cuffs for easy installation and removal from CPAP unit
- Exceptionally lightweight and flexible—40% lighter than standard tubing
- FlexiSlim™ tubing is 10% thinner than standard CPAP tubing
- Patient can freely move around in bed without creating mask leaks
- Tubing is designed to work with ResMed S9/S10 series and any devices with 15 mm tubing settings
- 90-day warranty

ITEM #	DESCRIPTION	QTY.
SLIM-FLEX6	FlexiSlim™ CPAP Tubing, 6 ft. Dark Gray with 22 mm Cuffs	each

BILLING TIPS (see page 19)

HCPCS CODES

A7037: CPAP Tubing

## CPAP REPLACEMENT HEADGEAR & ACCESSORIES

Patients may need to replace their headgear occasionally as it may become stretched or lose their elasticity. For optimal performance, provide your patients your patients with a personalized replacement plan to ensure optimal performance and help keep them compliant.

### CPAP Replacement Headgear



#ROS-FFHG



#ROS-NHG

BILLING TIPS (see page 19)

HCPCS CODES

A7035: Headgear

ITEM #	DESCRIPTION	QTY.
X1DEF2HS04	DreamEasy 2 Full Face Replacement Headgear, Gray	each
X1DEF2HS05	DreamEasy 2 Full Face Replacement Headgear Clips	pair
90377	Roscoe Medical® Headgear, for use with Roscoe® DreamEasy™ Full Face and Nasal CPAP Masks, Black	each
ROS-NHG	Universal Nasal Mask Headgear, Black	each
90395	Roscoe Medical® Replacement Headgear Clips, for use with Roscoe® DreamEasy™ Full Face and Nasal CPAP Masks	2/pk

## CPAP Replacement Chin Straps



#ROS-T09M

ITEM #	DESCRIPTION	QTY.
ROS-T09S	3-Point Chin Strap, Small (compare to Tiara/Cardinal #TMS-09S), Ruby Red	each
ROS-T09M	3-Point Chin Strap, Medium (compare to Tiara/Cardinal #TMS-09M), Ruby Red	each
ROS-T09L	3-Point Chin Strap, Large (compare to Tiara/Cardinal #TMS-09L), Ruby Red	each
ROS-T09ADJ	3-Point Adjustable Chin Strap,with adjustable 2 in. extension (compare to Tiara/Cardinal #TMS-09ADJ), Ruby Red	each
ROS-T09ADJXL	3-Point Adjustable Chin Strap, Extra Large, Ruby Red	each

## CPAP Mask Wipes

These 100% pure cotton wipes are specially formulated to naturally clean and protect your mask without the use of harmful chemicals. CPAP Mask Wipes effectively remove dirt, grease, oils and other residue from the patient’s mask keeping it fresh and clean. Wipes are made from 100% natural ingredients and contain no strong or unpleasant fragrances.

### PRODUCT FEATURES

- Contain only natural ingredients derived from coconut and other plant extracts
- Suitable for daily cleaning of CPAP masks, tubing and accessories
- Gentle on hands and face
- 100% pure cotton
- Alcohol and latex free
- Biodegradable



FGC13700

ITEM #	DESCRIPTION	QTY.
FGC13700 0000	Carex CPAP Mask Wipes, Unscented (62 wipes per tub)	each
FGC13800 0000	Carex CPAP Mask Wipes, Citrus Scent (62 wipes per tub)	each

## CPAP Tubing Hose Wraps

Roscoe Medical® Hose Wraps eliminate the most annoying parts of therapy—the cold medical tubing that connects the user’s system to the mask or interface. By adding our Roscoe® Hose Wraps to your patient’s system, you can help eliminate condensation and rainout that often disrupts sleep and therapy.

### PRODUCT FEATURES

- Lightweight and easy to travel with
- Ultra soft fleece material
- Easy to install and very patient friendly
- Fits any 6 ft. 22 mm CPAP tube
- Comes standard with easy grip closure
- Assists in controlling rainout
- Machine washable for easy cleaning
- Available in 2 colors: navy and burgundy



#HW-NAVY

ITEM #	DESCRIPTION	QTY.
HW-NAVY	Roscoe® Hose Wrap, Navy	each

## CPAP Adaptors, Port and Humidification Chambers



AG1642



HC325S



HC385S



10001-010067-01

ITEM #	DESCRIPTION	QTY.
10001-010067-01	iBreeze Replacement Humidifier	each
AG1642	Oxygen Bleed-In Port	each
HC325S	Humidification Chamber, for HC150 CPAP Systems, Single-Patient Reusable	each
HC385S	Standard Humidification Chamber for SleepStyle™ 200 CPAP Series	each

## Carex CPAP Pillow

The Carex CPAP Pillow can help reduce CPAP mask leaks, reduces mask pressure and helps increase comfort for a good night's rest. The larger cutouts on either side allow the user to sleep on either side without worrying about leaks.

### PRODUCT FEATURES

- Use with Full Face, Nasal or Direct Nasal interfaces
- Freedom to move without displacing mask
- Provides neck and head support
- Removable cover for easy cleaning



ITEM #	DESCRIPTION	QTY.
FGP11300 0000	Carex CPAP Pillow with Cover (box packaging)	each

## Contour CPAPMax Pillow

The CPAPMax Pillow can help reduce CPAP mask leaks, reduces mask pressure and helps increase comfort for a good night's rest. The larger cutouts on either side allow the user to sleep on either side without worrying about leaks.

### PRODUCT FEATURES

- Reduces mask interference, shifting and leaks
- Adjustable pillow height switched from regular to extra thick
- Minimizes facial pressure, sores and bruising
- Improve sleep comfort and CPAP compliance
- Ventilated foam is cool, breathes and reduces perspiration
- Ergonomic center cradles head and improves airway alignment
- Front edge matches the natural curve of neck and shoulders
- Pressure-free zones prevent interference between mask and pillow
- CoolTech 3D mesh cover
- 2-in-1 reversible memory foam or fiber top



ITEM #	DESCRIPTION	QTY.
PLW-CPMX	Contour CPAPMax Pillow	each

## CPAP HCPCS Codes

CODE	PRODUCT	REPLACEMENT SCHEDULE
A4604	Tubing with integrated heating element for use with positive airway pressure (PAP) device	1 per 3 months
A7027	Combination oral/nasal mask, used with continuous positive airway pressure (CPAP) device, each	1 per 3 months
A7028	Oral cushion for combination oral/nasal mask, replacement only, each	2 per month
A7029	Nasal pillows for combination oral/nasal mask, replacement only, pair	2 per month
A7030	Full face mask used with PAP device, each	1 per 3 months
A7031	Face mask interface, replacement for full face mask, each	1 per month
A7032	Replacement cushion for nasal application device, each	2 per month
A7033	Replacement pillows for nasal application device, pair	2 pairs per month
A7034	Nasal mask or cannula type application device, used with PAP device	1 per 3 months
A7035	Headgear used with PAP device	1 per 6 months
A7036	Chin strap used with PAP device	1 per 6 months
A7037	Tubing used with PAP device	1 per 3 months
A7038	Filter, disposable, used with PAP device	2 per month
A7039	Filter, non-disposable, used with PAP device	1 per 6 months
A7044	Oral interface used with PAP device	Contact payor
A7045	Exhalation port with or without swivel used with accessories for PAP devices, replacement only	Contact payor
A7046	Replacement water chamber for humidification, used with PAP device	1 per 6 months
A9279	Monitoring feature/device stand-alone or integrated, any type, includes all accessories, components and electronics, not otherwise classified	Contact payor
E0470	Respiratory assist device (RAD), bi-level pressure capability, without backup rate feature, used with noninvasive interface	Monthly rental rate
E0471	Respiratory assist device, bi-level pressure capability, with backup rate feature, used with noninvasive interface	Monthly rental rate
E0561	Humidifier, non-heated, used with PAP device	Contact payor
E0562	CPAP heated humidifier	Contact payor
E0601	CPAP device	Monthly rental rate

### \*HCPCS Disclaimer:

HCPCS codes provided by Compass Health® are intended as general guidelines. Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended nor should it be considered billing advice. Compass Health® does not guarantee coverage or reimbursement on any product.





6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

[www.compasshealthbrands.com](http://www.compasshealthbrands.com)

rev\_03\_251013