



DME CATALOG



- Bath Safety
- Mobility
- Seating and Positioning



Shower Chair without Back

The **ProBasics Shower Chair without Back** is designed for individuals who require seating support and stability when bathing or showering. The bench features a comfortable contoured textured seat that can support patient weights up to 300 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Blow-molded bench has a textured surface and is contoured with drainage holes to reduce slipping
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Two hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0245*

SPECIFICATIONS

Overall Footprint:	16.5"W x 19"D
Seat Dimensions:	19.5"W x 11.5"D
Adjustable Seat Height:	14.25" – 20.5"
Weight Capacity:	300 lb
Carton Dimensions:	23.5" x 17" x 19"
Shipping Weight:	19 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSC	ProBasics Shower Chair without Back	815067071920	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Shower Chair without Back

The **ProBasics Bariatric Shower Chair without Back** is designed for individuals who require seating support and stability when bathing or showering. The bench features a special cross brace design with a heavy duty blow-molded textured seat that can support patient weights up to 550 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Comfortable, blow-molded seat with integrated handles for additional support
- Dual cross brace design reinforced with aircraft-style rivets for extra stability
- High-strength aluminum frame to prevent corrosion
- Textured contoured surface with drainage holes to reduce slipping and increase comfort
- Wide-based tips help increase stability
- Legs are height adjustable in one-inch increments with heavy duty height adjustment pins to accommodate larger individuals
- Quick tool-free assembly
- HCPCS Code: E0245*

SPECIFICATIONS

Overall Footprint:	21"W x 17.5"D
Seat Dimensions:	20"W x 12"D
Adjustable Seat Height:	16.5" – 20.5"
Weight Capacity:	550 lb
Carton Dimensions:	18.5" x 22.5" x 17.5"
Shipping Weight:	12 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBC	ProBasics Bariatric Shower Chair w/o Back	815067071944	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Shower Chair with Back

The **ProBasics Shower Chair with Back** is designed for individuals who require seating support and stability when bathing or showering. The chair features a comfortable contoured textured seat and back that can support patient weights up to 300 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Blow-molded bench and back provide comfort and strength
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Two hole handles in seat provide security when sitting or standing
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0240*

SPECIFICATIONS

Overall Footprint:	16.5"W x 19"D
Seat Dimensions:	19.5"W x 11.5"D
Adjustable Seat Height:	14.25" – 20.5"
Weight Capacity:	300 lb
Carton Dimensions:	28" x 21" x 19"
Shipping Weight:	25 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSCWB	ProBasics Shower Chair with Back	815067071913	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Shower Chair with Back

The **ProBasics Bariatric Shower Chair with Back** is designed for individuals who require seating support when bathing or showering. The chair features a special cross brace design with a heavy duty blow-molded back and textured seat that can support patient weights up to 550 pounds. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Comfortable, blow-molded seat with integrated handles for additional support
- Textured contoured surface with drainage holes to reduce slipping and increase comfort
- Dual cross brace design reinforced with aircraft-style rivets for extra stability
- High-strength aluminum frame to prevent corrosion
- Legs are height adjustable in one-inch increments with heavy duty height adjustment pins to accommodate larger individuals
- Wide-based tips help increase stability
- Backrest is removable
- Quick tool-free assembly
- HCPCS Code: E0240*

SPECIFICATIONS

Overall Footprint:	21"W x 17.5"D
Seat Dimensions:	20"W x 12"D
Adjustable Seat Height:	16.5"– 20.5"
Weight Capacity:	550 lb
Carton Dimensions:	21" x 17" x 27"
Shipping Weight:	16.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBCWB	ProBasics Bariatric Shower Chair with Back	815067071937	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Deluxe Shower Chair with Back and Arms

The **ProBasics Deluxe Shower Chair with Back and Arms** is designed for individuals who require seating support and stability when bathing or showering. The chair features a comfortable contoured textured seat and back that can support patient weights up to 300 pounds. Arms and back are easy to remove. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Removable padded armrests provide comfort and added security when sitting or standing
- Blow-molded bench and back provide comfort and strength
- Quick tool-free assembly
- Wide-based crutch tips help increase stability
- Legs are height adjustable with one-inch snap button adjustments
- HCPCS Code: E0240*

SPECIFICATIONS

Overall Footprint:	18.5"W x 21"D
Seat Dimensions:	16"W x 13.5"D
Adjustable Seat Height:	17.5" – 21.5"
Weight Capacity:	300 lb
Carton Dimensions:	19" x 17" x 22"
Shipping Weight:	33 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSCWBA	ProBasics Deluxe Shower Chair w/ Back and Arms	815067071951	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Transfer Bench

The **ProBasics Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame and blow-molded bench, it features a reversible seat back for right or left handed transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, sturdy and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0247*

SPECIFICATIONS

Overall Footprint:	28.5"W x 17.5"D
Seat Dimensions:	25"W x 16"D
Adjustable Seat Height:	19" – 23"
Weight Capacity:	300 lb
Carton Dimensions:	28" x 17" x 11.5"
Shipping Weight:	23 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSTB	ProBasics Transfer Bench	815067071968	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Transfer Bench

The **ProBasics Bariatric Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame, it features added cross braces for stability and a heavy duty blow-molded seat. It converts easily for left or right hand transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, durable and corrosion-proof
- Durable blow-molded bench and back have a textured surface and are contoured with drainage holes to reduce slipping
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Added cross brace for stability
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0248*

SPECIFICATIONS

Overall Footprint:	32.5"W x 21.5"D
Seat Dimensions:	25"W x 15.5"D
Adjustable Seat Height:	18" – 24.25"
Weight Capacity:	500 lb
Carton Dimensions:	33.5" x 20.5" x 17"
Shipping Weight:	28 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBTB	ProBasics Bariatric Transfer Bench	815067071975	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Padded Transfer Bench

The **ProBasics Padded Transfer Bench** helps make transfers in and out of the bathtub safer and more comfortable. Constructed with an aluminum anodized frame and padded back and bench, it features a reversible back for right or left handed transfers. The angled legs have one-inch, easy-snap button adjustments and wide-based, non-skid feet for added sturdiness and safety. Tool-free assembly and limited lifetime warranty are standard.

FEATURES

- Aluminum frame is lightweight, sturdy and corrosion-proof
- Vinyl seat and back are padded with 1.25" open-cell foam for added comfort
- Wide-based tips and suction cups add sturdiness and prevent slipping
- Legs have one-inch easy snap-button adjustments
- Converts easily for left or right hand transfers
- Quick tool-free assembly
- HCPCS Code: E0247*

SPECIFICATIONS

Overall Footprint:	27" W x 16.5" D
Seat Dimensions:	25" W x 15.5" D
Adjustable Seat Height:	21" – 25"
Weight Capacity:	300 lb
Carton Dimensions:	30" x 17" x 14.5"
Shipping Weight:	25 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSTBP	ProBasics Padded Transfer Bench	815067071982	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Raised Toilet Seats

ProBasics Raised Toilet Seats are a must for anyone with bending or sitting difficulties. Available with or without arms and locking options to accommodate all needs. All toilet seats are constructed of heavy duty blow-molded plastic that is durable and easy to clean.



BSRTSL
Locking Seat



BSRTSLA
Locking Seat with Padded Arms

FEATURES

- Adjustable front locking mechanism tightens easily with one hand and fits most toilets. Prevents shifting of seat when transferring
- Elevates existing seat height by 4.5"
- Extra-large hygienic hole cutout
- Tool-free assembly
- Color: white
- HCPCS Code: E0244*

SPECIFICATIONS

Seat Width:	17"
Seat Depth:	17"
Width Between Arms:	19.25" (BSRTSLA only)
Seat Height:	BSRTSL: 4.5" BSRTSLA: 4.5"
Weight Capacity:	350 lb
Carton Dimensions:	BSRTSL: 17.5" x 7.5" x 18" BSRTSLA: 21" x 8" x 18"
Shipping Weight:	BSRTSL: 5.5 lb BSRTSLA: 7.3 lb
Warranty:	Limited lifetime on blow-molded seat One year on all parts

Item #	Description	UPC	UOM
BSRTSL	ProBasics Raised Toilet Seat with Lock	815067072002	1/cs
BSRTSLA	ProBasics Raised Toilet Seat with Lock and Padded Arms	815067072019	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Toilet Safety Frame

The **ProBasics Toilet Safety Frame** provides a safe, stable aid for individuals who struggle with the up and down demands of toileting. The frame is designed to accommodate both standard and elevated toilet seats. It is easy to clean and the aluminum frame mounts securely onto the bowl with an adjustable bracket. The toilet safety frame can safely support patient weights up to 300 pounds.

FEATURES

- Frame is made of anodized, non-corrosive one-inch aluminum tubing
- Arms are height and width adjustable
- Versatile design with tool-free assembly quickly secures to most toilets
- Contoured waterfall armrests provide additional comfort and support while providing an excellent gripping surface
- Powder-coated aluminum bracket easily attaches frame to toilet
- HCPCS Code: E0243*

SPECIFICATIONS

Adjustable Arm Width:	20" – 22"
Adjustable Arm Height:	26" – 29.5"
Weight Capacity:	300 lb
Carton Dimensions:	19.5" x 18" x 10"
Shipping Weight:	15 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSTF	ProBasics Toilet Safety Frame	815067071821	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Three-in-One Steel Commode

The **ProBasics Three-in-One Steel Commode** has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. It features a removable back and comes complete with commode bucket, splash guard and lid. The versatile design includes a deeper seat, armrests for added comfort, angled legs for stability and non-marring rubber tips.

FEATURES

- Can function as a free-standing commode, a toilet safety frame (with backrest removed), or as a raised toilet seat with arms
- Features a deeper seat for less mess and to allow more room for cleaning
- Durable steel-welded frame is attractive, strong and easy-to-clean
- Tool-free assembly, includes plastic snap on seat with lid
- Seat height readily adjusts in one-inch increments
- Plastic arms provide added comfort
- Comes complete with commode bucket with handle, cover and splash guard
- HCPCS Code: E0163*

SPECIFICATIONS

Width Between Arms:	17.5"
Outside Width:	24"
Seat Dimensions:	13.5"W x 16.5"D
Seat Opening:	8"W x 10"D
Adjustable Seat Height:	15.5" – 23.5"
Color:	Gray
Weight Capacity:	350 lb
Carton Dimensions:	25.5" x 25.5" x 25"
Shipping Weight:	57 lb
Warranty:	Limited lifetime on frame One year on parts

Item #	Description	UPC	UOM
BS31C	ProBasics Three-in-One Commode	815067071838	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Three-in-One Folding Commode

The **ProBasics Three-in-One Folding Commode** has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. This full-size steel commode folds compactly for storage and makes delivery easy. Clip-on seat can be removed for easy cleaning. Comes complete with a commode bucket, splash guard and cover.

FEATURES

- Durable steel-welded frame is attractive, strong and easy-to-clean
- Comfortable deep seat for ease of use
- Easily folds for transportation and storage
- Clip-on seat can be removed for easy cleaning
- Use as a bedside commode, toilet safety frame or raised toilet seat
- Seat height readily adjusts in one-inch increments
- Comes complete with seat, lid, commode bucket with handle, cover and splash guard
- Tool free assembly
- Color: gray frame, gray accessories
- HCPCS Code: E0163*

SPECIFICATIONS

Width Between Arms:	17.75"
Outside Width:	22"
Seat Dimensions:	13.5"W x 16.5"D
Seat Opening:	8"W x 10"D
Adjustable Seat Height:	16" – 22"
Weight Capacity:	350 lb
Carton Dimensions:	30" x 23" x 15"
Shipping Weight:	62.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSFC	ProBasics Three-in-One Folding Commode	815067071845	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Drop-Arm Commode

The **ProBasics Drop-Arm Commode** features a durable steel frame with an easy-release push-button arm mechanism allowing for safe lateral patient transfers to and from the commode. The anti-slip, soft Hypalon® armrests provide extra comfort and support. The frame can be configured to be used as a bedside commode, toilet safety frame or raised toilet seat. The tool free removable back makes configuring the commode easy for any application.

FEATURES

- Durable steel frame provides strength and stability
- Equipped with easy-to-release arm mechanism allowing for safe lateral transfers to and from the commode
- Push-button drop-arm allows easier side transfer
- Soft, slip-resistant Hypalon® covered armrests
- Seat height readily adjusts in one-inch increments
- Bucket slides easily into place
- Removable backrest
- Comes complete with 7.5 quart commode bucket with carry handle and splash guard
- Color: white frame, gray accessories
- HCPCS Code: E0163*

SPECIFICATIONS

Width Between Arms:	18"
Outside Width:	21"
Seat Dimensions:	13.5"W x 17"D
Seat Opening:	8"W x 10"D
Adjustable Seat Height:	17" – 22"
Weight Capacity:	300 lb
Carton Dimensions:	21.5" x 16.5" x 14"
Shipping Weight:	19.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSDAC	ProBasics Drop Arm Commode	815067071906	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Commode

ProBasics Bariatric Commode has a wider seat and steel frame to support a higher weight capacity and is ideal for larger patients. It has three versatile functions; it can be used as a toilet safety frame, a commode, and a raised toilet seat. It comes complete with commode bucket and seat lid. The versatile design includes a rectangular pan and pan holder, armrests for added comfort, angled legs for stability and non-marring rubber tips.

FEATURES

- Can function as a free-standing commode, a toilet safety frame (with bucket removed) or as a raised toilet seat with arms.
- Durable steel-welded frame is attractive, strong and easy-to-clean
- Tool-free assembly, includes plastic snap on seat with lid
- Plastic arms provide added comfort
- Seat height readily adjusts in one-inch increments
- Color: gray frame, white accessories
- HCPCS Code: E0168

SPECIFICATIONS

Seat height range:	15.75" – 20.87"
Width Between Arms	22"
Outside Width:	28.8"
Seat Opening:	8.39" x 10"
Pan Capacity:	10 quarts / 9.5 liters
Weight Capacity:	450 lb
Shipping Dimensions:	29.5" x 17.5" x 15.5"
Shipping Weight:	35.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSB31C	ProBasics Bariatric Commode	815067071852	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Drop-Arm Commode

The **ProBasics Bariatric Drop-Arm Commode** features a durable steel frame with an easy-release push button arm mechanism allowing for safe lateral patient transfers to and from the commode. Equipped with wider seat and reinforced steel frame to support patient weights up to 650 pounds, it comes complete with commode bucket and lid.

FEATURES

- Features a reinforced steel frame and an extra wide platform seat designed for comfort and stability
- Equipped with an easy-to-release arm mechanism allowing for safe lateral patient transfers to and from the commode
- Bucket easily slides into place
- Contoured armrests provide additional patient comfort
- Comes complete with a commode bucket and cover
- HCPCS Code: E0168*

SPECIFICATIONS

Width Between Arms:	24.5"
Outside Width:	28.5"
Seat Dimensions:	23"W x 18"D
Seat Opening:	8.25"W x 10"D
Adjustable Seat Height:	19.5" – 23.5"
Weight Capacity:	650 lb
Shipping Dimensions:	29.5" x 21" x 20"
Shipping Weight:	26 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
BSBDAC	ProBasics Bariatric Drop-Arm Commode	815067071869	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Commode Extra Wide Seat

ProBasics Bariatric Commode Extra Wide Seat features a 24-inch wide platform seat and reinforced steel frame to support patient weights up to 650 pounds and is ideal for larger patients. It comes complete with commode bucket and lid. The versatile design includes a commode bucket with cover, armrests for added comfort, angled legs for stability and non-marring rubber tips.

FEATURES

- Welded steel construction with 24-inch extra wide snap on platform seat designed for strength and stability
- Tool-free assembly
- Commode opening is contoured for increased comfort
- Seat height readily adjusts in one-inch increments
- Angled legs for stability with non-marring rubber tips
- Plastic armrests provide added comfort
- Comes complete with a commode bucket and cover
- Color: gray frame, gray accessories
- HCPCS Code: E0168

SPECIFICATIONS

Width b/w Arms:	26"
Outside Width:	27"
Seat Dimensions:	24"W x 18"D
Seat Opening:	9"W x 10.5"D
Adjustable Seat Height:	17" – 23"
Color:	Gray
Weight Capacity:	650 lb
Carton Dimensions:	32.5" x 22" x 26"
Shipping Weight:	41 lb
Warranty:	Limited lifetime on frame One year on parts

Item #	Description	UPC	UOM
BSB24C	ProBasics Bariatric Commode with Extra Wide Seat	815067071890	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Round Handle Cane

The **ProBasics Round Handle Cane** provides excellent support and is designed for patients who need additional assistance with balance. The round handle allows users to hook the cane over their arm while opening doors or performing other tasks. Supports patient weights up to 250 pounds.



Easy one-button height adjustment and anti-rattle ring



Round vinyl handle

FEATURES

- Easy one-button height adjustment with a locking anti-rattle ring allowing for quick, safe adjustment
- Adjusts in 1-inch increments
- Anodized aluminum finish with vinyl grip
- Curved handle allows cane to be retained on arm
- Easy-to-clean aluminum frame
- HCPCS* Code: E0100

SPECIFICATIONS

Handle Height Adjustment Range:	31" – 39"
Cane Diameter:	7/8"
Cane Tip:	3/4"
Weight Capacity:	250 lb
Product Weight:	13.44 oz
Carton Dimensions:	32" x 6.7" x 8.0"
Shipping Weight:	9.24 lb
Warranty:	Limited Lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNRSL	ProBasics Round Handle Cane	Silver	815067071265	10/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Folding Canes

ProBasics Folding Canes fold easily for storage or travel. The comfortable T-shape grip reduces hand stress and fatigue on the wrist for maximum comfort. Available in two attractive finishes with a limited lifetime warranty. Supports patient weights up to 300 pounds.



Folds compact, perfect for a purse or carry-on luggage



Contoured T-shape handle

FEATURES

- Folds in four convenient parts for easy storage
- Easy-to-use, one button height adjustment with a locking suppression ring allowing for safe fitting
- Comfortable, contoured T-shape handle
- Easy-to-clean aluminum frame
- Adjusts in 1-inch increments
- Attractive powder-coated finish available in black or bronze
- HCPCS* Code: E0100

SPECIFICATIONS

Adjustment Range:	32.5" – 36.5"
Cane Diameter:	7/8"
Cane Tip:	3/4"
Product Weight:	13 oz
Weight Capacity:	300 lb
Carton Dimensions:	20.5" x 6.3" x 12.6"
Shipping Weight:	16.6 lb
Warranty:	Limited lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNFBK	ProBasics Folding Cane	Black	815067071241	6/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Standard

Bariatric

Offset Canes

The **ProBasics Offset Cane** allows for easy height adjustment with a simple push of a button. The soft cushioned handle with convenient wrist strap provides comfort and grip support.



Easy one-button height adjustment



Soft cushioned handle with wrist strap

FEATURES

- Easy one-button height adjustment with a locking suppression ring allowing for quick, safe adjustment
- Bariatric model features reinforced tubing and a wider base for increased stability and support
- Adjusts in 1-inch increments
- Comfortable, contoured handle
- Easy-to-clean aluminum frame
- Attractive powder-coat finish
- Product covered under limited lifetime warranty
- HCPCS* Code: E0100

SPECIFICATIONS

Handle Height	Offset: 29" – 38"
Adjustment Range:	Bariatric Offset: 30" – 39"
Cane Diameter:	Offset: 7/8" Bariatric Offset: 7/8"
Cane Tip:	Offset: 3/4" Bariatric Offset: 3/4"
Weight Capacity:	Offset: 300 lb Bariatric Offset: 500 lb
Product Weight:	Offset: 12 oz Bariatric Offset: 2 lb 5oz
Carton Dimensions:	Offset: 33.1" x 7.9" x 7.3" Bariatric Offset: 30.7" x 11" x 9.1"
Shipping Weight:	Offset: 8.4 lb Bariatric Offset: 47.75 lb
Warranty:	Limited Lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNOFBK	ProBasics Offset Cane with Wrist Strap	Black	815067071135	10/cs
CNOFHD	ProBasics Bariatric Offset Cane with Wrist Strap	Black	815067071173	20/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Quad Canes

ProBasics Quad Canes feature a unique four-point design which allows for greater weight-bearing and stability than a single-point cane. The handle rotates allowing for a proper fit with right or left hand use. Small and large base models available in black with a limited lifetime warranty.

FEATURES

- Unique four-point design allows for greater weight-bearing and stability than a single-point cane
- Easy push button height adjustment in 1-inch increments with a locking, anti-rattle ring allowing for quick, safe adjustment
- Non-skid rubber tips on all four feet enhance durability for long term use
- Small base design allows for easier maneuverability while the large base design provides more stability
- Handle rotates for left or right hand use
- Quad canes stand up on their own when not in use
- Comfortable, contoured vinyl handle with sculpted grip reduces chance of slippage
- Lightweight aluminum upper shaft with welded steel base
- Tool-free left to right hand adjustments
- Bariatric model features a double-plated base for increased stability and support for patients up to 500 pounds
- HCPCS* Code: E0105

SPECIFICATIONS

Handle Height Adjustment Range:	Small & Large Base: 29"– 38" Bariatric: 27"– 36"
Cane Tip:	Small & Large Base: 3/4" Bariatric: 3/4"
Weight Capacity:	Adult: 300 lb Bariatric: 500 lb
Product Weight:	Small Base: 2.3 lb Large Base: 2.7 lb Bariatric Small Base: 3.2 lb Bariatric Large Base: 3.6 lb
Carton Dimensions:	Small Base: 25.6" x 10.6" x 8.7" Large Base: 25.6" x 13.8" x 11.8" Bariatric Small Base: 24.4" x 16.5" x 6.3" Bariatric Large Base: 24.2" x 12.2" x 8.9"
Shipping Weight:	Small Base: 14.6 lb Large Base: 17.8 lb Bariatric Small Base: 14 lb Bariatric Large Base: 8.5 lb
Warranty:	Limited Lifetime on frame 90 days on all parts

Item #	Description	Color	UPC	UOM
CNQSB	ProBasics Small Base Quad Cane	Black	815067071180	6/cs
CNQLB	ProBasics Large Base Quad Cane	Black	815067071197	6/cs
CNQHSBB	ProBasics Bariatric Small Base Quad Cane	Black	815067071227	4/cs
CNQHLBB	ProBasics Bariatric Large Base Quad Cane	Black	815067071234	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Aluminum Crutches

ProBasics Aluminum Crutches are sturdy, lightweight anodized aluminum crutches with push-button height adjustment. The extra thick underarm pads and grips offer additional comfort for use. The double extruded center tube is a solid one-piece design providing additional strength to weight-bearing area. Non-skid rubber tips provide excellent traction and absorb shock. Easy wing nut adjustments for hand grips. Available in youth, adult or tall sizes to accommodate all individuals. Supports patient weights up to 350 pounds.



Extra thick underarm pads



Easy push-button height adjustment

FEATURES

- Made of lightweight, anodized aluminum
- Reinforced one-piece center tube provides additional strength to weight-bearing area
- Push-button adjustment makes changing height easy
- Height adjustment in 1-inch increments
- Hand grip position readily adjusts with wing nuts
- Non-skid rubber tips provide excellent traction on many surfaces
- Extra thick underarm pads and grips offer additional comfort
- HCPCS Code: E0114*

SPECIFICATIONS

Patient Height Range:	Junior: 4'6" – 5'2" Adult: 5'2" – 5'10" Tall: 5'10" – 6'6"
Arm Pad to Hand Grip:	Junior: 9" – 13" Adult: 15" – 19" Tall: 16" – 20"
Arm Pad to Floor:	Junior: 37" – 46" Adult: 45" – 53" Tall: 53" – 61.5"
Weight Capacity:	350 lb
Product Weight:	Junior: 3.3 lb Adult: 3.6 lb Tall: 3.8 lb
Carton Dimensions:	Junior: 39" x 13" x 11" Adult: 47" x 13" x 11" Tall: 54" x 13" x 11"
Shipping Weight:	Junior: 29 lb Adult: 32 lb Tall: 34 lb
Warranty:	Limited lifetime on frame 90-day on all parts

Item #	Description	UPC	UOM
CRJ	ProBasics Aluminum Crutches, Junior	815067071289	8/cs
CRA	ProBasics Aluminum Crutches, Adult	815067071272	8/cs
CRAT	ProBasics Aluminum Crutches, Tall	815067071296	8/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Economy Two-Button Folding Steel Walker

ProBasics Economy Two-Button Folding Steel Walker provides maximum strength in a durable design, all backed by a limited lifetime warranty. Ample height adjustment options coupled with a two-button folding mechanism make it easy to use and accommodating to the most demanding patient needs. Available with or without 5" wheels and in junior or adult sizes. Supports patient weights up to 350 pounds.



Rear Glide Caps



Push Button Provides Audible Click When Locked



Height Adjustable

FEATURES

- 1-inch diameter steel frame construction provides superior strength
- Easy-to-use push button folding mechanism provides audible click when locked
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Front cross brace near top of walker allows for full stride
- Side bar adjusts with height to provide additional support
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage
- HCPCS* Code: E0135 (without wheels)
E0143 (with wheels)

SPECIFICATIONS

Width Inside Hand Grips:	17"
Width Inside Base Legs:	20"
Depth Opened at Base:	17"
Depth Folded:	4.5"
Wheel Size:	5"
Patient Height Range:	Adult: 5'2" – 6'4"
Weight Capacity:	350 lb
Product Weight:	WKSAN2B: 7 lb WKSAAW2B: 7.6 lb
Carton Dimensions:	WKSAN2B: 21.7" x 9.6" x 33" WKSAAW2B: 21.9" x 10.6" x 31"
Shipping Weight:	WKSAN2B: 32.6 lb WKSAAW2B: 36.2 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
WKSAN2B	ProBasics Economy Two-Button Steel Walker without Wheels, Adult	815067071050	4/cs
WKSAAW2B	ProBasics Economy Two-Button Steel Walker with 5" Wheels, Adult	815067070985	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Two-Button Folding Aluminum Walkers


ProBasics Two-Button Folding Aluminum Walkers provide maximum strength in a lightweight aluminum design, all backed by a limited lifetime warranty. Ample height adjustment options coupled with a two-button folding mechanism make ProBasics Folding Walkers easy to use and accommodating to the most demanding patient needs. Available with or without 5" wheels and in junior or adult sizes. Supports patient weights up to 300 pounds.



Rear Glide Caps Standard On Walker With Wheels



Push Button Provides Audible Click When Locked



Height Adjustable

FEATURES

- Attractive, lightweight aluminum frame
- Adjustable height settings
- Easy-to-use push button folding mechanism provides audible click when locked
- Front cross brace near top of walker allows for full stride
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage
- HCPCS* Code: E0135 (without wheels)
E0143 (with wheels)

SPECIFICATIONS

Width Inside Hand Grips:	17"
Width Inside Base Legs:	20"
Depth Opened at Base:	17"
Depth Folded:	4"
Patient Height Range:	Junior: 4'4" – 5'7" Adult: 5'2" – 6'4"
Weight Capacity:	300 lb
Product Weight:	WKAAN2B: 6.6 lb WKA AW2B: 7.5 lb WKA JW2B: 6.7 lb
Carton Dimensions:	WKAAN2B: 22" x 9.6" x 33" WKA AW2B: 22" x 9.6" x 31" WKA JW2B: 22" x 10.6" x 24"
Shipping Weight:	WKAAN2B: 26.6 lb WKA AW2B: 33 lb WKA JW2B: 29.2 lb
Warranty:	Limited Lifetime on frame One year on all parts

Item #	Description	UPC	UOM
WKAAN2B	ProBasics Two-Button Folding Aluminum Walkers, Adult	815067070954	4/cs
WKA AW2B	ProBasics Two-Button Folding Aluminum Walkers, Adult	815067070978	4/cs
WKA JW2B	ProBasics Two-Button Folding Aluminum Walkers, Junior	815067071036	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Sure Lever Release™ Folding Aluminum Walkers

The **ProBasics Sure Lever Release Folding Aluminum Walkers** are an ideal design for patients with manual dexterity limitations. The sure lever release folding mechanism allows users to close the walker without removing their hands from the walker handle. Ample height adjustment options coupled with a sure lever release folding mechanism make it easy to use and accommodating to the most demanding patient needs. Supports patient weights up to 350 pounds.



Sure Lever Release



Height Adjustable



Blue Flame Finish

FEATURES

- Attractive, lightweight aluminum frame
- Sure lever release folding mechanism allows users to close the walker without removing their hands from the walker handle
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Front cross brace near top of walker allows for full stride
- Adjustable height settings
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage
- HCPCS* Code: E0135 (without wheels)
E0143 (with wheels)

SPECIFICATIONS

Width Inside Hand Grips:	18"
Width Inside Base Legs:	20.5"
Depth Opened at Base:	17"
Depth Folded:	4.5"
Patient Height Range:	Junior: 4'4" – 5'7" Adult: 5'2" – 6'4"
Wheel Size:	5"
Weight Capacity:	350 lb
Product Weight:	WKAAWSL: 6.1 lb WKAAWSLB: 6.1 lb
Carton Dimensions:	WKAAWSL: 22.4" x 10.4" x 30.7" WKAAWSLB: 22.4" x 10.4" x 30.7"
Shipping Weight:	WKAAWSL: 27.4 lb WKAAWSLB: 27.4 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	UOM
WKAAWSL	ProBasics Sure Lever Release Folding Aluminum Walkers with 5" Wheels, Adult, Gray	815067071005	4/cs
WKAAWSLB	ProBasics Sure Lever Release Folding Aluminum Walkers with 5" Wheels, Adult, Blue Flame	815067071012	4/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Two-Button Folding Walkers

ProBasics Bariatric Two-Button Release Folding Walkers provide maximum strength and a rugged design, all backed by a limited lifetime warranty. Ample height adjustment options coupled with a two-button folding mechanism make ProBasics Folding Walkers easy to use and accommodating to the most demanding patient needs. Available with or without 5" wheels and in junior or adult sizes. Support patient weights up to 500 pounds.



Rear glide caps standard on walker with wheels



Push button provides audible click when locked



Height adjustable

FEATURES

- Attractive, lightweight aluminum frame
- Adjustable height settings
- Easy-to-use push button folding mechanism provides audible click when locked
- Front cross brace near top of walker allows for full stride
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl contoured hand grip for added comfort and long lasting wear
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces (come standard on models with 5" wheels)
- Folds flat for transport and storage
- HCPCS* Code: E0148 (without wheels)
E0149 (with wheels)

SPECIFICATIONS

Width Inside Hand Grips:	21.5"
Width Inside Base Legs:	32"
Depth Opened at Base:	20"
Depth Folded:	4.5"
Wheel Size:	5"
Patient Height Range:	Adult: 5'2" – 6'4"
Weight Capacity:	500 lb
Product Weight:	w/out wheels: 7.3 lb w/wheels: 10.6 lb
Carton Dimensions:	w/out wheels: 24" x 6" x 33" w/wheels: 26" x 6" x 31"
Shipping Weight:	w/out wheels: 16.7 lb w/wheels: 23.3 lb
Warranty:	Limited lifetime on frame One year on all parts


Item #	Description	UPC	UOM
WKABN2B	ProBasics Bariatric Two-Button Folding Walker without Wheels, Adult	815067071104	2/cs
WKABW2B	ProBasics Bariatric Two-Button Folding Walker with Wheels, Adult	815067071111	2/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Knocked Down Economy Steel Rollator w/ 6" Wheels

The **ProBasics Steel Rollator** is the perfect blend of style, comfort and affordability. Constructed of durable powder-coated steel, it offers great maneuverability and stability. It features easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Equipped with an easy folding mechanism for storage and transport. Supports patient weights up to 300 pounds. The new knock-down version features a quick tool-free first time assembly.



Water-Resistant
Nylon Pouch



Easy-to-Operate
Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Durable powder-coated steel frame
- Removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- Easy Tool Free setup - Installation guide included
- New and improved box size for increased product availability and reduced shipping costs
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy
Finish



Blue
Finish



SPECIFICATIONS

Overall Width:	23"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 9.5" (D)
Seat-to-Floor Height:	20.5"
Weight Capacity:	300 lb
Adjustable Handle Height:	30.5" – 36.5"
Wheel Size:	6"
Product Weight:	16.5 lb
Carton Dimensions:	23" x 4.75" x 22"
Shipping Weight:	19 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLS6BG-KD	ProBasics Economy Steel Rollator with 6-inch Wheels	Burgundy	815067079926	1/cs
RLS6BL-KD	ProBasics Economy Steel Rollator with 6-inch Wheels	Blue	815067079919	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Junior Aluminum Rollator w/ 6" Wheels

The **ProBasics Junior Rollator** is designed to fit individuals of shorter stature with a focus on style and affordability. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 250 pounds.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Lightweight aluminum frame with burgundy or blue finish
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy Finish



Blue Finish



SPECIFICATIONS

Overall Width:	22"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13" (W) x 13.5" (D)
Seat-to-Floor Height:	19"
Weight Capacity:	250 lb
Adjustable Handle Height:	29" – 34"
Wheel Size:	6"
Product Weight:	12 lb
Carton Dimensions:	22.6" x 7.3" x 31.3"
Shipping Weight:	14.3 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLAJ6BG	ProBasics Junior Rollator with 6-inch Wheels	Burgundy	815067071722	1/cs
RLAJ6BL	ProBasics Junior Rollator with 6-inch Wheels	Blue	815067071715	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Aluminum Rollator w/ 6" Wheels

The **ProBasics Aluminum Rollator w/6" Wheels** is designed to fit the majority of your patients' fundamental needs with a focus on style and affordability. Constructed of lightweight aluminum, it offers great maneuverability and stability. Equipped with easy-to-operate ergonomic hand brakes, a removable padded backrest and a water-resistant nylon pouch. Easily folds for storage and transport with built-in fold strap. Supports patient weights up to 300 pounds.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes



6-inch Wheels

FEATURES

- Lightweight aluminum frame in burgundy, blue and black
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 6" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

COLORS

Burgundy Finish  Blue Finish  Black Finish 

SPECIFICATIONS

Overall Width:	22.75"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13.75" (W) x 11.5" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	31" – 36"
Wheel Size:	6"
Product Weight:	13 lb
Carton Dimensions:	24.2" x 7" x 32.6"
Shipping Weight:	16 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
RLA6BG	ProBasics Standard Rollator with 6-inch Wheels	Burgundy	815067071630	1/cs
RLA6BL	ProBasics Standard Rollator with 6-inch Wheels	Blue	815067071647	1/cs
RLA6BK	ProBasics Standard Rollator with 6-inch Wheels	Black	815067071654	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Aluminum Rollator w/ 8" Wheels



The **ProBasics Aluminum Rollator** has large 8-inch wheels for improved maneuverability and is ideal for outdoor terrain. The padded seat and removable padded backrest provide a comfortable resting area. It features easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch, and a built-in fold strap for easy storage and transport. Supports patient weights up to 300 pounds.



Water-Resistant Nylon Pouch



Easy-to-Operate Ergonomic Hand Brakes

FEATURES

- Lightweight aluminum frame
- Two-button removable curved padded backrest
- Comfortable, flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0143
E0156 (seat)

SPECIFICATIONS

Overall Width:	22.75"
Width b/w Hand Grips:	18.25"
Seat Dimensions:	13.75" (W) x 11.5" (D)
Seat-to-Floor Height:	21"
Weight Capacity:	300 lb
Adjustable Handle Height:	31" – 36"
Wheel Size:	8"
Product Weight:	15.4 lb
Carton Dimensions:	24.8" x 9.1" x 35"
Shipping Weight:	18.3 lb
Warranty:	Limited lifetime on frame One year on all parts

COLORS

- Burgundy Finish  Blue Finish  Black Finish 

Item #	Description	Color	UPC	UOM
RLA8BG	ProBasics Rollator with 8-inch Wheels	Burgundy	815067071678	1/cs
RLA8BL	ProBasics Rollator with 8-inch Wheels	Blue	815067071661	1/cs
RLA8BK	ProBasics Rollator with 8-inch Wheels	Black	815067071685	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Bariatric Aluminum Rollator w/ 8" Wheels

The **ProBasics Bariatric Aluminum Rollator** is made with an extra-strength, reinforced aluminum frame and has large reinforced forks on 8-inch wheels for improved maneuverability and stability. The padded seat and removable padded backrest provide a comfortable resting area. It features easy-to-operate ergonomic hand brakes, a water-resistant nylon pouch, and a built-in fold strap for easy storage and transport. Supports patient weights up to 400 pounds



Extra Wide Seat and Water-Resistant Pouch



Easy-to-Operate Ergonomic Hand Brakes



Reinforced forks and 8-inch Wheels

FEATURES

- Lightweight aluminum frame
- Two-button removable curved padded backrest
- Comfortable, extra wide flip-up padded seat
- Easy-to-operate ergonomic hand brakes
- Height-adjustable for a custom fit
- Water-resistant nylon pouch attaches under seat for added privacy and security of personal items
- Equipped with a folding mechanism for easy storage and transport
- Non-marring 8" black casters with reinforced forks and soft grip tires, ideal for both indoor and outdoor use
- HCPCS* Code: E0149
E0156 (seat)

SPECIFICATIONS

Overall Width:	29"
Width b/w Hand Grips:	23.5"
Seat Dimensions:	19" (W) x 13" (D)
Seat-to-Floor Height:	22"
Weight Capacity:	400 lb
Adjustable Handle Height:	33.5" – 38.5"
Wheel Size:	8"
Product Weight:	18.7 lb
Carton Dimensions:	29.1" x 8.6" x 36.2"
Shipping Weight:	22.7 lb
Warranty:	Limited lifetime on frame One year on all parts

COLORS

Burgundy
Finish



Blue
Finish



Item #	Description	Color	UPC	UOM
RLAB8BG	ProBasics Bariatric Aluminum Rollator	Burgundy	815067071746	1/cs
RLAB8BL	ProBasics Bariatric Aluminum Rollator	Blue	815067071739	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.






Aspen 8" Height Adjustable Rollator

Experience comfort, stability, and convenience with the Aspen 8" Height Adjustable Rollator. Crafted from a lightweight yet durable aluminum frame, this rollator is built for years of worry-free use. With first-time tool-free assembly, setup is quick and hassle-free. Enjoy superior comfort with the extra-padded strap-style backrest, designed to support you during rest breaks. The rollator features five 1" handle height adjustments and four 1" seat height adjustments, making it easy to achieve the perfect ergonomic fit for a variety of users.

FEATURES

- Features an attractive, lightweight aluminum oval-tube frame that folds easily for storage and on the go travel – a perfect design for years of worry-free use
- 8" all terrain wheels are perfect for both indoor and outdoor use
- (4) 1" Seat height and (5) 1" handle height adjustments to custom fit a variety of users
- Extra padded flip-up seat offers excellent comfort and easy access to an under-seat storage bag
- New flip-up backrest that is flexible, wider, more comfortable, and height adjustable
- New braking system with covered calipers for easy adjustment
- Tool Free first-time assembly is fast and simple
- Weight capacity: 300 lbs
- Warranty: Limited lifetime on frame
- HCPCS* Code: E0143, E0156 (seat)

COLORS

Cobalt Blue  Brandywine  Charcoal 

SPECIFICATIONS

Weight Capacity:	300 lb
Overall Width:	23"
Seat Dimensions:	14.5" (W) x 14.5" (D)
Patient Height Range:	5'2" - 6'3"
Handle Height Range:	
Seat Height Setting	Handle Height Range
Seat-to-Floor Height: 20"	34" - 38"
Seat-to-Floor Height: 21"	35" - 39"
Seat-to-Floor Height: 22"	36" - 40"
Seat-to-Floor Height: 23"	37" - 41"
Wheel Size:	8"
Product Weight:	16.4 lbs
Carton Dimensions:	20" x 8" x 22"
Shipping Weight:	17 lbs
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Parts

* User body types vary, heights are approximate based on average person.

Item #	Description	UPC
RLAH08CB	Aspen Rollator, Cobalt Blue	815067012695
RLAH08BW	Aspen Rollator, Brandywine	815067012701
RLAH08BK	Aspen Rollator, Charcoal	815067012725

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Ascent Euro-Style Rollator

The Ascent is not only our newest model, it's the Flagship of our robust line of rollators. This lightweight, euro-style, 4-wheel rollator has both an adjustable height padded seat and backrest. The all new Ascent boasts a number of high-end features to ensure safety and longevity including dual sided, vented aluminum front forks, stylish 10" front and 8" rear casters with black non-marring tires and a unique top-to-bottom protected braking system. The Ascent is featured in three brilliant colors: Cobalt Blue, Brandywine, and Charcoal Black finish.



- Easy push button for handle height adjustments in 1" increments
- Aluminum brake cable conduit protects at the most frequent breakage point



- Innovative design hides brake lines for a clean look while protecting from wear and breakage
- Shielded adjustable braking mechanism offers excellent durability and performance



- Dual Sided Aluminum Front Forks with reinforced 10" Front & 8" rear casters for longevity and durability

FEATURES

- Lightweight Aluminum Frame for ease of use and longevity.
- Ergonomic handles with hidden brake cables; adjustable height in one-inch increments from 37.5" – 42.5" with an audible click as well as visual push button height settings.
- New and improved Dual Sided Aluminum Front Forks for increased stability.
- 18" x 10" thickly padded seat with a handle to quickly and easily fold and unfold side-to-side (rather than front-to back), allowing the unit to maneuver through narrow doorways.
- Dual Sided, vented aluminum front forks with 10-inch front and 8-inch rear heavy-duty wheels with black non-marring, soft-grip tires to traverse almost any terrain.
- Top to Bottom hidden braking system is fully adjustable, durable and stylish.

COLORS



SPECIFICATIONS

Weight Capacity:	300 lb
Overall Width:	23.5" (W)
Seat Dimensions:	18" (W) x 10" (D)
Seat Height:	Handle Height Range:
21"	37.5" - 40.5"
22"	38.5" - 41.5"
23"	39.5" - 42.5"
Patient Height Range:	5'6" - 6'5"
Adjustable Handle Height:	37.5" – 42.5"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	16.4 lb
Carton Dimensions:	20.7" x 8.85" x 26"
Shipping Weight:	17.5 lb
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Wear Parts

* User body types vary, heights are approximate based on average person.

Item #	Description	UPC	UOM
RLAS10CB	Ascent Rollator, Cobalt Blue	815067012732	each
RLAS10BW	Ascent Rollator, Brandywine	815067012640	each
RLAS10BK	Ascent Rollator, Textured Black	815067012756	each

* Map Pricing: \$259 Cannot be sold on Amazon

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



Ascent XR™ Euro-Style Rollator

The Ascent XR is a unique extra-rugged rollator that accommodates a larger range of users. It is featured in two colors: Cobalt Blue and Brandywine. The two major features of the Ascent XR are a larger weight capacity of 450 LBS and a rear shock absorbing suspension system.



- Easy push button for handle height adjustments in 1" increments



- Innovative design hides brake lines for a clean look while protecting from wear and breakage



- 10-inch Front Caster Reduces Resistance for Easy Navigation

- Aluminum brake cable conduit protects at the most frequent breakage point

- Shielded adjustable braking mechanism offers excellent durability and performance

FEATURES

- Lightweight, Aluminum frame with integrated suspension system is designed for durability and years of worry free use.
- 10-inch front and 8-inch rear heavy-duty wheels with black non-marring, soft-grip tires to traverse almost any terrain.
- Ergonomic handles – with top-down hidden braking system; Adjustable height in one-inch increments from 36" - 42.5" with audible "click" setting.
- Comes with a 24" x 13" thickly padded seat with a handle to fold and unfold side-to-side quickly and easily. This allows for the rollator to maneuver through narrow hallways and doors.
- Attractive front mesh tote attaches to the front of the frame for carrying personal belongings. It is easily attached and detached and can double as a shoulder bag when necessary.
- Folds up quickly and easily for transport or storage. Can easily fit in vehicles, storage closet, or public transport.

SPECIFICATIONS

Weight Capacity:	450 lb
Overall Width:	29" (W)
Seat Dimensions:	24" (W) x 13" (D)
Seat to Floor Height:	24"
Adjustable Handle Height:	36" – 42.5"
Wheel Size:	10" (front) 8" (rear)
Product Weight:	23 lb
Carton Dimensions:	29.7" x 11.2" x 37.6"
Shipping Weight:	30.4 lb
Limited Lifetime Warranty:	5 Year Warranty on Frame 1 Year on Wear Parts

COLORS

Cobalt Blue



Brandywine




Item #	Description	UPC	UOM
RLASXRCB	Ascent XR Rollator, Cobalt Blue	815067012657	each
RLASXRBW	Ascent XR Rollator, Brandywine	815067012688	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

Deluxe Knee Rollator

The **Deluxe Knee Rollator** offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Rollator is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required. The Knee Rollator is equipped with an easy folding mechanism that collapses using a thumb release lever and fits beneath the raised knee platform for ease of storage and transportation.



Adjustable Hand Brake

FEATURES

- Rollator can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Large 8-inch clear PU wheels with sealed bearings
- Adjustable hand brake ensures safety
- Padded knee platform with easy-adjustment mechanism
- Padded handlebar grips
- Limited lifetime warranty on frame, 90 days on all other parts

SPECIFICATIONS

Patient Height Range:	4'11" – 6'6"
Knee Platform Dimensions:	7" (W) x 14" (L)
Adjustable Handle Height:	31" – 41"
Knee Platform Height:	17.5" – 22"
Product Weight:	22.8 lbs.
Weight Capacity:	350 lbs.

COLORS

Black
Finish



White
Finish



Burgundy
Finish



Item #	Description	UPC	UOM
ROS-KSB	Deluxe Black Knee Rollator	092237601136	each
ROS-KSW	Deluxe White Knee Rollator	092237601471	each
ROS-KSBG	Deluxe Burgundy Knee Rollator	092237601143	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



K1 Wheelchair


STANDARD WHEELCHAIR



The **ProBasics K1 Manual Wheelchair** combines style and durability at an affordable price. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Supporting patient weights up to 250 pounds, the ProBasics K1 Standard Wheelchair is a popular choice among providers. Available with swing-away footrests or elevating legrests.



Flip-back armrests



Easy-to-operate wheel locks



Front riggings with heel loop

FEATURES

- Slide tube steel frame in black powder coat finish
- Single axle positioning
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- HCPCS* Code: K0001

K0195 OR E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16"
Seat-to-Floor Height ¹ :	19.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24 ½" 26 ½" 28 ½"
Overall Depth:	31 ¼" (without riggings)
Weight Capacity:	250 lb
Product Weight ² :	41– 43 lb
Carton Dimensions:	34" x 12" x 41"
Shipping Weight:	45 – 50 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC11616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070558	each
WC11816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070565	each
WC12016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070572	each
WC11616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070589	each
WC11816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070596	each
WC12016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070602	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



K1/K2 Wheelchair

STANDARD - HEMI HEIGHT



The **ProBasics K1/K2 Standard Wheelchair** is designed to meet a variety of patient profiles. With a patient weight capacity of 300 pounds, dual rear axle positions and multi-position caster forks, this wheelchair is designed at a standard hemi height making it one of the most versatile wheelchairs on the market today. Standard features include flip back removable padded armrests, easy-to-clean vinyl upholstery, chart pocket and heel loops on the front riggings. Available with swing-away footrests or elevating legrests.



Flip-back armrests



Easy-to-operate wheel locks



Front riggings with heel loop

FEATURES

- Slide tube, powder-coated silver vein steel frame
- Dual axle provides easy transition of seat height to hemi-level
- Flip-back, padded, upholstered armrests are removable
- Chart pocket on back for additional convenience
- Includes tool-free swing-away footrests with heel loops or elevating padded legrests with composite footplates
- Elevating legrests are equipped with a pneumatic piston for infinite elevating positions
- Easy push-to-lock wheel locks
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Composite handrims with 24" solid urethane tires
- 8" x 1" solid rubber front casters
- HCPCS* Code: K0002, K0001
K0195 OR E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16"
Seat-to-Floor Height ¹ :	19.5" Hemi 17.5"
Back Height:	18"
Back Style:	Fixed
Overall Width w/Arms:	24 ½" 26 ½" 28 ½"
Overall Depth:	31 ¼"
Weight Capacity:	300 lb
Product Weight ² :	41 – 43 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	45 – 50 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC21616DS	16" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070633	each
WC21816DS	18" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070640	each
WC22016DS	20" x 16"	Flip Back Desk-Length	Swing-Away Footrests	815067070657	each
WC21616DE	16" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070664	each
WC21816DE	18" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070671	each
WC22016DE	20" x 16"	Flip Back Desk-Length	Elevating Legrests	815067070688	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



K3 Lite Wheelchair

LIGHTWEIGHT WHEELCHAIR



The **ProBasics K3 Lite Wheelchair** is a standard lightweight, extremely versatile wheelchair.

The K3 Lite features a stylish, easy-to-clean nylon upholstery and comes standard with a variety of different features and benefits. The K3 Lite easily adapts to any patient's individual needs. It offers superior value and functionality for those users looking for a high-performance mobility product at an affordable price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



24" Rear Wheels

FEATURES

- High-Strength Steel Frame for years of worry free use
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14"
- 3-Position axles allow for removal of rear wheels for fast hemi-height adjustment
- Comfortable, lightweight and stylish easy-to-clean nylon upholstery
- 2 4" Rear Polyurethane tires and 8" composite casters for superior performance
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Includes tool-free adjustable front riggings
- Easy push-to-lock wheel locks (optional brake extensions included)
- HCPCS* Code: K0003

K0195 OR E0990 (elevating legrests only)
E0971 (anti-tippers)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16" – 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	19"
Back Style:	Fixed
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight ² :	35 – 38 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	49 - 55 lbs (Includes Legrests)
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC3L1616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012763	each
WC3L1816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012787	each
WC3L2016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067012800	each
WC3L1616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012770	each
WC3L1816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012794	each
WC3L2016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067012817	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



K4 Lite Wheelchair

HIGH STRENGTH LIGHTWEIGHT WHEELCHAIR

The **ProBasics K4 Lite Wheelchair** is a high-strength lightweight extremely versatile wheelchair. With an included seat and upholstery extension combined with all of its other unique features that come standard, the ProBasics K4 Lite easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance mobility product at an affordable price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Quick release rear axle

FEATURES

- 3-Position quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14"
- Adjustable height wrap-around back upholstery for maximum comfort
- Comfortable and stylish nylon upholstery
- Polyurethane tires and composite casters for superior performance
- Includes tool-free adjustable front riggings
- Adjustable anti-tippers included
- Easy push-to-lock wheel locks
- HCPCS* Code: K0195 or E0990 (elevating legrests only)
E0971 (anti-tippers)
E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16" – 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" – 17.5"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight ² :	32 – 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 – 39 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WCL41616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078547	each
WCL41816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078554	each
WCL42016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067078561	each
WCL41616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078516	each
WCL41816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078523	each
WCL42016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067078530	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



K4 Deluxe Wheelchair

HIGH STRENGTH K4 WHEELCHAIR



The **ProBasics K4 Deluxe Wheelchair** is a high-strength, feature rich wheelchair. This wheelchair comes standard with a 300 LB weight capacity, One -Button-Convertible Desk-To-Full Length Armrests, Quick Release Axle for Hemi-Height Adjustments. The K4 Deluxe easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance mobility product.



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Convertible armrests

FEATURES

- Convertible Armrests are both adjustable in height and length from Desk to Full in three settings
- Stylish rear Mag Wheels have brushed aluminum handrims and quick release rear axles for fast, too free hemi-height adjustment and for easy breakdown and transport
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- The chair comes standard with a wrap-around upholstery and an adjustable height and angle back (from 5 degrees to 20 degrees)
- Comfortable, slick nylon back upholstery is fluid resistant and easy-to-clean
- Adjustable anti-tippers included
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)
E0971 (anti-tippers)
E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20" 22"
Seat Depth:	16" - 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" - 18"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 LBS (16" / 18" / 20") 350 LBS (22")
Product Weight ² :	32 - 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 39 lb
Warranty:	Limited Lifetime on Frame, One Year on all Wear Parts

Item #	Seat	Arm Style	Leg Style	UOM
WC4DX1616CS	16" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX1816CS	18" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX2016CS	20" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX2216CS	22" x 16"	Height Adjustable, Flip Back, Convertible	Swing-Away Footrests	each
WC4DX1616CE	16" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX1816CE	18" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX2016CE	20" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each
WC4DX2216CE	22" x 16"	Height Adjustable, Flip Back, Convertible	Elevating Legrests	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.




Steel Transport Chair



The **ProBasics Steel Transport Chair** is a lightweight steel-construction design with a chip-proof and easy-to-maintain finish. The flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport. Comes complete with seat belt, tool-free swing-away composite footrests and heel loops. Fits patient weights up to 300 pounds.



Seat belt



Back-release hinge for easy folding



Folds flat for easy storage



8-inch wheels

FEATURES

- Lightweight steel construction (only 25 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish, silver vein carbon steel frame with permanent rear axle for strength and durability
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away composite footrests with heel loops included

SPECIFICATIONS

Overall Width:	21.25" (open) 9" (closed)
Width between Grips:	19"
Overall Height:	37"
Overall Depth:	38.75" with riggings 26.75" without riggings
Seat-to-Floor Height:	19"
Back Height:	16"
Seat Dimensions:	19" (W) x 16" (D)
Patient Weight Capacity:	300 lb
Product Weight:	25 lb
Carton Dimensions:	23" x 11" x 31"
Shipping Weight:	34 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	UPC	Seat Width	Colour	Wt. Cap.	UOM
TCS1916SV	ProBasics Steel Transport Chair	815067070480	19"	Silver Vein	300 lb	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Aluminum Transport Chair

The **ProBasics Aluminum Transport Chair** is a compact, lightweight design with an easy-to-maintain finish. The flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport. Comes complete with seat belt, tool-free swing-away composite footrests and heel loops. Supports patient weights up to 300 pounds.




Seat belt



Handy back-release hinge



Flared handles back to fold down flat



Push-to-lock rear wheel locks

FEATURES

- Lightweight aluminum construction (only 20 lbs)
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Stylish aluminum frame with permanent rear axle for strength and durability
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Push-to-lock rear wheel locks
- Composite 8" wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away composite footrests with heel loops included
- HCPCS Code: E1038*

SPECIFICATIONS

Overall Width:	22" open 9" closed
Width b/w Grips:	19.75
Overall Height:	37"
Overall Depth:	38.75" with riggings 26.75" without riggings
Seat-to-Floor Height:	19"
Back Height:	16"
Seat Dimensions:	19" (W) x 16" (D)
Weight Capacity:	300 lb
Product Weight:	20 lb
Carton Dimensions:	23" x 11" x 31"
Shipping Weight:	26 lb
Warranty:	Limited Lifetime on frame One year on all parts

COLORS



Item #	Description	Color	UPC	UOM
TCA1916BL	ProBasics Aluminum Transport Chair	Blue	815067070503	1/cs
TCA1916BG	ProBasics Aluminum Transport Chair	Burgundy	815067070497	1/cs
TCA1916PK	ProBasics Aluminum Transport Chair	Pink	815067070510	1/cs
TCA1916BK	ProBasics Aluminum Transport Chair	Black	815067072682	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Aluminum Transport with 12-inch Rear Wheels

The **ProBasics Aluminum Transport Chair with 12" rear wheels** has a lightweight aluminum frame and padded upholstery. The flared handles and handy back-release hinge allow the back to fold down for easy storage and transport. The companion activated hand brakes are ideal for locking the rear wheels during patient transfers, and provide for added safety and control. Comes complete with seat belt, swingaway footrests with composite footplates and heel loops. Supports patient weights up to 300 lb.



Flared handles for easy folding



Back-release hinge for easy folding



Folds flat for easy storage



12" rear wheels

FEATURES

- Lightweight aluminum construction (only 20 lb)
- 12" rear wheels for better performance on uneven surfaces
- Companion activated hand brakes are ideal for locking the rear wheels during transfers
- Flared handles and handy back-release hinge allow the back to fold down flat and compact for easy storage and transport
- Powder coated aluminum frame is lightweight and durable
- Padded full-length permanent armrests
- Padded nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Composite wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away composite footrests with heel loops included
- HCPCS Code: E1038*

SPECIFICATIONS

Overall Width:	21.5" (open) 10.25" (closed)
Width between Grips:	19.25"
Overall Height:	37"
Overall Depth:	39" with riggings 27.5" without riggings
Seat-to-Floor Height:	19"
Back Height:	16"
Seat Dimensions:	19" (W) x 16" (D)
Patient Weight Capacity:	300 lb
Product Weight:	20 lb
Carton Dimensions:	25" x 11" x 31"
Shipping Weight:	29.5 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
TCA191612BK	ProBasics Aluminum Transport with 12-inch Rear Wheels	Black	815067070527	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.




Bariatric Steel Transport Chair

The **ProBasics Bariatric Steel Transport Chair** with 12-inch rear wheels has dual reinforced cross braces and durable nylon upholstery. The handy back-release hinge allows the back to fold down for easy storage and transport. The companion activated hand brakes are ideal for locking the rear wheels during patient transfers, and provide for added safety and control. Comes complete with seat belt, swingaway footrests with aluminum footplates and heel loops. Supports patient weights up to 450 lb.



Companion activated hand brakes



Back-release hinge for easy folding



Back folds down quickly for easy storage



Solid rubber tires are lightweight & maintenance free

FEATURES

- Heavy-duty reinforced steel construction
- Stylish, silver vein steel frame with dual, reinforced cross braces for strength and durability
- 12" rear wheels for better performance on uneven surfaces
- Companion activated hand brakes are ideal for locking the rear wheels during transfers
- Handy back-release hinge allows the back to fold down for easy storage and transport
- Removable desk-length padded armrests
- Reinforced nylon upholstery is comfortable, lightweight and easy to clean
- Seat belt for added safety
- Nylon seat straps allow for ease of folding and lifting
- Composite wheels with solid rubber tires are lightweight and maintenance free
- Tool-free, swing-away aluminum footrests with heel loops come standard
- HCPCS Code: E1039*

SPECIFICATIONS

Overall Width:	27" (open) 11" (closed)
Width between Hand Grips:	20"
Overall Height:	37.75"
Overall Depth:	41.5" with riggings 31.5" without riggings
Seat-to-Floor Height:	20"
Back Height:	18"
Seat Dimensions:	22" (W) x 18" (D)
Patient Weight Capacity:	450 lb
Product Weight:	41 lb
Carton Dimensions:	30" x 12" x 33"
Shipping Weight:	66 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Description	Color	UPC	UOM
TCS221612SV	ProBasics Bariatric Transport with 12-inch Rear Wheels	Silver-Vein	8150670753	1/cs

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



K4 Transformer

TRANSFORMS FROM K4 HIGH STRENGTH
LIGHTWEIGHT WHEELCHAIR TO A TRANSPORT CHAIR

The perfect upgrade from a standard K0004 wheelchair, the **ProBasics TRANSFORMER™** is a dual purpose, high strength lightweight wheelchair that quickly transforms from a self-propelled wheelchair to a compact transport chair with 8" wheels all around. Great as a transitional chair for accompanied excursions while retaining self-propel capabilities. Also allows the slimming down of the chair to navigate narrow doorways.



24-inch quick release rear wheels



Transforms into a Transport Chair



Flip back, desk length, padded arms

FEATURES

- Lightweight aluminum frame
- 8-inch full swivel front casters mounted on aluminum forks
- 24-inch "quick release" rear wheels with composite handrims
- Bolted on 8-inch transport chair rear wheels
- Self-locking paddle release, flip-back padded desk length arms
- Lightweight black, woven nylon upholstery
- Fold down back with easy release paddles
- Complete with lightweight aluminum wheel locks for both 24-inch and 8-inch rear wheel modes
- Detachable flip-up rear anti-tippers
- Accessory pouch sewn into back
- Footrests include heel loops and composite footplates (optional elevating legrests sold separately - WCT4ELR)
- Folds for easy storage
- HCPCS* Code: K0004
K0195 or E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16"
Seat-to-Floor Height ¹ :	Adult 20" Hemi 17"
Back Height:	16"
Back Style:	Folds downs
Overall Width w/Arms:	18" – 29"
Overall Depth:	28" – 32" (without riggings)
Weight Capacity:	250 lb
Product Weight ² :	32 – 34 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 – 37 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WCT41616DS	16" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070817	each
WCT41816DS	18" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070824	each
WCT42016DS	20" x 16"	Flip Back, Desk-Length	Swing-Away Footrests	815067070831	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.



K7 Bariatric Wheelchair

HEAVY DUTY WHEELCHAIR

The **ProBasics K7 Heavy Duty Wheelchair** features a robust design to accommodate patients with wider seat and higher weight capacity requirements. With seat widths from 24 to 28 inches, a weight capacity of 450 - 600 pounds and a heavy-duty slide tube frame, the ProBasics K7 offers a rugged yet attractive option for the bariatric patient population.



Removable, desk length, padded arms



Heavy duty vinyl upholstery



Aluminum footplate

FEATURES

- Heavy Duty, slide tube silver vein steel frame
- 8" x 2" semi-pneumatic casters are adjustable in 3 positions
- Dual axle provides easy transition of seat height to hemi-level
- 24" rear urethane wheels with durable hand rims
- Removable desk length arms with padded armrests
- Heavy duty vinyl padded upholstery with chart pocket
- Tool-free footrests include heel loops and Aluminum footplates
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite elevating positions and aluminum footplates
- Easy push-to-lock wheel locks
- HCPCS* Code: K0007
K00195 or E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	22" 24" 26" 28"
Seat Depth:	18" 20"
Seat-to-Floor Height ¹ :	Adult 21" Hemi 19"
Back Height:	16" 18"
Back Style:	Fixed
Overall Width w/Arms:	30" 32" 34½" 37"
Overall Depth:	33" 34" (without riggings)
Weight Capacity:	450 – 600 lb
Product Weight ² :	51– 64 lb
Carton Dimensions:	38" x 13" x 35"
Shipping Weight:	58 – 70 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WC72218DS	22" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070855	each
WC72418DS	24" x 18"	Removable, Desk-Length	Swing-Away Footrests	450 lb	815067070862	each
WC72620DS	26" x 20"	Removable, Desk-Length	Swing-Away Footrests	500 lb	815067072057	each
WC72820DS	28" x 20"	Removable, Desk-Length	Swing-Away Footrests	600 lb	815067072064	each
WC72218DE	22" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070879	each
WC72418DE	24" x 18"	Removable, Desk-Length	Elevating Legrests	450 lb	815067070886	each
WC72620DE	26" x 20"	Removable, Desk-Length	Elevating Legrests	500 lb	815067072071	each
WC72820DE	28" x 20"	Removable, Desk-Length	Elevating Legrests	600 lb	815067072088	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.




Reclining Wheelchair



The **ProBasics Reclining Manual Wheelchair** provides full back and headrest support as patient reclines from 90 to 160 degrees. The 13-inch headrest and padded, full-length arms are removable. Black, padded vinyl upholstery resists mildew and bacteria. The frame has a chip-resistant silver vein finish. Mag-style wheels are lightweight and set back on frame to prevent patient from tipping when in the reclined position.



Reclines from 90° – 160°



Caregiver activated recline mechanism



13-inch headrest supports patient when reclined and is removable

FEATURES

- Slide tube steel frame with an attractive silver vein finish that is chip resistant and easy to maintain
- Padded Desk-length, detachable armrests provide added patient comfort
- Chart pocket on backrest provides additional convenience
- Tool-free elevating padded legrests are equipped with a pneumatic piston for infinite positions and composite footplates
- Composite handrims with mag-style wheels are lightweight and maintenance free
- Wheels are set back on frame to prevent tipping
- Rear anti-tippers come standard
- Headrest extension comes standard with head immobilizer
- 8" x 1" solid rubber front casters
- Easy push-to-lock wheel locks
- HCPCS* Code: K0001/K0006

E1226 (reclining back)

E0966 (headrest extension)

E0990 (elevating legrests only)

SPECIFICATIONS

Seat Width:	16" 18" 20" 22"
Seat Depth:	17"
Seat-to-Floor Height ¹ :	Adult 20"
Back Height:	19.5" 34"
Back Style:	Reclining
Overall Width w/Arms:	24½" 26½" 28½"
Overall Depth:	31¼" (without riggings)
Weight Capacity:	300 lb 400 lb (WCR2218E)
Product Weight ² :	51 – 60 lb
Carton Dimensions:	40" x 12" x 36"
Shipping Weight:	71 – 75 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	Weight Limit	UPC	UOM
WCR1616E	16" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070893	Each
WCR1816E	18" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070909	Each
WCR2018E	20" x 17"	Removable, Desk-Length	Elevating Legrests	300 lb	815067070916	Each
WCR2218E	22" x 17"	Removable, Desk-Length	Elevating Legrests	400 lb	815067070923	Each


* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.

Meridian General Use Foam Cushion

HCPCS
E2601



- Economical seating surface provides comfort and pressure reduction
- Designed to assist in the prevention of pressure ulcers
- Cushion contours the seating area
- Removable, waterproof nylon top with non-slip vinyl bottom
- Color: navy
- Weight bearing: See chart below
- One-year limited warranty
- TAA Compliant options available, ask your Compass Health Representative



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
WCF-163	16" x 16" x 3"	275 lbs.	each	E2601
WCF-183	18" x 16" x 3"	300 lbs.	each	E2601
WCF-203	20" x 16" x 3"	325 lbs.	each	E2601
WCF-2083	20" x 18" x 3"	350 lbs.	each	E2601
WCBF-223	22" x 16" x 3"	400 lbs.	each	E2602
WCBF-2283	22" x 18" x 3"	425 lbs.	each	E2602
WCBF-2483	24" x 18" x 3"	450 lbs.	each	E2602
GEL-182	18" x 16" x 2"		each	N/A

HCPCS
E2611




Meridian General Use Back Cushion

- Provides upper back and lumbar support
- Designed to reduce lower back pain and muscle tension
- Removable, waterproof nylon cover
- Stabilization board provides enhanced rigidity and support
- Non-skid backing
- Color: Black
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
BKF-1617	16" x 17"	NA	each	E2611
BKF-1817	18" x 17"	NA	each	E2611
BKF-2019	20" x 19"	NA	each	E2611
BKFB-2219	22" x 19"	NA	each	E2612

HCPCS
E2603



Meridian Optimum Comfort Gel Seat Cushion



- Unique channeled gel bladder prevents gel migration in high-pressure areas
- Gel channels provide relief over the entire cushion area
- Foam base provides additional support
- Inner foam core has a liner to protect foam's integrity
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-9B 16X16	16" x 16"	275 lbs.	each	E2603
HDC-9B 18X16	18" x 16"	300 lbs.	each	E2603
HDC-9B 18X18	18" x 18"	325 lbs.	each	E2603
HDC-9B 20X16	20" x 16"	350 lbs.	each	E2603
HDC-9B 20X18	20" x 18"	400 lbs.	each	E2603
HDC-9B 22X18	22" x 18"	425 lbs.	each	E2604
HDC-9B 24X18	24" x 18"	450 lbs.	each	E2604

HCPCS
E2607



Meridian Optimum Pommel Cushion



- Contoured base and pommel
- Four grades of foam
- Inner core has a liner to protect foam
- Cover is a 4 way stretch polyester/polyurethane
- Moisture repellent & low shear/friction
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-POM1616	16" x 16" x 3"	275 lbs.	each	E2607
HDC-POM1816	18" x 16" x 3"	275 lbs.	each	E2607
HDC-POM2016	20" x 16" x 3"	275 lbs.	each	E2607

HCPCS
E0190



Meridian Optimum Coccyx Cushion



- Unique design for relief of pressure on coccyx area
- Weight bearing: See chart below
- One-year limited warranty
- Proudly made in the USA



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-CX1616	16" x 16" x 3"	225 lbs.	each	E0190
HDC-CX1816	18" x 16" x 3"	225 lbs.	each	E0190
HDC-CX2016	20" x 16" x 3"	225 lbs.	each	E0190

Ordering & Shipping Information

ORDERING POLICY

A valid Compass Health account number is required when placing an order. Once an order for stock product has reached booked status (ready to ship) and is in the shipping process, changes to the order will not be accepted. If it is necessary to change or cancel an order for a non-stock item, Compass Health must be notified immediately, within the same business day. All orders placed with Compass Health must reach a minimum level of \$50.00 in product (before taxes and shipping charges). Credit cards are accepted: Visa, MasterCard and American Express.

FREIGHT POLICY

Compass Health ships all orders received by 3:00 p.m. EST same day. Stock Orders must have a booked (Ready-to-Ship) status. Orders placed on hold for credit reasons or price deviations will be delayed for shipment until all holds are cleared. All shipments from Compass Health are F. O. B. Origin. This means that the title of product transfers to the buyer at the time the product leaves the Compass Health shipping dock or manufacturing location.

Exceptions to Freight Policy:

- Published freight costs plus handling fees will be charged on shipments to Puerto Rico and Canada.
- Any order requesting expedited delivery will be charged at published list rates. Any order requiring special handling may be subject to additional charges.
- COD shipments will be subject to additional charges.
- Direct ship and clearance items will be charged at published list freight costs.

A customer may pick up any available stock product at a Roscoe Distribution Center with a minimum of four (4) hours advance notice, from 9:00 a.m. - 4:00 p.m., local distribution center time. If a product is not available in a customer's primary shipping location, Compass Health reserves the right to ship a product from its choice of location.

A customer may obtain tracking and delivery confirmation from the carrier, or from a member of our Customer Service Team at (800) 376-7263 or at www.compasshealthbrands.com.

WARRANTY

We warrant each item to be free of defects in material and workmanship under normal use. This warranty is in lieu of all other warranties expressed, implied or previously printed and does not apply in the event of negligence, abuse, misuse or normal wear and tear. We are not responsible for injury that may result from abuse, misuse or user negligence and our liability is limited to original selling dealer or distributor. Unauthorized service or modification of product invalidates warranty. If a covered defect is found, Compass Health at its option will do one of the following: 1) Replace any defective or nonconforming part or product. 2) Issue a credit or refund at the original selling price.

MERCHANDISE INSPECTION

The customer is required to inspect all merchandise upon delivery. Any visible damage or carton shortage should be noted on the carriers delivery receipt. Claims for content shortage should be made to Compass Health Customer Service Department within five (5) days of delivery or any such claim is waived. Concealed damage requires immediate inspection by the carrier or any claim is waived. The customer is responsible for filing of all freight and concealed damage claims with the carrier.

REFUSED SHIPMENTS/RECONSIGNMENT

It is the responsibility of the customer to contact Compass Health within two (2) business days of the attempted delivery to inform Roscoe that a shipment has been refused. Refusing a shipment in full or part does not guarantee a credit to the consignee. If the consignee refuses the shipment, the carrier will return the product back to the original shipping location.

The refused shipment will undergo one or all of the following:

- RA assignment in accordance to RA assignment policy
- Returned to Consignee when applicable

Allow for thirty (30) days to process any RA associated with a refused shipment order. If a reconsignment of shipment is necessary, and it is found the fault of Compass Health, then Roscoe will pay reconsignment charges and have the freight delivered to the original specified location. If a reconsignment is necessary, and it is found to be the fault of the consignee, then the consignee must accept and receive delivery of the shipment. Once the shipment is received, the consignee can negotiate a re-delivery of the product at the consignee's expense.

RETURN GOODS POLICY

Compass Health will accept returns within 30 days from the original ship date. 25% restocking fee will be subtracted from the original cost and a credit for the difference will be issued. Returns must be shipped with freight charges prepaid. Shipping costs are not refundable.

All products must be returned in original packaging and must contain all components, accessories and user manuals. In the event that any components are missing you are responsible for the cost of the replacement component and the 25% restocking fee. All returns must be approved with a Return Authorization Number (RA). Please call a member of our customer service team at (800) 376-7263 to obtain a RA number.

Provide the following information when calling:

- Item number
- Original order number
- Product serial number/lot number
- Reason for return

The RA number must be clearly marked on the returned carton and is valid for 10 business days from the date of issue. Returned merchandise must be in the same unit of measure as originally purchased.

Return labels or call tags can be issued by our customer service department to return merchandise. Associated fees and return freight charges will apply. All returns of drop-shipped items are subject to a restocking fee as well as inbound and outbound freight charges.

Returns will not be accepted on items that are:

- Missing their serial number
- Specials order items
- Returned more than 30 days after delivery
- Returned without notification

Freight Overages, Shortages and Damaged Merchandise: All freight overages, shortages and/or damages must be reported immediately upon receipt of shipment. Visible damage should be noted on the freight company's receipt at the time of delivery. All packing and shipping documentation must be retained for damage claims. Contact a member of our customer service team at (800) 376-7263 for final disposition.



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

www.compasshealthbrands.com

rev_05_2025