


15303 BARIATRIC BED 54" (M SERIAL #)



1. Head Board - 15355H
Head & Foot Board (1 of ea)- 15355
2. Head Frame- 15315
3. Foot Frame- 15325
4. Foot Board- 15355F
5. Head Board Panel- 15303P1005
Foot Board Panel- 15303P1006
6. Head Board Frame- 15303P1051
Foot Board Frame- 15303P1052
Panel Frame Screw- 15005PS1051
7. Head Motor- 15302MH
8. Head Motor Actuating Tube- 15300PH1014
- 9a. Hand Pendant Junction Box- 15302P1000
- 9b. Hand Pendant Ext. Cord – Head- [3]
- 9c. Hand Pendant Ext. Cord – Foot- [3]
10. Head or Foot Board Gear Box- 15005P1028
11. Emergency Crank Head/Foot- 15005P1023
12. Hi/Lo Motor- 15300HL
13. Adjustable Leg Head 565mm- 15300P1015L [1]
Adjustable Leg Foot 505mm- 15300PF1015L [1]
14. Caster Receptacle- 15300EC
15. Caster w/Lock- 15300CL
Caster w/o Lock- 15300C
16. Foot Motor- 15300MF
17. Q/R Spring Knob- 15005P-1028
18. 3 Port Motor Control Box Screw Type- 15005JB [2]
3 Port Motor Control Box Q/R Type- 15005JBN [2]
19. Hand Pendant Screw On- 15005HC
Hand Pendant Quick Release- 15005HCN
20. Power Cord- 15005PC
21. Foot Motor Actuating Tube- 15300P1014
22. Attachment Pin & Hardware- 15300P1016
23. Frame Support Rail- 15300FSR
24. Hi/Lo Motor- 15300HL
25. Head Motor Tube- 15005P1025
26. Foot Motor Tube- 15005PF1025
27. Emergency Crank Hi/Lo- 15544C

[1] Includes cable & plastic guides, does not include part 14.

[2] When stock is gone use part number 15005JBU.

[3] No Longer Available