


15003PN



- 1. Head Frame - 15545-01
- 2. Foot Frame - 15545-02
- 3. Hi/Low Shaft - 15545-03
- 4. Frame Fastener Clip & Strap- 15560-12
- 5. Head Board - 15546-01
- 6. Foot Board - 15546-02
- 7. Hand Crank, 15571, 1/ea- 15561-10

- 8. Locking Caster- 15044CL
- 9. Non-Locking Caster- 15044CNL
- 10. Manual Crank Head- 15555-01
- 11. Manual Crank Foot- 15555-02
- 12. Head Tube Extension- 15555-03
- 13. Attaching Pin- 15005P-1016
- 14. Crank Mounting Knobs, (Pair)- 15555-05

09/23/2021