


15005LPN



- 1. Head Frame- 15545-01
- 2. Foot Frame- 15545-02
- 3. Hi/Low Shaft - 15545-03
- 4. Frame Fastener Clip & Strap- 15560-12
- 5. Head Board- 15546-01
- 6. Foot Board- 15546-02
- 7. Hand Crank- 15561-10
- 8. Caster with Lock- 15044CL
- 9. Caster w/o Lock- 15044CNL
- 10. Head Actuating Bar- 15552-01
- 11. Foot Actuating Bar- 15552-02
- 12. Motor Guard- 15552-03

- 13. Motor Guard Mntg Knobs (Pair)- 15552-04
- 14. Emergency Motor Crank- 15552-05
- 15. Cotter Pin- 15005P1016
- 16. Motor Q/R Knob- 15005P-1028
- 17. Head Motor- 15005MH-D
- 18. Foot Motor- 15005MF-D
- 19. Control Box- 15005JBN-D
- 20. Hand Pendant- 15005HC-D
- 21. HiLo Motor- 15005HL-D
- 22. HiLo Motor Mount- 15033MBC
- 23. Hi/Lo Motor Bracket Knob (ea)- 15033MBC-KNOB

11/08/2021