


15235 ELECTRIC HOMECARE BED DELTA LOW



- | | |
|--|--|
| 1. Head Spring- 15541 | 19. Motor Mount Cover- 15030-7RMM [4]
Motor Mount Cover- 15030P-1006 [5] |
| 2. Foot Spring- 15549 | 20. Motor Cover- 15030MC |
| 3. Head Board- 15547H [1] | 21. Motor Strap (2 ea)- 15030-MS |
| 4. Head Board Panel- 15547H-HP-1 | 22. Caster w/Lock- 15044CL
Caster w/o Lock- 15044CNL
Caster Set (2 of ea)- 15044CSET |
| 5. Head Board Panel Frame- 15547H-PF | 23. Caster Receptacle- 15030P-1001 |
| 6. Head Board Gear Box- 15544HGB | 24. Frame Linkage- 15030-FRLINK |
| 7. Foot Board- 15547F [1] | 25. Motor Coupler- Not Available [6]
Motor Coupler- 15033MHLN-C-N [7] |
| 8. Foot Board Panel- 15547F-FP-1 | 26a. Motor Bracket- 15033MBN [6] |
| 9. Foot Board Panel Frame- 15547F-PF | 26b. Motor Bracket- 15033MBN-DELTA [7] |
| 10. Foot Board Gear Box- 15544FGB [2] | 26c. Motor Bracket Mount- 15544GBMB |
| 11. Emergency Crank- 15547-CRANK | 27. Extension Spring Head or Foot- 15005-SPRING
Extension Spring Side- 15005-SPRING-1 |
| 12. Hi/Lo Motor/Gear Box- 15033MHLN-MGBA [3]
Gear Box Converter- 15033MHLN-GB | 28. Spring Deck Bolt- 15030-FRLINK [8] |
| 13. Hi/Lo Shaft w/Collar- [9] | 29. Sleep Surface Hook- 15005P-1040 |
| 14. Hi/Lo Shaft Collar- 1503C | |
| 15. Head/Foot Motor- 15037MN | |
| 16. Hand Pendant- 15033HC-G | |
| 17. Head Actuating Bar- 15030PN-1035 | |
| 18. Foot Actuating Bar- 15030PN-1036 | |

- [1] Head & Foot Board - 15547H-F
- [2] Without manual crank.
- [3] Before s/n 1S120329900 must order parts 27b. & 27c.
- [4] For 7R or 9R serial number.
- [5] Pre Delta with Okin-Ilcon motor.
- [6] Before serial number 1S120329900
- [7] After serial number 1S120329900
- [8] 2 bolts included with this part.
- [9] Discontinued. 15005S is sub. Does not Include collar.

[Diagrams](#)

[Home](#)

15235_rev. 08_03_2021