


**15005L LOW FULL ELECTRIC HOMECARE BED**

- |                                             |                                                 |
|---------------------------------------------|-------------------------------------------------|
| 1. Head Spring- 15008                       | 13. Hand Pendant- 15005HCN-L                    |
| 2. Foot Spring- 15013-L                     | 14. Power Cord- 15005PC-L                       |
| 3. Foot Board- 15547F                       | 15. Foot Motor- 15005MF-L                       |
| 4. Head Board- 15547H                       | 15a. Foot Motor Tube- 15005PF1025-L             |
| Head & Foot Board (1 ea)- 15547H-F          | 16. Sleep Surface Hook- 15005P-1040             |
| 5. Head Board Panel- 15547H-HP-1            | 17. Hi/Lo Shaft- 15005S                         |
| Foot Board Panel- 15547F-FP-1               | 18. Hi/Lo Shaft Adjustment Button- 15005SPN     |
| 6. Foot Board Panel Frame- 15547F-PF        | 19. Head Board Gear Box- 15544HGB               |
| Head Board Panel Frame- 15447H-PF           | 20. Extension Spring Head of Foot- 15005-SPRING |
| 6a. Panel Frame Screw- 15005PS1051          | Extension Spring Side- 15005-SPRING-1           |
| 7. Head Motor- 15005MH-L                    | 21. Frame Linkage- 15030-FRLINK                 |
| 7a. Head Motor Tube- 15005P1025-L           | 22. Head or Foot Motor Cable- 15005HFC-L        |
| 8. Foot Gear Box- 15544FGB                  | 23. Q/R Spring Knob- 15005P-1028                |
| 9. Emergency Hi/Lo Crank- 15544C            | 24a. Hi/Lo Motor- 15005HL-L                     |
| 10. Caster Receptacle- 15030P-1001          | 24b. Hi/Lo Motor Mount Bracket- 15005P1029-L    |
| 11. Caster w/Lock- 15044CL                  | 24c. Hi/Lo Motor Cable- 15005HLC-L              |
| Caster w/o Lock- 15044CNL                   | 24d. Hi/Lo Motor Coupler- 15005HLS              |
| Caster Set (2 of ea)- 15044CSET             | 25. Motor Attachment Pin- 15005P1016            |
| 12. Motor Control Box- 15005JBN-L           | 26. Spring Deck Bolt- 15030-FRLINK *            |
| Motor Control Box Mtg. Bracket- 15005JBMB-L | 27. Emergency Head/Foot Crank- 15005EC-L        |

\* 2 bolts included with this part.