

NRS185007



1. Seat with Insert- NRS185006-01

2. Seat Clip- NRS185006-02

3. Back- NRS185006-03

4. Arm Pads (pr) - NRS185006-04

5. Left Armrest Bracket- NRS185006-05

6. Right Armrest Bracket- NRS185006-06

7. Armrest Tube- [1]

8. Caster with Leg- NRS185006-08

9. Frame Clips (4 ea) - NRS185006-09

10. Swing Away Footrest- NRS185006-11

11. Footplate (pr) - NRS185006-12

12. Angle Adj. Footplate Mnt @ Hdwr- NRS185006-13

13. Down Tube & Plastic Insert- NRS185006-14

14. Footrest Height Adj. Knob & Bolt- NRS185006-15

15. Footrest Pivot Saddle & Hdwr. - NRS185006-16

16. Footrest Frame Latch- NRS185006-17

17. Back Locking Knob (pr)- NRS185006-21

Substitute (ea)- TRID-2A44

18. Ant-Rattle Insert- NRS185006-22

19. Frame Plugs (6 ea) – [1]

20. Pail- 11106-1

21. Seat Belt- STDS851-SA [2]

[1] Stock is gone the part is no longer available.

[2] When stock is gone the part is no longer available,
Footrest have an adjustment range from 13.5" to 18.5"

01/18/2023