

11130



1. Back	11130B
2. Armrest	11130A
3. Leg	11130L
4. Seat	11130S
5. Push Button	12002PB
6. Tips	RTL10389GB
7. Commode Pail	11120B
8. Pail Bucket	
3 Screws on ea side	11135BRK **
4 Screws on ea side	11135BRKN *
* Serial # 1A1104000001 or later	
Hardware Kit	11130HDW

** When Stock is gone, the part is no longer available.

Commode Hole- 8.5" X 10"

08/29/18

[Diagrams](#)

[Home](#)