




SUNFIRE PLUS EC



1. Headrest 18" [3]
Before Ser # 2W301G0610000042
Before Ser # (W076306)
Headrest 18"- SEC2010N
After Ser # 2W301G0610000042
After Ser # (W076306)
Headrest 20"- SF2010 [2]
2. Seat Assembly w/o Arms
18" SEC9008-18
20" [1]
3. Armrest Left- SF9007L
Armrest Right- SF9007R
Arm Pad- SF9007GT
4. Joystick- S350169701
5. Joystick Arm- [3]
6. Joystick Screw Set- SI69001

7. Joystick Adj. Knob- [3]
8. Seat Plate- SEC9009
9. Seat Post- SF9010
10. Frame [3]
11. Articulating Bar- [3]
12. Footrest- SPCGW05
13. Fork Bearings- [3]
14. Fork- [3]
15. Axle & Spacer- SEC9028 [2]
16. Caster Gray- SEC9031
17. Drive Wheel- SEC9060
18. Wire Harness- SF9038GT [2]
19. Fender Decorations
19. Red- SEC9050R [2]
19. Blue- [3]
19. Gold- [3]

20. Main Shroud- SEC9050
21. Battery Cover- SEC9051
22. Controller- S350169739
23. Recline Hinge Cover- [3]
24. Charger- [3]
25. Motor & Gear Box Lt- SEC9015L [4]
25. Motor & Gear Box Rt- SEC9015R [4]
26. Rear Axle Bolt- SEC9061
27. Anti-Tippers (pr)- SEC9045
28. Batteries (2 needed)- LRA402208
29. Battery Strap- [3]
30. Hub Cap- [3]
31. Caster Bearings (pr)- SEC9031B
32. Motor Key Way- [3]

[1] Stock is gone use part # ECI6004.

[2] When stock is gone the part is no longer available.

[3] Stock is gone the part is no longer available.

[4] Any motor before serial # 2W301B119000869 must be replaced in pairs.

02/10/2021