




SUNFIRE PLUS GT




1. Headrest- SF2010 [1]
2. Seat Assembly 18"- SF9008GT
20"- ECI6004
3. Armrest Left- SF9007L
Right- SF9007R
- 3a. Armrest Pad- SF9007GT
4. Joystick- S350169701
5. Joystick Screw Set- SI69001
6. Joystick Arm Knob- SP3CG2W02
7. Seat Post- SF9010
8. Footrest- SP3CG2W05
9. Caster- SF9031
10. Drive Wheel- SF9060
11. Wire Harness- SF9038GT
12. Controller- S350169734 [1]
13. Charger- SP3CG2W14
Charger AC Cord- SF9056
14. Motor/Gear Box Left- SF9015LGT
Motor/Gear Box Right- SF9015RGT
15. Anti-Tippers- SF9045
16. Batteries (2 needed)- 402208



To attach ELRs

- A. ECELR-P
- B. SFELR
- C. LK3JELR



[1] When stock is gone the part is no longer available.

02/11/2021

[Diagrams](#)

[Home](#)