


[1] When stock is gone, parts will be no longer available.

1. Controller- SF9001 **[1]**
2. Armrest Assy. Right- **[2]**
3. Armrest Assy. Left- **[3]**
4. Headrest- SF2010 **[1]**
5. Seat Post- SF9010
6. Anti-Tipper- SF9045
7. Motor Right- SF9015R **[1]**
Motor Brush- SF5012
8. Drive Wheel- SF9060
9. Caster- SF9031
10. Wire Harness- SF9038 **[1]**
11. Foot Rest- SF9035
12. Battery (2 needed)- LRA402208
13. Charging Power Cable- SF9056
14. "L" Mounting Bracket- SF9018



[1] When stock is gone the part is no longer available.

[2] Stock is gone use SF9007R

[3] Stock is gone use SF9007L

02/11/2021