


WILDCAT 450

1. Arm Pad- LPW180011 [1]
2. Drive Wheel- LPW180009 [1]
3. Motor/Brake/Gear Box Assembly
3. Left- LPW182002-01 [1]
3. Right- LPW182002-02 [1]
4. Caster Fork- LPW182007 [1]
5. Footrest Right- LPW180018R [1]
6. Battery (2 needed)- LRB402308
7. Joystick- LPW182001
7. Joystick Overlay- DKREMD01
7. Joystick Cable- LPW182001-02 [1]
8. Lower Battery Box- LPW182006-02 [1]
9. Seat Belt- LDR410017
10. Charger- LRA402210



[1] When stock is gone the part is no longer available.

02/18/2021

[Diagrams](#)

[Home](#)