




- 1. Arm Pad- LPW180011 [1]
- 2. Drive Wheel- LPW180009 [1]
- 3. Motor/Brake/Gear Box Left- LPW180002-01 [1]
Motor/Brake/Gear Box Rt- LPW180002-01
- 4. Caster- LDR410007
- 5. Caster Bearing- STDS2A4314

- 6. Battery-LRA402208
- 7. Joystick- LPW180001
- 8. Battery Case Bottom- LPW180006-02 [1]
- 9. Seat Belt- LDR410017
- 10. Charger- LRA402209
- 11. Right Footrest- LPW180018R

[1] When stock is gone the part is no longer available.

02/18/2021