

# AP8000



specifications	
Pump Power	110 VAC, 60 Hz
Pump Airflow	4 LPM
Pump Cycle Time	5 Minutes
Pump Type	Compressor
Pump Alarms	No
Mattress Dimensions	80" x 36" x 8"
Mattress Cover Material	Fluid-Resistant Nylon
Weight Capacity	350 lbs.
Warranty	1 Year Limited

- The AP8000 alternating pressure mattress system is the most economical 8" therapeutic support surface in the market offering optimum pressure redistribution for the prevention and treatment of pressure ulcers.
- Quiet user friendly pump delivers air volume of 4 LPM.
- The removable top cover is low shear, anti-microbial and fluid resistant.
- 8" Alternating cell depth.
- CPR valve allows for rapid deflation.
- 350 lb. weight capacity.

1.	Pump Only	AP8120
2.	Mattress Complete	AP8125
3.	Cover	AP8140
4.	Bladder	AP8150

Part Photo Search

[Diagrams](#)

[Home](#)