


Daytona 3



Wheel

- | | |
|-----------------------------|------------------------------------|
| 1. Wire Harness - S31048 | 9. Tiller Bearing Set - S31021 |
| 2. Handle Bar - S31002 | 10. Controller - S31010 |
| 3. Throttle - S12-034-00300 | 11. Motor Brush - SF5012 |
| 4. Speed Pot - S31001 | 12. Axle Key - S31030-KEY |
| 5. Hand Grip - S31003 | Motor Coupler - S31097GT |
| 6. Battery Gauge - S31208 | 13. Rear Wheel - S31035 |
| 7. Key Switch - S31006 | 14. Battery (2 needed) - LRM412123 |
| Key - S3KEY | 15. Battery Charger - LRA402210 |
| 8. Headlight - S31022 | |

12/07/20

Diagrams

Home