


Daytona 3GT



1. Hand Grip - S31003
2. a. Tiller Cover Black - S31209BK
b. Speed Pot - S31001
c. Horn Button - S31012
3. Throttle - S12-034-00300
4. Head Light - S31022
5. Key Switch - S31006
Key - S3KEY
6. Tiller Adjustment Knob - S31028
7. Tiller Bearing Assembly - S31021
8. Front Shroud Red - S31015GTR
9. Front Wheel - S31038GT
10. Front Frame - S31099GT
11. Battery (2 needed) - LRM412123
12. Battery Charger - LRA402210

13. Motor Brush - SF5012
Transaxle - S31031GT
Motor/Transaxle - S31030GT [1]
Motor Coupler - S31097GT
Axle Key Way - S31030-KEY
14. Rear Wheel - S31035
15. a. Rear Shroud Red - S31018GTR
Rear Shroud Blue - S31018GTB
b. Rear Bumper - S31019GT
16. Seat Post - S31053
17. Adjustment Knob - S31052
18. Charger Port - S31147
19. Rear Axle Bolt - S350016910
20. Electro-Magnetic Brake - S31033BGT
21. Anti-Tipper Wheel - S31026

[1] When stock is gone the part is no longer available.

12/07/20

[Diagrams](#)

[Home](#)