

PHOENIX 4


1. Seat Assembly- S31045N
2. Armrest (subs) Left- SF9007L
Right- SF9007R
Arm Pad For Sub- SF9007GT
3. Seat Plate- ~~S31039~~ **NLA**
4. T Handle Frame Lock- S350169716
5. Tension Knob- S350169715
6. Handgrip- S31003
7. **Upper Tiller Cover**
Digital Battery Gauge- S31004GT
Horn Button- S31012
8. Speed Dial Assembly- S31001
Speed Dial Knob- LRM412103-11
9. Headlight Assembly- S31022
Head Light Switch- S31122
10. Throttle Old Style- ~~S31000~~ *
New Style- S12-034-00300

***Stock is gone replace with new style.**

11. Basket Mount- BC35051
12. Key Switch- S31006
Key- S3KEY
13. Seat Post- S31049
14. Seat Pin- S350169717
15. Wire Harness- S31048
16. Steering Axle- 350102 **[1]**
17. Tiller Adjustment Knob- S31028
18. Steering Post- S350101 **[1]**
19. Tie Rod- S350103
20. Front Bumper- S350104
21. Tiller Swivel Assembly- S350109
22. Main Shroud- ~~S350100~~ **NLA**
Decoration Blue Left- S3501BL
Right- S3501BR **[1]**
Red Left- S3501RL
Right- S3501RR
23. Foot Shroud- S350107
Decoration Blue Left- S3501BL-1
Right- S3501BR-1
Red Left- S3501RL-1
Right- S3501RR-1
24. Tiller Bearing Assembly- S31021
25. Charger- LRM302111
26. Battery (2 needed)- LRM402109
Battery Pack- S350169707
Case Only- S350169706
Battery Contacts- S31060

27. Rear Shroud- S350111
Decoration Blue Left- S35011B-L
Right- S35011B-R
Red Left- S35011R-L
Right- S35011R-R
28. Anti-Tipper Wheel- S31026
29. Rear Frame- S31130-4
30. Front Frame- S350108
31. Rear Wheel- S31132B
32. Rear Axle Nut- S350016910
33. Motor & Brake- S31031
Motor, Brake & Transaxle- S31125
34. Arm Adjustment Knob- S31052
35. Controller- S31008GT
36. Frame Contact Front- S31056
Rear- S31057 **[1]**
37. Front Wheel- S350100
Front Wheel Bearing- S310100-B

38. Transaxle Mounting Bracket- ~~S350112B~~ **NLA**
39. Battery Gauge Analog- S31208
Digital- S31208B
40. Basket- S35003-P
41. Seat- BC35045
Must discard seat plate and use existing seat plate.
42. Armrest Left- BC35040L
Right- BC35040R
Ampad- BC953014
43. Seat- S35045
44. Armrest Left- BC35040L
Right- BC35040R
Ampad- BC953014
45. Freewheel Knob- S31031K
46. Rear Lifting Handle- S35011FH



[1] When stock is gone the stock is no longer available.

01/14/2021