


SPITFIRE EX2 (4 WHEEL)



1. Control Panel Assembly- LRM414100
2. Lower Control Panel- LRM414102
- 2a. Tiller Charger Port- LRM414103
3. Ctrl. Panel Decorations Blue- LRM412103-03B
Ctrl. Panel Decorations Red- LRM412103-03R
4. Speed Control- LRM412103-10
Speed Control Knob- LRM412103-11
- 5a. Horn & Horn Switch- LRM412103-08
- 5b. Headlight Switch- LRM414103-16
6. Battery Gauge- LRM412103-LED
7. Key Switch- LRM402323-02-HL
Key- LRM112110-02
8. PC Board- LRM402104-02
9. Hand Grip Cover- LRM412101-11
10. Throttle Assembly- LRM302102-02
Reverse Throw Throttle- LRM302102-02RT
11. Front Compartment Door- LRM412123-05
12. Basket- LRA402226
13. Basket Bracket- LRM412101-08
14. Tiller Boot- LRM402200-06

29. Seat Assembly- L
Back Cover- LRM
Seat Cover- LRM
30. Seat Plate- LDR42
a. Seat Swivel Knob
b. Seat Accessory K
31. Seat Post- LRM41
32. Seat Post Pin- S3
33. 12ah Battery (2 ne
12ah Battery Cas
12ah Battery Cas
21ah Battery (2 n
21ah Battery Cas
21ah Battery Cas
34. Battery Box Hand
35. 20 Amp Circuit Br
Substitute- LRM30
36. Battery Box Charg
37. Front Frame- LRM
38. Anti-Tipper Wheel

-
- 15. Steering Post- LRM412131-01
 - 16. Tiller Bearing Set- LRB402348
 - 17. Tiller Adjustment Knob- LRA402201-23
 - 18. Front Wheel- LRM412112-02
 - 19. Hub Cap- LRM412202-HC
 - 20. Decoration Set (inc. pt. #3) Blue- LRM412117B
 Decoration Set (inc. pt. #3) Red- LRM412117R
 - 21. Front Shroud Only- LRM412118
 - 22. Front Shroud Trim- LRM414108
 - 22a. Head Light Assembly- LRM414107
 - 23. Tie Rod (ea)- LRM412130
 Tie Rod Spacer- LRM412130-01
 - 24. Steering Knuckle Left- LRM414105
 Steering Knuckle Right- LRM414104
 - 25. Articulating Beam- LRM414106
 - 26. Floor Mat Left- LRM412118L
 Floor Mat Right- LRM412118R
 - 27. Armrest Pad- LRM302115-01
 - 28. Armrest Assembly Left- LRM302115L-02
 Armrest Assembly Right- LRM302115R-02
 Armrest Knob- S31052

- 39. Rear Shroud- LRM412112-02
- 40. Rear Wheel- LRM412112-02
- 41. Battery Case Clamp- LRM302115-01
- 42. Motor/Gearbox/Braze- LRM302115-01
 Motor- LRM302115-01
 Gearbox- LRM412112-02
 Brake- LRM412112-02
- 43. Rear Frame- LRM412112-02
- 44. Controller- LRM412112-02
- 45. Charger- LRM302115-01
- 46. Frame Latch- LRM412112-02
- 47. Wire Harness- LRM412112-02
- 48. Front Wheel Bearing- LRM412112-02
- 49. Battery Contact Pad- LRM302115-01
- 50. Tiller Lock- LRM412112-02
- 51. Lower Tiller Swivel- LRM412112-02
- 52. Seat Post Screw- LRM412112-02
- 53. Contact Block Set- LRM302115-01
- 54. Tiller- LRM414111

[Diagrams](#)

[Home](#)

Spitfire EX2 (4 wheel)