

10200



1. Push button - 12002PB
2. Leg Only - 10200L
3. Tips - RTL10386GB or RTL10386G-8

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10200-4	17"	19"	21.75"	16.5"	39"	32"	6.5 lbs.	28 lbs.	350 lbs.

08/18/16