

10202-2



1. Push Bottom - 12002PB
2. Leg - 10200L
3. Tips - RTL10386GB or RTL10386G-8

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10202-2	17"	20.5"	23"	18.5"	39"	32"	6 lbs.	14 lbs.	300 lbs.

01/24/18