

1. Push Button - 12002PB
2. Tips - RTL10386GB or RTL10386G-8
3. 3" Walker Wheels Pr. - 10111
4. Leg Only - 10200L



SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10207-2	17"	20"	22"	15.5"	32"	25"	7 lbs.	17 lbs.	350 lbs.

08/19/16