



1. Straight handle pr - 10219S
2. Knob - 10219CL
3. Leg/wheel left - 10219CWL
Leg/wheel right - 10219CWR
4. Front leg - 10219FL
5. Rear leg - 10219RL
6. Tips - RTL10386GB or
RTL10386G-8

Hardware pack - 10219HRDWR
 One directional front wheels -
 10219WK

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Depth (opened at base)	Handle Height Range	Height (highest setting)	Height (lowest Setting)	Product Weight	Carton Shipping Weight	Weight Limit
10219	Adjustable	24"	25"	29" - 37"	37"	29"	8 lbs.	11 lbs.	175 lbs.