


10227



1. Push Button- 12002PB
2. Leg with 5" Wheels (pr)- 10109
3. Leg with Mushroom Cap (ea)- 10200L-MCB
4. Mushroom Cap (pr)- 10107MC-B
5. Leg with Rubber Tip (ea)- 10200L
6. Rubber Tip- 10325B



SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10227-4	18"	20.5"	25.5"	19"	33.5"	26.5"	6.5 lbs.	30 lbs.	350 lbs.

10/27/17

[Diagrams](#)
[Home](#)