



10233

ORIGINAL SEAT

1. Seat [1]
2. Rear Leg & Glide Brake- 10233RL
3. Walker Tip Grey- RTL10386GB
Black- 10325B
4. Leg & Wheel- 10109
5. Seat Hinge Pin [1]
6. Adjustment button- 12002PB
7. Basket- 10230B

[1] Part is no longer available.



SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Seat Width	Seat Depth	Seat Height Highest	Seat Height Lowest	Handle Ht. (highest setting)	Handle Ht. (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10233	18.5"	21.5"	27"	24"	19"	11"	23.5"	19.5"	39"	32"	11 lbs.	10.5 lbs.	300 lbs.

06/21/2021

[Diagrams](#)

[Home](#)