



10243

ORIGINAL SEAT

1. Hand Brake- 10230BRAKE
2. Cable- 950340129901
3. Rear Wheel- 10243RW
4. Rear Leg Assemble Left- 10243LAL
Right- 10243LAR
5. Front Leg w/Fork- 10243F
6. Front Wheel- 10243W
7. Caster Lock Pin- 10243CPA
8. Adjustment Button- 12002PB
9. Seat- **No Longer Available**
10. Seat Hinge Pin- **No Longer Available**
11. Basket- 10230B



SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Seat Width	Seat Depth	Seat Height Highest	Seat Height Lowest	Handle Ht. (highest setting)	Handle Ht. (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10243	19"	21.5"	27"	26.5"	19"	11"	26"	23"	39"	35"	16 lbs.	19 lbs.	300 lbs.

10/09/18

[Diagrams](#)[Home](#)