

10248N



1. Push Button - 12002PB
2. Tips - RTL10386GB or RTL10386G-8
3. Height Adj. Bar (Pr)- 10253CB
(Includes Orange Clip)

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10248N	17"	19"	22"	17"	38"	28"	6.5 lbs.	7.5 lbs.	300 lbs.

10/17/16