

10253

UNIVERSAL ADULT/JR WALKER



1. Height Adj. Bar (pr)- 10253CB
(includes orange clips)
2. Rear Legs (pr)- 10253-RL
3. Front Legs & Wheels (pr)- 10128
4. Adjustment Button- 12002PB
5. Glide Caps (pr)- 10107MC-G

SPECIFICATIONS

Width	17" inside handgrip/18.75" inside back legs
Overall Width	24"
Depth (Opened at Base)	17.75"
Height	28.75" - 38.5"
Weight	7 lbs.
Carton Shipping Weight (10253-1/10253-4)	8/30.5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

10/17/16

[Diagrams](#)

[Home](#)