

990X KNEE WALKER



- | | |
|---|---|
| 1. Handlebar- 990X-01 [1] | 13. Hub Cap- [2] |
| 2. Brake Lever- 990X-02 [1] | 14. Wheel Bearings- STDS1012 |
| 3. Handgrip- 9503H79002 | 15. Lower Seat Post Assembly- 990X-08 [1] |
| 4. Handlebar Adj Clamp- 9505B79006 | 16. Water Bottle w/ Mount- 990X-09 [1] |
| 5. Handlebar Column- 9503H79005 | 17. Knee Pad Post Clamp- 990X-10 [1] |
| 6. Brake Cable Guide- 9503H79008 | 18. Rear Wheel Assembly Left- [2] |
| 7. Brake Cable- 990X-04 [1] | Right- 990X-12R [1] |
| 8. Steering Column- 9503H79007 | 19. Brake Caliper- 990X-13 [1] |
| 9. Column to Fork Locking Bolt- 9503H79012 | 20. Brake Disc- 990X-14 [1] |
| 10. Stem Bearings- 9503H79013 | 21. Rear Shroud- 990X-15 [1] |
| 11. Front Wheel Assembly Left- 990X-07L [1] | 22. Upper Kneepad Post Assembly- [2] |
| Right- 990X-07R [1] | Upper Kneepad Post Assembly- 990X-11N [1] |
| 12. Axle Bolt- 990X-06 [1] | Serial # 3H1803001101 & higher. |
| | 23. Knee Pad- 990X-16 [1] |

[1] When stock is gone the part is no longer available

[1] When stock is gone the part is no longer available.
[2] The part is no longer available.

7/19/2021

Dimensions

Diagrams

Home