


CRUISER II (4S SERIAL #)



- | | |
|---|--|
| 1. Handgrip- STDS3J2403 | 24. Release Lever- STDS3J2419 |
| 2. Rear Wheel- STDS3J2424 | 25. Adj. Ht Arm Assy. Full Left- STDS3J24AFL |
| 3. Hand Rim- STDS3J2400 | Adj. Ht Arm Assy. Full Right- STDS3J24AFR |
| 4. Hand Rim Hardware (9 ea)- STDS004609-9 | Adj. Ht Arm Assy. Desk Left- STDS3J24AL |
| 5. Rear Seat Saddle- STDS3J2414 | Adj. Ht Arm Assy. Desk Right- STDS3J24AR |
| 6. Rear Wheel Bearing- STDS3J2413 | 25a. Upper Adj. Arm Full Lt- STDS3J24AFL-UFL * |
| 7. Rear Axle- STDS3J2425 | Upper Adj. Arm Full Rt- STDS3J24AFR-UFR * |
| 8. Hub Cap- STDS3J2421 | Upper Adj. Arm Desk Lt- STDS3J24AL-UDL * |
| 9. Push to Lock Wheel Lock Left- STDS4S221L | Upper Adj. Arm Desk Rt- STDS3J24AL-UDR * |
| Push to Lock Wheel Lock Right- STDS4S221R | * Includes arm pad. |
| Pull to Lock Wheel Lock Left- STDSK234-PTL | 26. Adj. Height Release Left- STDS2A0046 |
| Pull to Lock Wheel Lock Right- STDSK234-PTR | Adj. Height Release Right- STDS2A0044 |
| Wheel Lock Tip- STDS8M4603 | 27. Arm Pad Desk- STDS3J2419 |
| 10. Front Seat Saddle- STDS3J2412 | Arm Pad Full- STDS3J2419 |
| 11. Hanger Pins (4 ea)- STDS4S1069 | 28. Back Upholstery 16"- STDS4S2406-AD-V |
| 12. Bumper- STDS3J2450 | Back Upholstery 18"- STDS4S2408-AD-V |
| 13. Footplate & Heel Loop Left- STDS3J2410L | Back Upholstery 20"- STDS4S2402-AD-V |

Footplate & Heel Loop Right- STDS3J2410R
Heel Loop & Peg- STDS831
Heel Loop Peg- STDS3J24SF-HP
14. Footrest Adj. Knob- STDS3J4502
15. Caster- STDS4S2432
16. Caster Axle- STDS8M4631
17. Caster Bearing- STDS4S3602
18. Fork- STDS8M4633
19. Fork Stem Bearing- STDS8M4601
20. Fork Stem Cap- STDS8M4651
21. Side Frame Left- STDS5D24L1
Side Frame Right- STDS5D24R1
22. Cross Frame 16"- STDS3J24X6-AD
Cross Frame 18"- STDS3J24X8-AD
Cross Frame 20"- STDS3J242X-AD
23. Seat Upholstery 16"- STDS4S2416-AD-V
Seat Upholstery 18"- STDS4S2418-AD-V
Seat Upholstery 20"- STDS4S2420-AD-V

29. Flip Back Arm Brkt- STDS3J2409
30. Cross Frame Bolt- STDS001061-K3
31. Frame Cap- STDS3J2422
32. Cross Frame Assist 16"- STDS3J2436
Cross Frame Assist 18"- STDS3J2438
Cross Frame Assist 20"- STDS3J2440
33. Rear Arm Latch- STDS5D2499
34. Seat Extension Screw (pr)- STDS5D24XF-S
Seat Extension & Screw (pr)- STDS5D24XF-ET

Footrest (pr)- STDS3J24SF
Elevating Leg Rest (pr)- LK3JELR
Anti Tippers- STDS824
Seat Belt- STDS850-PB
Brake Extensions- STDS601-K2

[Diagrams](#)

[Home](#)