


## Sentra ECHD 20" x 18" for s/n Y and 4S

1. Back Upholstery- STDS1S2630 20"
2. Handgrip- STDS871
3. Back Upholstery Screw (8 ea)- STDS5D5504
4. Rear Wheel- STDS1S1000HD
5. Hub Cap- STDS1000C
6. Rear Axle Bolt- STDS1053
7. Rear Axle Bearing- STDS1S1012
8. Hand rim- STDS005008
10. Hand Rim Bolt & Spacer- STDS1016-8
11. Cross Frame- 20" x 18"  
STDS4S2620-HD
12. Cross Brace bolt- STDS001061
13. Brake Tip- STDS8M4603
14. Brake- STDS4019R (right), L (left)
15. Front Axle Bolt- STDS8M4631
16. Axle Bearings- STDS1011HD
17. Caster- STDS1075-2
18. Front Fork- STDS1S3045
19. Fork Stem Bearings- STDS8M4601
20. Fork Stem Dust Cap- STDS1051
21. Seat Upholstery- 20" x 18"  
STDS2328BLK
22. Seat Upholstery Screws- STDS001062
23. Arm Assm.- Desk STDSDDAECR or L  
Full STDSDFAECR or L  
Adj. Desk- STDSADDASPR or L  
Adj. Full- STDSADFASPR or L
24. Arm Pad- Desk STDSDSKVB  
Full STDSFLVB
25. Arm Rest Panel- Desk STDS002044, Full STDS002043
26. Arm Release Button- STDSSB
27. Slide Tube- STDS005200
28. Footrest- STDSF-TF
29. Footplate- STDS00004A
30. Rubber Bumper- STDS1026



06/01/18

Diagrams

Home