

JBSTD18FB-FF & JBSTD22FB-FF

JET BLUE WHEELCHAIRS

**Only parts listed below
are available.**

1. Back Upholstery- STDS2AJB 02 (22")
2. Hand Grip- STDS871
3. Fork- STDS2AJB14
4. Frame Cap- STDS3J2422
5. Fork Stem Cap- STDS1051
6. Fork Nut- STDS001045
7. Stem Bearings- STDS1060
8. Front Caster- STDS1002
9. Axle bearings- STDS1011
10. Axle bolt- STDS1002-1
11. Hand Rim Screws & Spacers- STDS1073
12. Hand Rim- STDS004012
13. Rear Axle Bolt- STDS1005
14. Hub Cap- STDS1000C
15. Rear Wheel- STDS1006
16. Rear Axle Bearing- STDS1012
17. Footrest Set- STDF



04/22/2021

[Diagrams](#)

[Home](#)