




15901C - PCB301 LONG TERM CARE BED



1. Head Deck Grey- SP01-PCB300-0104AE
Yellow- [3]

2. Foot Deck Grey- SP01-PCB300-01-55E
Yellow- [3]

3. Thigh Deck Grey- SP01-PCB300-0105AE
Yellow- [3]

4. Head Caster Assy Left Grey- SP01-PCB300-01-58E
Yellow- SP01-PCB300-01-58Y [2]
Right Grey- SP01-PCB300-01-60E
Yellow- SP01-PCB300-01-60Y [2]

5. Foot Brake Assy Left Grey- SP01-PCB300-01-59E
Yellow- [3]
Right Grey- SP01-PCB300-01-61E
Yellow- [3]

6. Head/Foot/Thigh Deck Hdw- SP01-PCB300-015301
7. Head Foot Bolts Male- SP01-UF16108030E
Female- SP01-UN90108090E

8. Power Cord- SP01-CA65738

Grey

Yellow

9. Strain Relief Bracket- SP01-PCB300-01501E
Strain Relief Screws- SP01-PCB300-01500E

10. Wall Guard Gray- SP01-PCB300-010701
Yellow- [3]

11. Head/Foot End Mattress Retainer
Gray- SP01-PCB300-010403
Yellow- SP01-PCB300-01040Y [2]

12. Head Board Mtg Brkt
a. 2 Hole Gray- SP01-PCB300-01-62E
b. 3 Hole Yellow- SP01-PCB300-01-62Y [2]

13. Side Mattress Retainer Gray- SP01-PCB300-015302
Yellow- [3]

14. Foot Board Mtg Brkt
a. 2 Hole Gray- SP01-PCB300-01-63E
b. 3 Hole Yellow- SP01-PCB300-01-63Y [2]

15. Foot End Ratchet Assy- SP01-PCB300-01-57E

16. Locking Caster- SP01-TL500-37E

17. Non-Locking Caster- SP01-TL500-29E

18. Hand Pendant Gray- SP01-72298
Yellow- SP01-PCB601-03E-DM
Extendant Cord Pendant- SP01-40662

19. Control Box- SP01-66702

20. Hi/Lo Actuator- SP01-72487

21. Foot Actuator- SP01-72497

22. Head Actuator- SP01-72492
Cotter Pin for Actuators- SP01-PCB300-015002

23. Fixed Assist Bar- P300-03AE [1]
Rotating Assist Bar- P300-07AE [1]
Rotating Assist Rail- P300-08AE [1]

Side/Assist Rail Pin Assy- SP01-PCB300-52-01E
Side/Assist Rail Pivot Bolt- SP01-ML29160430E01

24. Caster Bales Gray- SP01-PCB300-011202
Yellow- [3]

25. Hardware Pack- SP01-PCB300-01-64E

26. Transport Bracket- SP01-PCB300-06E

[1] Accessory not a part.

[2] When stock is gone the part is no longer available.

[3] Part is no longer available.

07/16/20

[Diagrams](#)

[Home](#)