


Only parts listed are available.



1. Handgrip Left- 950330039901

Right- 950330049901

2. Brake (top mount)- 9503HBRAKE [1]

3. Cable Only- 950000199901 [1]

4. Anti-Rattle- 950000869901

10. Lower Pivot Bolt- 9500002051 [1]

11. Fork Stem Bearing-STDS1060

12. Dust Cap- STDS1051

13. Upper Pivot- 95000010001 [1]

- 5. Height Adjustment Knob
 - 1 Piece- 10208KN
 - 2 Piece- 10208K
- 6. Cable Clip- 10208C
- 7. Cable Fastener- 9500004101
- 8. Rear Axle Bolt Cap- STDS001000
- 9. Wheel Assembly- 9500H30001

- 13. Upper Pivot- 950000619901 [1]
- 14. Cable Adjuster- 9500001232
- 15. Tote Bag Black- 808
- 16. Tray- 816T
- 17. Wheel Bearing- STDS1011

[1] When stock is gone the part is no longer available.

07/12/2021

[Diagrams](#)

[Home](#)