

ALL CONTROLLERS & PROGRAMMERS POWER CHAIRS

UNIT	ELECTRONICS	PROGRAMMER
CIRRUS PLUS EC	OLD OR NEW VSI	SP-5
CIRRUS PLUS HD	NEW VSI	SP-5
COBALT	DYNAMIC JOYSTICK	DC DX-HHP
COBALT	DYNAMIC SHARK	DC DX-HHP
COBALT	NEW VSI	SP-5
DENALI	VR2	SP-5
GEO	VR2	SP-5
GLADIATOR	CURTIS 1227-2402	BC31000
IMAGE (ORIGINAL S# DMAN)	DYNAMIC SHARK	DC DX-HHP
IMAGE (ORIGINAL S# DMAN)	VR2	SP-5
IMAGE (S# 2W)	VR2	SP-5
IMAGE EC	VR2 GREY CONNECTION	SP-5
IMAGE EC	VR2 YELLOW CONNECTION	SP-5
IMAGE GT	VR2	SP-5
INTREPID	DYNAMIC JOYSTICK	DC DX-HHP
INTREPID	NEW VSI	SP-5
MEDALIST	DYNAMIC JOYSTICK	DC DX-HHP
MEDALIST	P & G JOYSTICK	SP-5
MEDALIST 450	DYNAMIC JOYSTICK	DC DX-HHP
MEDALIST 450	P & G JOYSTICK	SP-5
RENEGADE	DYNAMIC JOYSTICK	DC DX-HHP
RENEGADE	NEW VSI	SP-5
SUNFIRE GENERAL	VR2	SP-5
SUNFIRE GLADIATOR	VR2	SP-5
SUNFIRE LS	DYNAMIC SHARK	DC DX-HHP
SUNFIRE PLUS	OLD VSI	SP-5
SUNFIRE PLUS EC	VR2	SP-5

TITAN	DYNAMIC SHARK	DC DX-HHP
TITAN LTE	DYNAMIC LINX	DLX-HKEY01-A *
TITAN AXS	DYNAMIC LINX	DLX-HKEY01-A *
TRIDENT	VR2	SP-5
TRIDENT HD	DYNAMIC LINX	DLX-HKEY01-A *
TRIDENT HD	VR2	SP-5
WILDCAT	OLD OR NEW VSI	SP-5
WILDCAT 450	VR2	SP-5

* DC DX-HHP can also be used.

[Home](#)