


# Pure expressions™

ELECTRIC BREAST PUMPS  
TECHNOLOGY

# Pure expressions™

ELECTRIC BREAST PUMPS




## True2Nature™ Expression

- Stimulation and Expression modes imitate the natural breastfeeding experience



## Digital LCD Display

- Displays 9 different suction strengths and allows easy adjustment for each breast individually\*
- 30-minute timer and memory function which auto-saves last pumping session



## Anti-Backflow Technology\*

- The closed-loop, anti-back-flow system guarantees no milk backup into the tubing or motor.



**BP001**  
Breast Pump Carry Bag



**BP002**  
Insulated Cooler Bag\*\*



**BP004**  
6oz Wide Neck Storage Bottle with Cap



**BP005/BP0013**  
27mm & 31mm Pump Kit without Insert



**BP006/BP0014**  
24mm & 27mm Silicone Insert\*\*



**BP007**  
Valve Set with Membrane



**BP008**  
Membrane 4pcs per Bag



**BP009**  
Tubing Set



**BP0011**  
6V Power Adapter



**BP0012**  
6V External Battery Pack\*



**BP0015**  
Feeding Nipple Set\*\*



**BP0016**  
Bottle Base

# Pure expressions™

ELECTRIC BREAST PUMPS



©2018 Medical Depot, Inc. dba Drive Medical Design and Manufacturing®, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical Design and Manufacturing are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies.

**FLEXIBILITY TO CUSTOMIZE**  
suction strength for each breast



# Pure expressions™

ELECTRIC BREAST PUMPS



\* On Models RTLBP2000 only

\*\* Sold Separately

\* See inside for RTLBP0200 accessory options

## DELUXE DUAL CHANNEL ELECTRIC BREAST PUMP



Drive Item #	Description	UOM
RTLBP2000	Dual Channel Electric Breast Pump	1/ea

- LCD Display provides status of programs and functions
- Independent and custom suction for each breast with 9 adjustable suction strengths
- 2 Phase Expression; Stimulation and Expression Mode
- Easy to assemble and clean; quiet for discreet pumping
- Closed Loop Anti-backflow System
- Single or Double Pumping options
- BPA and Latex Free
- A/C Power Supply Included; Internal battery power source; 4 AA alkaline batteries (not included)
- One Year Limited Warranty

### Specifications

Product Dimensions	4" (W) x 6.25" (L) x 2.75" (H)
Product Weight	.9 lbs
Retail Box Dimensions	RTLBP2000 12.2" (W) x 7.5" (L) x 11" (H)
Retail Box Weight	RTLBP2000 3 lbs 4 oz
Batteries	4 AA Alkaline
Contains Latex	No
Warranty	One Year Limited

### Accessories

BP001	Pure Expressions Carry Bag	1/ea
BP002	Pure Expressions Insulated Cooler Bag	1/ea
BP004	Pure Expressions 6oz. Storage Bottle	1/ea
BP005	27MM Pump Kit w/out Insert	1/pr
BP006	24MM Silicone Insert	1/pr
BP007	Valve Set - w/membrane	1/pr
BP008	Membrane (4pcs)	2/pr
BP009	Tubing Set	1/pr
BP0011	6V Power Adapter	1/ea
BP0012	6V External Battery Pack	1/pr
BP0013	31MM Pump Kit w/out Insert	1/ea
BP0014	27MM Silicone Cushion Insert	1/ea
BP0015	Feeding Nipple Set	1/ea
BP0016	Base (base to hold bottles)	1/ea

## DUAL CHANNEL ELECTRIC BREAST PUMP



Drive Item #	Description	UOM
RTLBP0200	Dual Channel Electric Breast Pump	1/ea

- LCD Display provides status of programs and functions
- Control suction and speed to suit your comfort level with a total of 28 pumping state options
- 2 Phase Expression; Stimulation and Expression Mode
- Easy to assemble and clean; quiet for discreet pumping
- Passive prevention of milk backflow
- Single or Double Pumping Options
- BPA and Latex Free
- Piston Pump System imitates baby's suckling rhythm
- Includes A/C Power Supply & External battery pack; 6 AA alkaline batteries (not included)
- One Year Limited Warranty

### Specifications

Product Dimensions	7" (W) x 4.75" (L) x 3.5" (H)
Product Weight	.93 lbs.
Retail Box Dimensions	10.35" (W) x 6.02" (L) x 8.66" (H)
Retail Box Weight	3 lbs.
Batteries	6 AA Alkaline
Contains Latex	No
Warranty	One Year Limited

### Accessories

BP001	Pure Expressions Carry Bag	1/ea
BP002	Pure Expressions Insulated Cooler Bag	1/ea
MQ9401	24MM Breast Pump Set	1/pr
MQ9402	27MM Breast Pump Set	1/pr
MQ9403	Silicone Insert	1/pr
MQ9404	Valve Set - w/membrane	1/pr
MQ9405	Membrane (4pcs)	2/pr
MQ9406	4 oz Storage Bottle	1/ea
MQ9407	Tubing Set	1/ea
MQ9408	Power Adapter	1/ea
MQ9409	Feeding Nipple Set	1/ea
MQ9130C	Car Charger	1/ea

## PURE EXPRESSIONS PRODUCT SUITE COMPARISON CHART



Model	Pure Expressions RTLBP2000	Pure Expressions RTLBP0200	Pure Expressions RTLBP1000
Item Description	Double Electric Breast Pump	Economy Double Electric Breast Pump	Single Electric Breast Pump
Adjustable Vacuum	Yes	Yes	Yes
BPA Free	Yes	Yes	Yes
Latex Free	Yes	Yes	Yes
Double Electric Pumping	Yes	Yes	Yes
Single Electric Pumping	Yes	Yes	Yes
Closed System	Yes	No	Yes
Expression Phases	2	2	2
Vacuum Chamber	2	1	1
Vacuum Strength	30 mm Hg ~ 270 mm Hg	30 mm Hg ~ 285 mm Hg	30 mm Hg ~ 270 mm Hg
Suction Level	9	5	9
Independent Suction Strength for each breast	Yes	No	No
Clock/ Timer	Yes (30 min countdown)	Yes (30 min countdown)	Yes (30 min countdown)
Memory Function	Yes (settings last used)	Yes (settings last used)	Yes (settings last used)
LCD Display	Yes	Yes	Yes
Portability	Yes	Yes	Yes
Pump Dimensions	4" x 6.5" x 8"	3.5" x 5" x 7.2"	2.2" x 4.2" x 5"
Pump Weight	0.82 lbs (13.1oz)	0.92 lbs (14.74 oz)	0.54 lbs (8.64 oz)
Weight with box	3 lbs 4oz	2.70 lbs	1.92 lbs
AC Adapter	Yes	Yes	Yes
Battery Option	External battery pack (4 AA batteries)	External battery pack (6 AA batteries)	Internal Battery Pack (4 AA batteries)
Breast Shield Sizes	(1) 27mm, 24 mm	(1) 24mm	(1) 27mm, 24 mm
Silicone Cushion Insert	No- (2) 24mm	No	No- (1) 24mm
Breastmilk Containers	2 (4 oz)	2 (4 oz)	1 (4 oz)
Carrying Bag	Yes- Carry bag w/zipper	No	Yes-drawstring bag
Insulated Cooler Bag	No	No	No
Tubing Length	36	36	36

Drive DeVilbiss Healthcare provides a comprehensive line of quality breast pump products for the nursing mother. According to the Affordable Care Act (ACA) most healthcare plans must provide breastfeeding equipment and counseling for pregnant and nursing women. These services may be provided before and after birth. For more information please contact your healthcare insurance provider.