

NEBULIZER TO FILTER LISTED BY MODEL NUMBER

<u>MODEL #</u>	<u>NAME</u>	<u>FILTER #</u>
4056	Eagle	18090F
7090	Aeroneb Go	No Filter
7091	Aeroneb Go	No Filter
18002	Power Neb 2	18025F
18016	Hercules Beetle	No Filter
18018	Voyager Pro	No Filter
18070	Pacific Elite	18090F
18071	Pacific Elite	18090F
18080	Power Neb Ultra	18090F
18081	Power Neb Ultra	18090F
46502	PulmoMate	5650D-601 or 5650D-602
18040-P	Checker	18090F
18040-Y	Checker	18090F
18041-P	Checker	18090F
18041-Y	Checker	18090F
18062	Pacifica 2	18090F
18090-BE	Beagle	18090F
18090-FS	Fish	18090F
18090-PB	Panda	18090F
18090-PG	Penguin	18090F
18091-BE	Beagle	18090F
18091-FS	Fish	18090F
18091-PB	Panda	18090F
18091-PG	Penguin	18090F
2655D	Sunrise Nebulizer	2655D-601
3050-2	Sport Neb	3050-2-FILTER
3650D	Pulmo-Aide LT	5650D-601 or 5650D-602
3655D	Pulmo-Aide Compact	3655D-601 or 3655D-602
3655DX	PulmoNeb	3655D-601 or 3655D-602
3655LT	PulmoNeb LT	3655LT-601 or

		3655LT-602
3655LTR	PulmoNeb LT	3655D-601 or 3655D-602
3655TX	PulmoNeb	3655D-601 or 3655D-602
4650D	PulmoMate	5650D-601 or 5650D-602
5610D	Pulmo-Aide	6710D-607 *** 5650D-601 or 5650D-602
5650D	Pulmo-Aide	5650D-601 or 5650D-602
5650H	Pulmo-Aide	6910D-605 6710D-607 ***
6910D	Traveler	7025D-609
8650D	Heavy Duty Compressor	8650D-606 8650D-607
MQ0911	Fire Truck	MQ0054
MQ5300	Airial Voyager	MQ0054
MQ5500	Airial Voyager	MQ0054
MQ5501	Airial Voyager	MQ0054
MQ5600	MEDNEB	MQ0054
MQ5700	MEDNEB+	MQ0054
MQ5800	Compact Compressor	MQ0057 *
MQ5800	Compact Compressor	MQ92692 **
MQ5900	Compact Compressor	MQ0057 *
MQ5900	Compact Compressor	MQ92692 **
MQ6000	Compact Style	MQ0054
MQ6001	Duck	MQ0054
MQ6002	Penguin	MQ0054
MQ6003	Panda	MQ0054
MQ6004	Panda	MQ0054
MQ6005	Panda Pink	MQ0054
MQ6006	Penguin	MQ0054
MQ6008	Airial Express	MQ0054
MQ7000	Townhouse	MQ0054

- * Filter located bottom of unit, cover has 5 little holes in it.
- ** Old style, cylinder shape with small cap.
- *** 6 filters w/filter cap.

[Home](#)