



# 2026

## Power Mobility Catalog

Enhancing quality  
of life is our  
driving force

**drive**   
MEDICAL

## Smooth and Stylish Power Mobility

Our power mobility solutions utilize cutting-edge technologies blended with thoughtful design to ensure user comfort, convenience, and safety. From ergonomic hand controls to sensitive tiller adjustment, every feature is designed down to the last detail to ensure a smooth and effortless ride. Padded armrests and customizable swivel, reclining, and fold-down seating options complete the experience. Many options even offer color customization, to reflect your own personal style! If you're looking for the ultimate in sophisticated power mobility, turn to Drive Medical.



## CarbonTouch™ Power Chair

Item # CFPC17FS



- Includes a built-in USB charger to keep your devices powered on the go
- Comes standard with phone and cup holder
- Folds in seconds with a single touch for easier transport
- Easily switch between powered mode and manual mode for uninterrupted mobility
- Removable joystick can be mounted on either armrest for user preference
- Two charging modes — charge directly through the controller port or use off-board charging
- Equipped with safety belt and anti-tippers
- Breathable mesh seat and backrest offer ergonomic support with enhanced airflow for everyday comfort
- Flip-back armrest with depth adjustability

[Shop Now](#)



## Scout Compact 3-Wheel Scooter

Item # SFSCOUT3



- Quick and easy disassembly
- Interchangeable color panels
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy freewheel operation

[Shop Now](#)



## Scout LT 4-Wheel Scooter

Item #	Battery
SCOUTLT4	12 AH
SCOUTLT4-EXT	20 AH



- A fun, easy to disassemble 4-wheel scooter for the budget conscious shopper
- Delta Tiller for one-handed operation
- Adjustable armrests
- Generously-sized front basket for storage of personal items
- Includes easily interchangeable red and blue shrouds
- Up to 14 miles per charge
- Non-marking, solid flat-free tires
- Padded Delta tiller
- Battery Size: 12 AH - SCOUTLT4
- Battery Size: 20 AH - SCOUTLT4-EXT

[Shop Now](#)



## Phoenix LT 4-Wheel Scooter

Item #	Seat Size
PHOENIXLT4	18"
PHOENIXLT4-20	20"



- 4-wheel design and anti-tip wheels assist with overall stability for a more secure experience
- Height-adjustable seat with flip-back and adjustable-width arms
- Interchangeable color panels
- 350 lbs weight capacity
- Seat swivels 180 degrees for easy transfer
- Headlight for increased visibility
- Breaks down to 5 pieces for easy transport
- Comes with rugged 9"x3" wheels
- Padded Delta tiller

[Shop Now](#)



## Prowler 4-Wheel Mobility Scooter

Item # PROWLER3410MG22CS

- Adjustable delta tiller with ergonomic throttle control for better grip and additional steering comfort
- Height-adjustable swivel seat features fold-down backrest making storage easier.
- Reclining seat back and adjustable headrest
- Seat position can be adjusted forward or backward
- Pneumatic tires and full suspension allow for a smooth ride
- Large plastic carry basket to take items on-the-go
- Programmable controller to adjust driving parameters based on user's needs
- Rearview mirrors , anti-tip wheels, and front headlight for additional safety
- Additional hand-operated drum brake

[Shop Now](#)



## Panther 4-Wheel Heavy Duty Scooter

Item # PANTHER22CS

- Front and rear suspension on 14" pneumatic tires  
One-touch infinite tiller adjustment
- Built-in USB port for charging your mobile devices
- Extra-bright LED lighting for superb visibility
- Ergonomic delta control handles for ease of steering and throttle control
- Large front plastic basket with carry handle
- Captain's seat is height adjustable and features a fold-down, reclining backrest with adjustable headrest
- LED headlights, parking lights, and turn signals
- Dual rear view mirrors
- Anti-tip wheels with easy freewheel mode

[Shop Now](#)



## ZooMe 3-Wheel Recreational Scooter

Item # ZOOME3

- Powerful hub driven motor with speeds up to 15 mph
- Folding tiller for easy transportation or for storage
- Two keys included
- Hand brake with locking mechanism for parking on inclines
- Forward and reverse controls are easily changed with the switch of a button
- Fast and slow speed settings easily changed with the switch of a button
- Front basket for storage
- Pneumatic tires for the smoothest ride possible
- Extra bright LED headlight
- Anti-tip wheel



[Shop Now](#)

## ZooMe Auto-Flex Folding Travel Scooter

Item # FLEX-AUTO

- Fully electronic folding and unfolding
- Aluminum alloy frame is lightweight
- Airline safe lithium battery
- High-gloss shrouds for a protective finish
- Easy throttle control
- Delta tiller for users with limited dexterity
- Auto-folding seat with padded, adjustable armrests
- Flat-free, non-marking tires for worry-free travel
- Easy-to-adjust tiller places the scooter's controls in the ideal position
- Secondary manual-release levers and toggle switch
- Comes with two ignition keys and one key fob
- Front and rear reflectors and anti-tip wheels



[Shop Now](#)

## Rear Basket Set, LType,

Item # SF8020

- Roomy rear basket boasts strong steel construction
- Ideal for carrying large or heavy items that do not fit in front baskets
- Dimensions: 11.75"(L) x 16.75"(W) x 11.25"(H)

[Shop Now](#)



## Cup Holder

Item # AZ0060

- Plastic construction
- Take your favorite hot or cold beverage with you anywhere
- Easily attaches to your product's armrest
- Fits all Drive Medical Power Scooters except Spitfire Series

[Shop Now](#)



# Power Mobility

## Oxygen Holder Set

Item # SF8010

- Allows you to carry your oxygen tank with you
- Securely holds "D" and "E" size oxygen cylinders
- Sturdy metal frame

[Shop Now](#)





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  [www.drivemedical.com](http://www.drivemedical.com)

CANADA OFFICE  256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7  866-752-1363

©2026 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.