


THERAPEUTIC SUPPORT SURFACES CATALOG



GROUP I & II Support Surface Qualifiers CMS Coding Eligibility

GROUP I QUALIFIER

IS THE PATIENT COMPLETELY IMMOBILE?

NO

YES

Does the patient have limited mobility **OR** any stage pressure ulcer on the trunk or pelvis **AND** one of the following?

1. Altered sensory perception
2. Compromised circulatory status
3. Fecal or urinary incontinence
4. Impaired nutritional status

NO

No
Reimbursement
Available

YES

PATIENT QUALIFIES FOR A GROUP I SUPPORT SURFACE

Code	Description
E0180 E0181 E0182	Power pressure relieving mattress overlays (APP's)
E0184	Foam Mattress
E0185	Gel Overlay

GROUP II QUALIFIER

DOES THE PATIENT HAVE?

- Multiple Stage II pressure ulcers located on the trunk or pelvis; **AND**
- Patient has been on a comprehensive ulcer treatment program for at least the past month which has included the use of an appropriate Group I surface; **AND**
- The ulcers have worsened or remained the same over the past month; **OR**
- Large or multiple stage III or IV pressure ulcer(s) on the trunk or pelvis; **OR**
- Recent myocutaneous flap or skin graft for a pressure ulcer on the trunk or pelvis (surgery within the last 60 days); **AND**
- The patient has been on a group 2 or 3 support surface immediately prior to a recent discharge from a hospital or nursing facility (discharge within the past 30 days).

NO

No
Reimbursement
Available

YES

PATIENT QUALIFIES FOR A GROUP II SUPPORT SURFACE

Code	Description
E0277	Alternating Mattress Replacement System

Medicare Criteria: This information is not intended to be, nor should it be considered, billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare programs, and should consult an attorney or other advisor to discuss specific situations in further detail.

THERAPEUTIC SUPPORT SURFACES

Homecare Mattresses	4
Pressure Redistribution Foam Mattresses	5
Gel Overlays	6
Alternating Pressure Pump & Pad Systems	6
Alternating Pressure and Low Air Loss Mattress Systems	8
Perimeter Covers	12
Fall Prevention (Rail Pads & Floor Mats)	13
Wheelchair Cushions	14

 **WARNING**

Whenever a support surface and bed rails are used as part of any hospital or home care bed system, the potential for bed rail entrapment exists. To reduce the risk, a patient must always be assessed regarding their risk of entrapment before any product is used. Routine patient monitoring is also required once a product is in use. To learn more about the risk of bed rail entrapment, contact your Roscoe Medical sales representative for a copy of our "Bed Rail Safety Guide" or visit the U.S. FDA's website at <http://www.fda.gov> and search for "Bed Rail Entrapment".



Fiber Core Homecare Mattress



- 100% virgin densified polyester fiber core is more economical than foam
- Lighter weight than innerspring mattress, reduces shipping costs
- Breathable polypropylene cover keeps patient cool and comfortable
- Mattress and cover are odor free and non-toxic
- Polypropylene cover is fluid, stain and mold resistant, as well as 100% recyclable
- Meets Federal Fire Code 16 CFR 1633
- One-year limited warranty



Item #	Description	Size	Wt. Cap.	UOM
PBHBMATFB	Fiber Core Mattress with Polypropylene Cover	80" x 36" x 6"	250 lbs.	each



Firm Support Innerspring Mattress



- For use in home care, nursing homes and hospitals where firm support is needed
- Easy-to-clean vinyl cover is durable and fluid-resistant
- High quality 312 bonnell coil provides resilient support base
- 2oz densified polyester fiber topper adds firmness and support
- Heavy-duty resinated pad with a specialized mesh woven back
- Meets 16 CFR 1633 flammability standard
- One-year limited warranty



Item #	Description	Size	Wt. Cap.	UOM
PBHBMATIS	Firm Support Innerspring Mattress	80" x 36" x 6"	300 lbs.	each
PBHBMAT84	Firm Support Innerspring Mattress	84" x 36" x 6"	300 lbs	each

**HCPCS
E0184**



SATINFOAM Group 1 Compressed Foam Mattress




GROUP 1

- Premium quality foam homecare mattress
- Easy-to-clean, sheer-resistant nylon cover is durable and fluid-resistant
- Non-slip underside keeps mattress firmly in place
- Super compressed packaging; box footprint is only 14" x 15" x 38"
- Meets 16 CFR 1633 flammability standard
- One-year limited warranty



Item #	Description	Size	Wt. Cap.	UOM
PBSATINFOAM3680	Group 1 Foam Mattress	80" x 36" x 6"	300 lbs.	each

**HCPCS
E0184**



Aruba-200 Group 1 Contour Compressed Mattress



GROUP 1

- The sculptured foam surface is designed to provide pressure relief and comfort in three areas: the head and shoulders, the torso and the lower extremities
- The cross-cut design maximizes airflow and provides comfort, support and pressure redistribution
- Deluxe nylon cover is durable, zippered, fluid-resistant and anti-microbial
- Super compressed packaging; box footprint is only 14" x 15" x 38"
- The nylon top cover reduces friction, heat and contamination. The bottom cover is reinforced vinyl with non-slip protection for safety
- Meets 16 CFR 1633 flammability standard
- Two year warranty on foam, one year on cover



Item #	Description	Size	Wt. Cap.	UOM
ARUBA-200	Aruba Foam Mattress	80" x 36" x 6"	350 lbs.	each

**HCPCS
E0184**



Aruba-42 Group 1 Bariatric Compressed Foam Mattress



GROUP 1

- Compressed Group 1 Foam to help reduce and prevent pressure sores
- Multi-density foam layers for patient comfort
- Smaller footprint for easier delivery and shipping savings
- Supports those up to 600 lbs.
- Fluid-resistant zippered cover for easier cleaning
- One-year limited warranty

Item #	Description	Size	Wt. Cap.	UOM
ARUBA-42	Aruba Bariatric Foam Mattress	80" x 42" x 6"	600 lbs.	each



Gel-Lite Gel/Foam Overlay

GROUP 1

- Independent gel tubes prevent migration of gel in high pressure areas
- Raises mattress height for patients
- Triple seal polyurethane gel tube construction ensures seal integrity
- High-density foam perimeter on all sides helps prevent patient roll-out
- Conforms to body contours and allows air circulation to reduce heat
- Convoluted foam topper provides additional comfort
- Cover is vapor permeable and moisture resistant; low shear/low friction
- Weight bearing capacity: 250 lbs
- Meets flammability standards TB 117-2013
- One-year limited warranty



Item #	Description	Size	Wt. Cap.	UOM
GO-78343-C	Gel-Lite Overlay	78" x 34" x 3"	250 lbs.	each



Meridian SatinAir™ APP System

GROUP 1

- Interlocking bubble cells exchange air every 5 minutes
- Heavy duty 16 gauge PVC
- End flaps secure pad to mattress
- 130 individual bubble cells for maximum comfort and effectiveness
- High capacity output 5 liters/minute
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range
- Cycle time: 5 minutes
- Illuminated on/off switch
- Built-in bed mount hooks
- HCPCS Codes: E0181 (pad and pump)
- HCPCS Codes: A4640 (replacement pad)
- One-year limited warranty on pump, 90 days limited warranty on mattress pad



Item #	Description	Size	Wt. Cap.	UOM
MER-APP2	Meridian SatinAir APP System	78" x 34" x 2.5"	300 lbs.	each
APP-BPD	Basic Replacement Pad	78" x 34" x 2.5"	300 lbs.	10/cs
APP-DLXP	Deluxe Replacement Pad with Flaps	78" x 34" x 2.5"	300 lbs.	10/cs

A photograph of an elderly man with grey hair, wearing a light blue hospital gown, sitting in a hospital bed. He is smiling broadly and looking towards the right. A woman with short grey hair, wearing a bright pink t-shirt, stands next to him, also smiling and looking towards the camera. The background shows a hospital room with medical equipment and a white curtain.

CHOOSING THE RIGHT SURFACE

Pressure ulcers are high on the list of "never events" developed by the Center for Medicare and Medicare (CMS) stating that these events should never happen. We carry a broad line of systems with superior quality and exceptional value to meet your patients' individual needs.

HPCPS
E0277



Meridian SatinAir™ 8” APM

GROUP 2

ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM

- Includes mattress, cover and pump
- Pillow function maintains the air in cells at the head for patient comfort
- Air cells can be easily removed and replaced
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- 8 liter-per-minute pump in an easy-to-use design
- Two modes of therapy - static or alternating
- Variable pressure setting for maximum comfort
- Support patients up to 350 lbs.
- One-year limited warranty on mattress and pump



Item #	Description	Size	Wt. Cap.	UOM
MER-SAPM2	Meridian SatinAir™ 8” APM	80" x 36" x 8"	350 lbs.	each
X1SAPM201	Replacement Mattress	80" x 36" x 8"	350 lbs.	each
PBSATCOVER	Replacement Cover	80" x 36" x 8"	N/A	each

HPCPS
E0277



Meridian SatinAir™ 5” + 3” Foam Base

GROUP 2

5” ALTERNATING PRESSURE AND LOW AIR LOSS MATTRESS SYSTEM WITH 3” FOAM BASE

- Includes mattress, cover and pump
- 3-inch foam base prevents “bottoming out” in the event of a power outage
- Pillow function maintains the air in cells at the head for patient comfort
- Air cells can be easily removed and replaced
- CPR quick deflate at head of mattress
- Quilted nylon zipper cover is low shear, fluid-resistant, and vapor permeable
- 8 liter-per-minute pump in an easy-to-use design
- Two modes of therapy - static or alternating
- Variable pressure setting for maximum comfort
- Pump features both audio and visual alarms
- Support patients up to 350 lbs.
- One-year limited warranty on mattress and pump



Item #	Description	Size	Wt. Cap.	UOM
MER-SAPM53	Meridian SatinAir™ 5-3	80" x 36" x 8"	350 lbs.	each
X1SAPM5301	Replacement Mattress	80" x 36" x 8"	350 lbs.	each
X1SAPM5302	Replacement Cover	80" x 36" x 8"	N/A	each



GROUP 2

Ultra-Care™ 4800

STANDARD ALTERNATING PRESSURE/LOW AIR LOSS MATTRESS SYSTEM

HCPCS
E0277



ULTRACARE



- Includes mattress, cover and pump
- 8" alternating pressure mattress with low air loss
- Quick disconnect cells for easy cell replacement
- Air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- 8 LPM pump in an easy-to-use design
- 10 minute cycle time
- Front panel switch provides two modes of therapy - static or alternating
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range: 20 mmHg – 70 mmHg
- Supports patients up to 350 lbs.
- Two-year limited warranty on mattress and pump
- Optional upgrade to digital pump and 450 lb. weight capacity (Order 4800E)



Item #	Description	Size	Wt. Cap.	UOM
4800	Ultra-Care 4800 System	80" x 36" x 8"	350 lbs.	each
4800E	Ultra-Care 4800 System with 12 LPM Pump	80" x 36" x 8"	450 lbs.	each
4802 DC	Replacement Mattress	80" x 36" x 8"	450 lbs.	each
4802-C-88	Replacement Cover	80" x 36" x 8"	N/A	each

Magic-Air® Pump

The Quick-Fill Pump is designed to quickly inflate and deflate Low Air Loss Mattresses reducing your patient set-up time by 65%, and allowing you time to attend to your patient's other needs. This quick 51 cubic feet per minute flow pump allows your technician to get the mattress up and therapeutic in minutes.



- Light weight, compact and all-steel construction
- Pumps air at a rate of 51 CFM
- Quieter as compared to standard plastic air pumps
- Weight: 3 lbs.
- Dimensions: 13" length x 5.5" wide height
- 12 foot cord
- Works with Quick Fill Valve 90153 (sold separately)
- One-year limited warranty

Item #	Description	UOM
109-516811	Magic-Air® Pump	each
90153	Quick Fill Valve	each



GROUP 2

**HPCPS
E0277**



ULTRACARE
○○○○○○○○○○○○○○○○



Ultra-Care™ 4500

ALTERNATING PRESSURE/LOW-AIR LOSS MATTRESS SYSTEM WITH 3-INCH INTEGRATED FOAM BASE

- System includes mattress, cover and pump
- Features a 5" alternating pressure upper surface combined with an integrated 3" foam base to provide a full eight inches of therapeutic support with bottom-out protection in the event of a power failure
- Separate compartment for foam base inside polyurethane cover prevents soiling
- Quick disconnect cells for easy cell replacement
- Air cells made with nylon/polyurethane for increased durability and reliability
- CPR quick deflate located at the head of the mattress
- Quilted cover, moisture resistant, vapor permeable nylon; low shear/low friction
- 8 LPM pump in an easy-to-use design
- 10 minute cycle time
- Front panel switch provides two modes of therapy - static or alternating
- Variable pressure setting allows for maximum comfort
- Adjustable pressure range: 20 mmHg – 70 mmHg
- Two-year limited warranty on mattress and pump
- Optional upgrade to digital pump and 425 lb. weight capacity (Order 4500E)



Item #	Description	Size	Wt. Cap.	UOM
MER-4500	Ultra-Care 4500 Excel	80" x 36" x 8"	325 lbs.	each
4500E	Ultra-Care Excel System with 12 LPM Pump	80" x 36" x 8"	425 lbs.	each
4502	Replacement Mattress	80" x 36" x 8"	425 lbs.	each
4502-C-88	Replacement Cover	80" x 36" x 8"	N/A	each

Standard Alternating Pressure Pump



GROUP 2

The Manual Pump for APMs offers both alternating and static functions to provide patients with optimal treatment of pressure ulcers and maximum comfort. It flows at 8 liters-per-minute in a simple, easy-to-use design.



- Front panel switch provides two modes of therapy - static or alternating - for maximum comfort
- LED backlit switches with dial on panel for simple set up and operation
- Includes retractable hooks for bed use
- 8 LPM pump
- Supports patients up to 350 lbs
- Two-year limited warranty



Item #	Description	UOM
4581	Standard APM Pump	each



Safe-T-Guard™ Cover

PERIMETER COVER TO REDUCE RISK OF FALLS AND ENTRAPMENT

- Reduces the risk of falls and entrapment
- Fits most standard-size hospital mattresses as well as standard eight inch alternating air mattresses
- Secures with three quick-release straps (positioned top, bottom and middle) that connect underneath the mattress
- The quilted top is water resistant and has a covered foam border that rises approximately four inches above the mattress surface and is sloped inward to reduce the risk of roll outs
- Includes cover and 4 foam inserts with separation to allow for patient egress
- Nylon, polyurethane construction to reduce heat build-up
- One-year limited warranty

SGU UNIVERSAL DESIGN



- Secures to mattress with three quick-release straps positioned top, bottom and middle for a snug fit
- Fits most standard hospital bed mattresses 36" W x 80" L x 6" H

SG8 DESIGN (FOR 8" APM SYSTEMS)



- Fits most 8" air mattress systems
- Low shear, nylon zippered cover



LOOKING FOR SOMETHING SPECIFIC?

Your Compass Health Representative is a knowledgeable, trained professional. If you can't find what you need, assistance is just a phone call away. Call **1-800-376-7263** or visit our online catalog at **www.compasshealthbrands.com**.

Item #	Description	L x W	Qty.
SGU_FG	Universal Safe-T-Guard Cover with QR Straps	80" x 35.5"	each
SG8_FG	Safe-T-Guard Zipper Cover for 8" Air	80" x 35.5"	each
SGU42_FG	Safe-T-Guard Cover for 42" Mattress	80" x 42"	each



Full



Half

Bed Rail Pads

- Protects against bumps, bruises
- Vinyl cover wipes clean for easy disinfecting
- Pads do not interfere with the operation of the bed rails
- Attach with hook & loop straps
- One-year limited warranty



Item #	Description	Size	UOM
30038	Full-Length Bed Rail Pads	48" x 15" x 1"	pair
BR50	Half-Rail Bed Pads	31" x 16" x 1"	pair



High Density Bedside Foam Fall Mat

- Help reduce the possibility of injuries due to falls
- Easy-to-clean vinyl cover
- High density, high compression polyurethane foam
- Flame retardant meets with California Technical Bulletin #117
- Bi-Fold design held together with hook & loop fasteners
- One-year limited warranty



Item #	Description	Size	UOM
30036	High Density Foam Fall Mat	66" x 24" x 2"	each



Meridian Fall Mat

- PVC Material offers tear resistance, durability and is easy to clean and disinfect
- Non-skid bottom helps keep the mat safely in place
- Anti-fatigue properties provide additional comfort for caregivers
- For use when the patient is in the bed
- Ships in brown box - 20.9" X 8.7" X 8.7"
- One-year limited warranty




Item #	Description	Size	UOM
FM7220	Meridian Medical Fall Mat	72" x 20" x 3/4"	each

WHEELCHAIR CUSHIONS

If a patient has a wheelchair for which they meet Medicare coverage criteria, they also qualify for a General Use Seat and Back Cushion (see page 18).



**HPCPS
E2601**




Meridian General Use Foam Cushion

- Economical seating surface provides comfort and pressure reduction
- Designed to assist in the prevention of pressure ulcers
- Cushion contours the seating area
- Removable, waterproof nylon top with non-slip vinyl bottom
- Color: navy
- Weight bearing: See chart below
- One-year limited warranty
- TAA Compliant options available, ask your Compass Health Representative



ITEM #	SIZE	Wt. Cap.	UOM	HPCPS
WCF-163	16" x 16" x 3"	275 lbs.	each	E2601
WCF-183	18" x 16" x 3"	300 lbs.	each	E2601
WCF-203	20" x 16" x 3"	325 lbs.	each	E2601
WCF-2083	20" x 18" x 3"	350 lbs.	each	E2601
WCBF-223	22" x 16" x 3"	400 lbs.	each	E2602
WCBF-2283	22" x 18" x 3"	425 lbs.	each	E2602
WCBF-2483	24" x 18" x 3"	450 lbs.	each	E2602
GEL-182	18" x 16" x 2"		each	N/A

**HCPCS
E2603**



Meridian Optimum Comfort Gel Seat Cushion



- Unique channeled gel bladder prevents gel migration in high-pressure areas
- Gel channels provide relief over the entire cushion area
- Foam base provides additional support
- Inner foam core has a liner to protect foam's integrity
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-9B 16X16	16" x 16"	275 lbs.	each	E2603
HDC-9B 18X16	18" x 16"	300 lbs.	each	E2603
HDC-9B 18X18	18" x 18"	325 lbs.	each	E2603
HDC-9B 20X16	20" x 16"	350 lbs.	each	E2603
HDC-9B 20X18	20" x 18"	400 lbs.	each	E2603
HDC-9B 22X18	22" x 18"	425 lbs.	each	E2604
HDC-9B 24X18	24" x 18"	450 lbs.	each	E2604

**HCPCS
E2611**



Meridian General Use Back Cushion

- Provides upper back and lumbar support
- Designed to reduce lower back pain and muscle tension
- Removable, waterproof nylon cover
- Stabilization board provides enhanced rigidity and support
- Non-skid backing
- Color: Black
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
BKF-1617	16" x 17"	NA	each	E2611
BKF-1817	18" x 17"	NA	each	E2611
BKF-2019	20" x 19"	NA	each	E2611
BKFB-2219	22" x 19"	NA	each	E2612

HCPCS
E2607



Meridian Optimum Pommel Cushion



- Contoured base and pommel
- Four grades of foam
- Inner core has a liner to protect foam
- Cover is a 4 way stretch polyester/polyurethane
- Moisture repellent & low shear/friction
- Weight bearing: See chart below
- Proudly made in the USA
- One-year limited warranty



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-POM1616	16" x 16" x 3"	275 lbs.	each	E2607
HDC-POM1816	18" x 16" x 3"	275 lbs.	each	E2607
HDC-POM2016	20" x 16" x 3"	275 lbs.	each	E2607

HCPCS
E0190



Meridian Optimum Coccyx Cushion



- Unique design for relief of pressure on coccyx area
- Weight bearing: See chart below
- One-year limited warranty
- Proudly made in the USA



ITEM #	SIZE	Wt. Cap.	UOM	HCPCS
HDC-CX1616	16" x 16" x 3"	225 lbs.	each	E0190
HDC-CX1816	18" x 16" x 3"	225 lbs.	each	E0190
HDC-CX2016	20" x 16" x 3"	225 lbs.	each	E0190

Medicare Wheelchair Seat Cushion CMS Coding Eligibility



Medicare Criteria: This information is not intended to be, nor should it be considered, billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare programs, and should consult an attorney or other advisor to discuss specific situations in further detail.



6753 Engle Road • Middleburg Hts., OH 44130 Toll Free: (800) 376-7263 • Fax: 440-572-4261

www.compasshealthbrands.com

rev_04_2025