



# 2026

## Walkers Catalog

Enhancing quality  
of life is our  
driving force

**drive**<sup>®</sup>   
MEDICAL

## Independent Mobility

Support your mobility and safely enhance your independence with a walker from Drive Medical! We offer a variety of styles and features—including wheeled walkers, knee walkers, and walkers with built-in seating. Many of our walkers are designed to be operated by users with limited manual dexterity, and a wide variety are lightweight and feature a compact fold that makes them portable and easy to transport. Explore our products and find the perfect walker to complement your busy lifestyle!

### Table of Contents

Walkers	3-13
Accessories	14-23



# Walkers

## AIR-GO™ WALKER



### Walk on Air™

The Air-Go Walker redefines mobility—delivering comfort, ease, and freedom beyond traditional walkers. It's comfortable on the hands, stable underfoot, and ready when it's time to rest. With oversized, cushioned Airgonomic™ handles, water-dispersing treaded wheels, and a padded mesh seat, it gives users both freedom of movement and confidence in every step. For you, it delivers showroom differentiation, stronger customer loyalty, and reliable margins with MAP pricing you can count on.

### The support of a walker meets the convenience of a rollator



#### Comfortable Hold, Every Time

Oversized Airgonomic™ handles cushion hands and are easy to grasp — even for users with arthritis or reduced dexterity.



#### Easily Rest

Enjoy the freedom to sit and rest when needed on the ¾-length open-and-go padded mesh seat—designed shorter to let users walk naturally—while four height adjustments ensure a comfortable, customized fit.



#### Compact, Stable, Travel-Ready

The scissor-fold frame collapses in seconds for storage and transport, and the 6" water-dispersing treaded wheels with rear ski glides provide smooth, stable movement across various terrain.



[Shop Now](#)

#### Accessories



#### Storage Bag

Item # RTL102662-BAG



#### Cup Holder

Item # RTL102662-CUP



#### Cane Holder

Item # RTL30100-CANE



#### Product Specs

Height Range	Adjustable for users 4'9" - 6'7"
Product Dimensions	25.6" (W) x 25.2" - 26.9" (L) x 33.1" - 36" (H)
Seat Height	20" - 23"
Caster Dimensions	6" (D) x 1" (W)
Folded Dimensions	7.5" (W) x 26.2" (L) x 20.5" (H)
Product Weight	8 lb
Weight Capacity	300 lb
Warranty	Limited Lifetime
Retail Box Dimensions	21.5" (L) x 7.1" (W) x 28.5" (H)

Item #	Description	UPC
MD30100BK	Air-Go Walker, Onyx <span style="color: black;">●</span>	822383031071
MD30100BL	Air-Go Walker, Sapphire <span style="color: blue;">●</span>	822383031088

# Walkers

## Deluxe Folding Walker, Two Button



Item #	Description	Packaging
10200-1	Adult	Each
10200-4	Adult	Case

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight

[Shop Now](#)



## Folding Walker, Two Button with 5" Wheels



Item #	Description	Packaging
10210-1	Adult	Each
10210-4ASM	Adult	Case

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter construction ensures maximum strength while remaining lightweight
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Includes wheels
- Polypropylene (PP) and PVC wheels allow walker to roll easily and smoothly over irregular surfaces

[Shop Now](#)



# Walkers

## Folding Walker, Two Button with 5" Wheels

Item #	Description	Packaging
10244-4	Adult	Case

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter steel construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Comes with wheels



[Shop Now](#)

## Deluxe Folding Walker, Two Button

Item #	Description	Packaging
10201-4	Youth	Case

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight



[Shop Now](#)

# Walkers

## Deluxe Folding Walker, Two Button with 5" Wheels



Item #	Description	Packaging
10211-4ASM	Youth	Case

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Includes wheels
- Polypropylene (PP) and PVC wheels allow walker to roll easily and smoothly over irregular surfaces

[Shop Now](#)



## Folding Junior Walker, Two Button with 5" Wheels



Item #	Description	Packaging
10247-4	Youth	Case

- Easy push-button release mechanisms
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter steel construction ensures maximum strength
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Vinyl-contoured hand grip
- Comes with wheels

[Shop Now](#)



# Walkers

## Deluxe, Trigger Release Folding Walker

Item #	Description	Packaging
10224-4	Adult	Case

- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity
- Trigger release can be reached with hands of all sizes
- Durable, Composite Trigger Release resists cracking
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- Patented

[Shop Now](#)



## Deluxe, Trigger Release Folding Walker with 5" Wheels

Item #	Description	Description
10226-1	Adult	Single
10226-4	Adult	Case

- Sturdy aluminum construction ensures maximum strength while remaining lightweight
- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle
- Ribbed, contoured handgrips relieve pressure on hands
- Easily adjust the height in 1" increments for a customized fit
- 5" wheels are ideal for maneuvering both indoors and outdoors
- Rear glide cap allows for walker to slide easily
- Folds for easy transport with detachable sides
- Quick and tool-free assembly

[Shop Now](#)



# Walkers

## Deluxe, Trigger Release Folding Walker



Item #	Description	Packaging
10225-4	Youth	Case

- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle and is ideal for individuals with limited finger dexterity (Figure A) Trigger release can be reached with hands of all sizes
- Durable, Composite Trigger Release resists cracking
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight
- Patented

[Shop Now](#)



## Deluxe, Trigger Release Folding Walker with 5" Wheels



Item #	Description	Packaging
10227-4	Youth	Case

- Sturdy aluminum construction ensures maximum strength while remaining lightweight
- Trigger release feature for opening and closing allows users to close the walker without removing their hands from the handle
- Ribbed, contoured handgrips relieve pressure on hands
- Easily adjust the height in 1" increments for a customized fit
- 5" wheels are ideal for maneuvering both indoors and outdoors
- Rear glide cap allows for walker to slide easily
- Folds for easy transport with detachable sides
- Quick and tool-free assembly

[Shop Now](#)



## Bariatric Aluminum Folding Walker, Two Button



Item #	Description	Description
10220-1	Adult	Single

- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lb
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip



[Shop Now](#)

## Bariatric Aluminum Folding Walker, Two Button



Item #	Description	Description
10220-1WW	Adult	Single

- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- Vinyl-contoured hand grip



[Shop Now](#)

# Walkers

## Deluxe Folding Travel Walker with 5" Wheels

Item #	Description	Packaging
RTL10263KDR	Adult	Single

- Easily folds tool-free to ½ the size of a traditional walker and can be stored in a "Ready Set Go" Pouch
- Legs are attached with strong elastic cords so legs can be folded
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces
- Two side-storage pouches with large pockets for personal items are held in place by hook and loop fastener tabs to both side frames



[Shop Now](#)

## Clever-Lite LS, Adult Walker

Item #	Description	Packaging
10233	Adult	Single

- Rear glide push-down brakes
- Sturdy 1" diameter anodized, extruded aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker by raising seat to step inside the frame
- Easily folds with patented dual-lever side-paddle release (Figure B)
- Soft, flexible backrest provides comfort and security while seated
- 5" fixed front wheels



[Shop Now](#)

# Walkers

## Clever-Lite Walker, Adult, with 8" Casters



Item #	Description	Packaging
10243	Adult	Single

- Front wheels can be set in either "swivel" or "fixed" position
- Special loop lock made of internal aluminum casting operates easily and ensures safety
- Sturdy 1" diameter anodized aluminum construction
- Comes with flip-up seat with built-in carry handle
- Allows individual to be seated or use as traditional walker
- Folds with patented side-paddle release



[Shop Now](#)

## Two Wheeled Walker with Seat



Item #	Description	Packaging
1239RD	Adult	Single

- Features comfortable and convenient padded seat and backrest
- Quickly folds when not in use
- Rear leg tips act as brakes when pressed down
- Height adjustable handles
- Large 5" swivel wheels for added maneuverability
- Contoured handles for a comfortable grip



[Shop Now](#)

# Walkers

## Side Walker

Item #	Description	Packaging
10240-1	Adult	Single



- Designed for individuals with the use of only one hand or arm
- Lighter than a walker and more stable than a cane
- Folds easily with one hand
- Ideal for users who do not require a walker but need a wider base for support

[Shop Now](#)



## Steerable Knee Walker

Item #
796



- Knee walker can be steered for increased maneuverability
- Ideal for individuals recovering from foot surgery, breaks, sprains, amputation and ulcers of the foot
- New, two-piece padded leg pad provides comfort
- Deluxe braking system
- Removable front basket
- Durable, steel frame with silver vein finish
- Tool free height adjustment
- 7.5" casters are ideal for indoor/outdoor use

[Shop Now](#)



## Nitro® Glide Knee Walker

### Item #

791BL  
791RD

### Description

Blue  
Red



- Width expanding front axels (16", 18", 20") for adjustment of turning radius
- Easy-to-replace brake cables enclosed in frame tubing
- Tool-free height adjustment of handle ( 30"- 39") and knee pad (17" - 23").
- Offset, reversible knee pad for left or right leg use
- Includes removable storage basket with handle
- Lightweight and foldable frame for convenient storage and transport
- Simple-to-operate hand brakes control dual rear braking system

[Shop Now](#)



## Economy Folding Knee Walker

### Item #

RTL799



- Easily steerable, while tight turning radius
- Dual rear on-wheel braking system provides immediate stopping ability and control while navigating
- Push-button parking brake incorporated into handbrake for when you need to rest safely
- Padded leg support provides comfort, stability and is easily compatible with left or right leg
- Quickly and easily adjust the height of the kneepad or handlebar to cater to your individual needs, no tools required
- Foldable design allows for convenient storage and portability
- Strong, durable frame supports users up to 350 lb
- 8" wheels are ideal for indoor and outdoor use

[Shop Now](#)



# M+D Air-Go™ Walker Accessories

## Storage Bag

**Item #**

RTL102662-BAG

- Spacious storage bag can be used on the M+D by Drive Air-Go™ Walker

[Shop Now](#)



## Cane Holder

**Item #**

RTL30100-CANE

- Coming Soon!
- Fits on the M+D by Drive Air-Go™ Walker with built-in accessory mount

[Shop Now](#)



## Cup Holder

**Item #**

RTL102662-CUP

- Fits on the M+D by Drive Air-Go™ Walker with built-in accessory mount
- Securely holds cups, water bottles, and travel mugs all at an arms reach

[Shop Now](#)



## NexStride® by Drive®

### Item #

FG-1000

- Audible Metronome: NexStride®'s metronome emits an audible signal, engaging a different pathway in your brain to enable you to move forward.
- Visual Laser Line: In addition to audible external signaling, NexStride® uses visual external signaling by way of a laser signal that shows you where to step, giving your brain a target to focus on and a new way to engage with the body (this can also help correct stride issues).
- Regain Your Independence: Use NexStride walking cane, rollator, or walker accessories in your home to navigate between rooms and use the bathroom unassisted. Take NexStride on walks through the neighborhood or to meet up with friends and continue to enjoy the freedom and independence that comes with walking on your own.

[Shop Now](#)



## WalkWise® by Drive®

### Item #

RTLWWA100

- Peace of Mind: Stay confidently connected to your loved one from anywhere.
- 24/7 Monitoring: See if they're active and consistently using their mobility device to help reduce fall risk.
- Real-Time Notifications: Get simple alerts if the device tips or goes inactive for extended periods (monthly subscription required).
- Unobtrusive & Easy: No cameras, wearables, charging, or smartphone needed, it runs automatically in the background.
- Activity Tracking: Built for walker and rollator users to capture distance and daily movement that fitness trackers miss.

[Shop Now](#)



## Universal 5" Walker Wheels

### Item #

10109

- Converts folding walker into wheeled walker
- Fits most manufacturers' walkers
- Allows for 8 height adjustments
- Polypropylene (PP) and PVC wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C) allowing use on all surfaces
- 350 lb Weight Limit

[Shop Now](#)

## Tall Extension Legs with Wheels, Combo Pack (Adds 4")

### Item #

10108WC

- Allows for height adjustments from 36" to 43" (walker handgrip to floor) (Adds 4" to height of walker)
- 5" rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- For use with all Drive and most leading manufacturers' walkers
- Includes rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- Weight Capacity 350 lb

[Shop Now](#)



## 5" Bariatric Walker Wheels with Two Sets of Rear Glides

**Item #**

10118SV

- Converts heavy-duty walker into wheeled walker
- Allows for 8 height adjustments
- 10118T Adds 4"
- Dual rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- For use with items # 10220-2 and 10223-2

[Shop Now](#)

## Swivel Wheel with Lock & Two Sets of Rear Glides

**Item #**

10115

- Can be locked for fixed wheel or unlocked for a swivel wheel
- Allows for 8 adjustments
- Rubber wheels allow walker to roll easily and smoothly over irregular surfaces
- Comes with rear glide caps, (item #10107) and glide covers (item #10107C), allowing use on all surfaces
- Weight Capacity 350 lb

[Shop Now](#)



## Rear Glide Brakes with 3" Wheels

**Item #**

10106

- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfort
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place
- 10105CB (Crutch Attachment) sold separately

[Shop Now](#)



## Platform Walker/Crutch Attachment

**Item #**

10105-1

**Description**

Single

10105-2

Case

- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfort
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place

[Shop Now](#)



## Bariatric Platform Walker/ Crutch Attachment

### Item #

10105HD-2

- Designed for those who are unable to grip the walker
- Cradles forearm with soft vinyl padding for comfortable weight bearing
- Large top bracket allows the platform attachment to be placed over the hand grip of the walker
- Adjustable strap securely holds arm in place
- Attachment for Underarm Crutch for item #10105 must be ordered separately (Cannot be used with Bariatric model)
- Bariatric Weight Capacity 500 lb

[Shop Now](#)

## Glide Ski

### Item #

### Description

10110	White
10110B	Black

- For use with all 1" folding walkers
- Easily added to walker
- Durable composite construction ensures long life
- Ideal on most surfaces
- For use with all Drive, M+D by Drive and most leading manufacturers' walkers

[Shop Now](#)



## Walker Glide Ski

**Item #**

RTL100017

- Allows walker to glide easily and smoothly over most surfaces
- Easily slides onto 1 1/8" tubing of walker
- Durable construction for long life
- Tool-free installation
- For use with most leading manufacturers' walkers

[Shop Now](#)



## Walker Glides

**Item #**

781-010

- For use with all 1" folding walkers
- Easily added to walker
- Durable composite construction ensures long life
- Ideal on most surfaces

[Shop Now](#)



## Tennis Ball Glides with Replaceable Glide Pads

Item #

10121

- Provides a quiet, smooth and durable glide experience when used with a walker
- Easy and safe to install
- Patent Pending

[Shop Now](#)

## Universal Walker Tray with Cup Holder

Item #

10124

- Allows personal items to be carried from room to room
- Contains one cup holder
- Fits most manufacturers' walkers
- Made of easy-to-clean durable plastic
- Easy-to-install
- Size: 23" (L) x 17" (W) x 1.5" (H)

[Shop Now](#)



## Walker Tray with Cup Holders

**Item #**  
10125

- Allows personal items to be carried from room to room
- Contains two cup holders
- Fits most manufacturers' 2-button walkers
- Easily folds down when not in use
- Made of easy-to-clean, durable plastic
- Easy-to-install
- Size: 16" (L) x 12" (W) x 2" (H)

[Shop Now](#)



## E-Z Walker Caddy

**Item #**  
RTL10131

- Allows personal items to be carried from room to room
- Includes a removable top that conceals personal items
- Fits most manufacturers' 2-button walkers
- Easily attaches with side clips
- Made of easy-to-clean durable plastic
- Easy-to-install
- Size: 15½" (L) x 12½" (W) x 5" (H)

[Shop Now](#)



## Walker Basket

### Item #

10200B

- For use with all 1" folding walkers
- 10200B comes with plastic insert tray with cup holder (10200B-1)

[Shop Now](#)

## Mobility Bags

### Item #

RTL10254BFL

### Description

Blue Floral

RTL10254GR

Army Green

RTL10254TFL

Tropical Floral

- Dual attachment options: back strap attachment with secure clip or adjustable hanging straps
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Four front exterior pockets and key ring zipper for added security of belongings
- Side mesh pocket holds water bottle
- Large, interior zip pocket
- Fully-lined pockets and interior
- Treated with water resistant technology
- Spot clean with wet cloth

[Shop Now](#)



## Multi-Use Accessory Bag

Item #	Description
RTL10255NV	Navy
RTL10255PK	Pink

- Back strap attachment with secure, adjustable clip
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Front pocket with magnetic closure for added security of belongings
- Reinforced back foundation maintains bag shape for added support
- Interior small pockets and zipper pocket
- Fully-lined pockets and interior
- Removable, optional cross-body carrying strap
- Treated with water resistant technology
- Spot clean with wet cloth

[Shop Now](#)



## Walker Accessory Bag

Item #
RTL10256NV

- Adjustable bag loop fasteners attach to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Three front pockets hold essentials such as keys, snacks or cell phone
- Two back pockets allow items to be accessed from user's side
- Lined interior and pockets
- Treated with water resistant technology
- Spot clean with wet cloth

[Shop Now](#)





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  [www.drivemedical.com](http://www.drivemedical.com)

CANADA OFFICE  256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7  866-752-1363

©2026 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.