

2026

Wheelchair & Transport Chair Catalog

Enhancing quality
of life is our
driving force



drive[®] 
MEDICAL


Wheelchair & Transport Chairs

Safe & Comfortable Mobility

Maintaining mobility is important when it comes to protecting independence and supporting wellbeing—but doing so safely and comfortably is also crucial. Drive Medical has put decades of clinical research and care into developing and building wheelchairs and transport chairs that put users and caregivers in the driver's seat. From easy-to-use wheel locks and anti-tippers, to padded seats, adjustable armrests, and elevated leg rests, our products are loaded with customizable features that both protect users and provide them with a comfortable, customizable ride. Our wheelchairs and transport chairs are also lightweight and fold compactly for easy storage and transport, to support and enhance the busy lives of users and those who care for them.

Table of Contents

Wheelchairs	3-19
Transport Chairs	20-25
Accessories	26-40



Wheelchairs



AeroFly™ Ultra-Lite Wheelchair

Ultra-lightweight yet ultra-capable, the AeroFly™ Ultra-Lite Wheelchair is engineered for users who live life on the move. With a sleek, modern design and a frame that weighs just 13.3 lb.*, it's among the lightest wheelchairs in its class. It's competitively priced and incredibly easy to ship, store, and deliver.

Item #	Description	UPC	UOM
AF18FDABK-SF	AeroFly Ultra-Lite Wheelchair, 18", Black	822383030135	1/ea
AF18FDAGY-SF	AeroFly Ultra-Lite Wheelchair, 18", Gray	822383030142	1/ea
AF20FDABK-SF	AeroFly Ultra-Lite Wheelchair, 20", Black	822383033204	1/ea
AF20FDAGY-SF	AeroFly Ultra-Lite Wheelchair, 20", Gray	822383033211	1/ea

Product Specs	18"	20"
Frame Material	Magnesium Alloy	
Frame Weight (w/o rear wheels and anti-tippers)	13.3 lb.	14 lb.
Product Weight (includes rear wheels)	19 lb.	20 lb.
Seat Dimensions	18" (W) x 16" (D)	20" (W) x 16" (D)
Seat-to-Floor Height	19"	
Weight Capacity	250 lb.	300 lb.
Warranty	Limited Lifetime	

Accessories



Carrying Bag
Item #AFCARRYBAG



Cup Holder
Item #STDSI040S



Gray Cushions, 18"
Item #AEROFly-06-18 (Back)
Item #AEROFly-28-18 (Seat)



Orange Cushions, 18"
Item #AEROFly-03-18 (Back)
Item #AEROFly-25-18 (Seat)

Effortless Patient Adoption

Offers stress-free portability your customers will love thanks to lightweight yet durable magnesium alloy construction, fold-down backrest, quick-release rear wheels, and built-in carry handles.

Fully Equipped to Keep Customers Safer

Safety comes standard with included wheeled anti-tippers, plus a dual braking system, treaded tires for improved traction and stability, and European-style footrests with integrated heel loops to keep feet secure.

Protects Your Margins

Competitively priced against similarly equipped ultra-lightweight wheelchairs to maximize your returns with MAP pricing to ensure you won't be undercut.

Stands Out to Sell More

Differentiate your offerings and draw customers to your showroom with the sleek and stylish patent-pending frame design in a choice of two premium finishes.

**Transport weight, 18" model – excludes rear wheels, anti-tippers, and footrests. Total product weight is 19 lbs with rear wheels attached.*

[Shop Now](#)

Wheelchairs

Nitro Cruise Wheelchair With Swing-Away Footrests



Item #	Seat Width	Color
NW16FBADARD-SF	16"	Red Frame
NW16FBADABK-SF	16"	Black Frame
NW18FBADARD-SF	18"	Red Frame
NW18FBADABK-SF	18"	Black Frame
NW20FBADARD-SF	20"	Red Frame
NW20FBADABK-SF	20"	Black Frame

- Modern, lightweight aluminum frame
- Lightweight spoke wheels with treaded tires for added traction and control with quick release rear axles
- Adjustable seat depth, height adjustable armrests, fold down backrest and angle adjustable footplates
- Hemi-Height adjustable
- Includes anti-tippers for added safety
- Includes removable, machine washable padded back cushion with removable foam adjustment blocks
- Includes mounting bracket for cup holder (102662-CUP), cup holder sold separately



Coming Soon

Nitro Cruise Wheelchair With Elevating Leg Rests



Item #	Seat Width	Color
NW16FBADARD-ELR	16"	Red Frame
NW16FBADABK-ELR	16"	Black Frame
NW18FBADARD-ELR	18"	Red Frame
NW18FBADABK-ELR	18"	Black Frame
NW20FBADARD-ELR	20"	Red Frame
NW20FBADABK-ELR	20"	Black Frame

- Modern, lightweight aluminum frame
- Lightweight spoke wheels with treaded tires for added traction and control with quick release rear axles
- Adjustable seat depth, height adjustable armrests, fold down backrest and angle adjustable footplates
- Hemi-Height adjustable
- Includes anti-tippers for added safety
- Includes removable, machine washable padded back cushion with removable foam adjustment blocks
- Includes mounting bracket for cup holder (102662-CUP), cup holder sold separately



Coming Soon

Wheelchairs

Rebel Wheelchair With Swing-Away Footrests

Item #	Seat Width	Description
RTLREB18DDA-SF	18"	Flip-Back Desk Arm

- Quick-release rear wheels allow for easy storage or transport
- Fold-down back and folding frame allows transport in almost any vehicle
- Easily removable swing-away footrests
- Tool-free removable front casters
- Desk-style armrests make for easy transfers and allows for comfortable seating at a table or desk
- Lightweight and extremely maneuverable
- Large accommodating seating surface

[Shop Now](#)



Blue Streak Wheelchair With Swing-Away Footrests

Item #	Seat Width	Description
BLS16FBD-SF	16"	Flip-Back Desk Arm
BLS18FBD-SF	18"	Flip-Back Desk Arm
BLS20FBD-SF	20"	Flip-Back Desk Arm

- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- solid rubber tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap
- Comes with swing-away footrests

[Shop Now](#)



Wheelchairs

Blue Streak Wheelchair With Elevating Leg Rests



Item #	Seat Width	Description
BLS16FBD-ELR	16"	Flip-Back Desk Arm
BLS18FBD-ELR	18"	Flip-Back Desk Arm
BLS20FBD-ELR	20"	Flip-Back Desk Arm

- Blue Powder-coated frame, black cross brace
- Flip-back, padded desk-length arm
- Nylon upholstery is durable, lightweight, attractive and easy to clean
- solid rubber tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Comes with calf strap
- Comes with elevating leg rests

[Shop Now](#)



Silver Sport 1 Wheelchair With Swing-Away Footrests



Item #	Seat Width	Description
SSP118FA-SF	18"	Full Arm

- Easy to clean powder coated silver vein steel frame
- Padded fixed armrests
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- polyurethane tires mounted on composite wheels provide durability and low maintenance
- Comes with push-to-lock wheel locks
- Carry pocket on backrest provides additional convenience
- Comes with swing-away footrests
- Folding back-canes for transport and storage

[Shop Now](#)



Wheelchairs

Silver Sport 2 Wheelchair With Swing-Away Footrests



Item #	Seat Width	Description
SSP216FA-SF	16"	Full Arm
SSP216DDA-SF	16"	Detachable Desk Arm
SSP216DFA-SF	16"	Detachable Full Arm
SSP218FA-SF	18"	Full Arm
SSP218DDA-SF	18"	Detachable Desk Arm
SSP218DFA-SF	18"	Detachable Full Arm
SSP220DDA-SF	20"	Detachable Desk Arm
SSP220DFA-SF	20"	Detachable Full Arm

- Powder-coated silver vein steel frame
- Durable, lightweight, and easy-to-clean vinyl upholstery
- Polyurethane tires mounted on composite wheels
- Dual axle transitions seat height from 19.5"-17.5"
- 8" front casters are adjustable in three positions
- Padded armrests for comfort
- Carry pocket on backrest provides storage
- Footrests include plastic foot plates with heel loops
- Push-to-lock wheel locks



[Shop Now](#)

Silver Sport 2 Wheelchair With Elevating Leg Rests



Item #	Seat Width	Description
SSP216FA-ELR	16"	Full Arm
SSP216DDA-ELR	16"	Detachable Desk Arm
SSP216DFA-ELR	16"	Dedicated Full Arm
SSP218FA-ELR	18"	Full Arm
SSP218DDA-ELR	18"	Detachable Desk Arm
SSP218DFA-ELR	18"	Detachable Full Arm
SSP220DDA-ELR	20"	Detachable Desk Arm
SSP220DFA-ELR	20"	Detachable Full Arm

- Powder-coated silver vein steel frame
- Durable, lightweight, and easy-to-clean vinyl upholstery
- Polyurethane tires mounted on composite wheels
- Dual axle transitions seat height from 19.5"-17.5"
- 8" front casters are adjustable in three positions
- Padded armrests for comfort
- Carry pocket on backrest provides storage
- Footrests include plastic foot plates with heel loops
- Push-to-lock wheel locks



[Shop Now](#)



Silver Sport Recliner with Vinyl Upholstery and Elevating Leg Rests

Item #	Seat Width	Description
SSP16RBDFAV	16"	Detachable Full Arm
SSP16RBDDA	16"	Detachable Desk Arm
SSP18RBDFAV	18"	Detachable Full Arm
SSP18RBDDAV	18"	Detachable Desk Arm
SSP20RBDFAV	20"	Detachable Full Arm
SSP20RBDDAV	20"	Detachable Desk Arm

- Vinyl upholstery is comfortable and easy-to-clean
- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Arms include support extension for reclining position
- Carbon Steel frame with silver vein finish
- Composite Mag-style wheels with 8" front casters
- Padded armrests provide added patient comfort
- Wheels are set back on frame to prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Headrest extension with cushioned head immobilizer Comes with push to lock wheel locks

[Shop Now](#)





Cruiser III Wheelchair With Swing-Away Footrests

Item #	Seat Width	Description
K316DDA-SF	16"	Detachable Desk Arm
K316DFA-SF	16"	Detachable Full Arm
K316ADDA-SF	16"	Flip Back, Detachable, Adj. Desk Arm
K316ADFA-SF	16"	Flip Back, Detachable, Adj. Full Arm
K318DDA-SF	18"	Detachable Desk Arm
K318DFA-SF	18"	Detachable Full Arm
K318ADDA-SF	18"	Flip Back, Detachable, Adj. Desk Arm
K318ADFA-SF	18"	Flip Back, Detachable, Adj. Full Arm
K320DDA-SF	20"	Detachable Desk Arm
K320DFA-SF	20"	Detachable Full Arm
K320ADDA-SF	20"	Flip Back, Detachable, Adj. Desk Arm
K320ADFA-SF	20"	Flip Back, Detachable, Adj. Full Arm

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18"
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Weighs under 36 lb (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer
- Frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- Dual axle provides easy transition of seat height to hemi-level
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

[Shop Now](#)





Cruiser III Wheelchair With Elevating Legrests

Item #	Seat Width	Description
K316DDA-ELR	16"	Detachable Desk Arm
K316DFA-ELR	16"	Detachable Full Arm
K316ADDA-ELR	16"	Flip Back, Detachable, Adj. Desk Arm
K316ADFA-ELR	16"	Flip Back, Detachable, Adj. Full Arm
K318DDA-ELR	18"	Detachable Desk Arm
K318DFA-ELR	18"	Detachable Full Arm
K318ADDA-ELR	18"	Flip Back, Detachable, Adj. Desk Arm
K318ADFA-ELR	18"	Flip Back, Detachable, Adj. Full Arm
K320DDA-ELR	20"	Detachable Desk Arm
K320DFA-ELR	20"	Detachable Full Arm
K320ADDA-ELR	20"	Flip Back, Detachable, Adj. Desk Arm
K320ADFA-ELR	20"	Flip Back, Detachable, Adj. Full Arm

- Built-in seat rail extensions and extendable upholstery easily adjust seat depth from 16" to 18"
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Weighs under 36 lb (excluding front riggings)
- Carbon steel frame with silver vein finish
- Removable flip-back arms allow for easy transfer
- Frame style eliminates seat guides and allows for custom back inserts and accessories
- Nylon upholstery is durable, lightweight, attractive and easy-to-clean
- Composite, Mag-style wheels are lightweight and maintenance free
- Dual axle provides easy transition of seat height to hemi-level
- 8" front casters are adjustable in three positions
- Padded armrests provide additional comfort
- Tool-free adjustable length riggings
- Comes with push-to-lock wheel locks

[Shop Now](#)





Cruiser X4 Wheelchair With Swing-Away Footrests

Item #	Seat Width	Description
CX416ADDA-SF	16"	Flip Back, Detachable, Adj. Height Desk Arm
CX416ADFA-SF	16"	Flip Back, Detachable, Adj. Height Full Arm
CX418ADDA-SF	18"	Flip Back, Detachable, Adj. Height Desk Arm
CX418ADFA-SF	18"	Flip Back, Detachable, Adj. Height Full Arm
CX420ADDA-SF	20"	Flip Back, Detachable, Adj. Height Desk Arm
CX420ADFA-SF	20"	Flip Back, Detachable, Adj. Height Full Arm

- Weighs 34 lbs or less (excluding front riggings).
- Height-adjustable, flip-back, and removable armrests
- Back height is adjustable.
- Aluminum and carbon steel frame
- Flame-retardant nylon upholstery resists mildew and bacteria.
- Seat depth adjustable from 16" to 18".
- Lightweight, maintenance-free mag-style wheels with quick-release axles.
- Dual over-center support for added frame rigidity.
- Swing-away footrests with push-to-lock wheel locks.
- Durable, lightweight composite footplates
- Adjustable 8" front casters
- Easy transition of seat height

[Shop Now](#)



Wheelchairs



Cruiser X4 Wheelchair With Elevating Legrests

Item #	Seat Width	Description
CX416ADDA-ELR	16"	Flip Back, Detachable, Adj. Height Desk Arm
CX416ADFA-ELR	16"	Flip Back, Detachable, Adj. Height Full Arm
CX418ADDA-ELR	18"	Flip Back, Detachable, Adj. Height Desk Arm
CX418ADFA-ELR	18"	Flip Back, Detachable, Adj. Height Full Arm
CX420ADDA-ELR	20"	Flip Back, Detachable, Adj. Height Desk Arm
CX420ADFA-ELR	20"	Flip Back, Detachable, Adj. Height Full Arm

- Weighs 34 lbs or less (excluding front riggings).
- Height-adjustable, flip-back, and removable armrests
- Back height is adjustable.
- Aluminum and carbon steel frame
- Flame-retardant nylon upholstery resists mildew and bacteria.
- Seat depth adjustable from 16" to 18".
- Lightweight, maintenance-free mag-style wheels with quick-release axles.
- Dual over-center support for added frame rigidity.
- Swing-away footrests with push-to-lock wheel locks.
- Durable, lightweight composite footplates
- Adjustable 8" front casters
- Easy transition of seat height

[Shop Now](#)



Wheelchairs

Viper Plus GT Wheelchair with Swing-Away Footrests



Item #	Seat Width	Description
PLA416FBUARAD-SF	16"	Universal Armrests
PLA418FBUARAD-SF	18"	Universal Armrests
PLA420FBUARAD-SF	20"	Universal Armrests
PLA422FBUARAD-SF	22"	Universal Armrests

- Convert the armrests into desk arm, half arm, and full arm
- Adjustable flip-back, removable, padded arms
- Easily adjust seat depth from 16" to 18"
- Dual Axle allows for multiple height positions
- Crace eliminates need for seat guides
- Adjustable angle back easily adjusts from 5° – 20°
- Durable, heavy-gauge flame-retardant nylon upholstery
- Quick release, 24" composite Mag-style wheels
- Adjustable angle caster forks have three height adjustments and angle adjustability
- Adjustable anti-tipper with flip-up wheels
- Comes with push-to-lock wheel locks



[Shop Now](#)

Viper Plus GT Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
PLA416FBUARAD-ELR	16"	Universal Armrests
PLA418FBUARAD-ELR	18"	Universal Armrests
PLA420FBUARAD-ELR	20"	Universal Armrests
PLA422FBUARAD-ELR	22"	Universal Armrests

- Convert the armrests into desk arm, half arm, and full arm
- Adjustable flip-back, removable, padded arms
- Easily adjust seat depth from 16" to 18"
- Dual Axle allows for multiple height positions
- Crace eliminates need for seat guides
- Adjustable angle back easily adjusts from 5° – 20°
- Durable, heavy-gauge flame-retardant nylon upholstery
- Quick release, 24" composite Mag-style wheels
- Adjustable angle caster forks have three height adjustments and angle adjustability
- Adjustable anti-tipper with flip-up wheels
- Comes with push-to-lock wheel locks



[Shop Now](#)

Wheelchairs

Viper Plus Reclining Wheelchair



Item #	Seat Width	Description
PLA416RBDDA	16"	Detachable Desk Arm
PLA416RBDFA	16"	Detachable Full Arm
PLA418RBDDA	18"	Detachable Desk Arm
PLA418RBDFA	18"	Detachable Full Arm
PLA420RBDDA	20"	Detachable Desk Arm
PLA420RBDFA	20"	Detachable Full Arm

- Easily adjusts seat depth from 16" or 18"
- Attendant wheel locks
- Spreader bar acts as a back stabilizer and a convenient push handle for the attendant
- Offset rear wheel position adds additional stability and equally distributes the patient's weight
- Hydraulic reclining mechanism adjusts up to 180°
- Durable, flame-retardant nylon upholstery
- Carbon steel frame with double-coated chrome
- Flip-back style padded armrests for easy transfers
- Easy adjustment of seat-to-floor height to hemi-level
- Front caster forks are adjustable in three positions

[Shop Now](#)



Pediatric Viper Wheelchair with Swing-Away Footrests



Item #	Seat Width	Description
L412DDA-SF	12"	Detachable Desk Arm
L412DDA-SF	14"	Detachable Desk Arm

- Frame style eliminates seat rail guides Back height adjustable from 18"-20"
- Carbon steel black powder-coated frame
- Durable, flame retardant nylon upholstery
- Composite Mag-style wheels with aluminum rims
- Polyurethane rear tires and 8" front casters
- Nylon, padded armrests provide added user comfort
- Flip-back arms allow for easy transfers
- The front mounted release lever can be flipped up or down
- Precision sealed wheel bearings in front and rear
- Dual axle provides easy adjustment of seat height
- Composite foot plates with heel loops
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks

[Shop Now](#)



Wheelchairs

Pediatric Viper Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
L412DDA-ELR	12"	Detachable Desk Arm
L414DDA-ELR	14"	Detachable Desk Arm

- Frame style eliminates seat rail guides Back height adjustable from 18"-20"
- Carbon steel black powder-coated frame
- Durable, flame retardant nylon upholstery
- Composite Mag-style wheels with aluminum rims
- Polyurethane rear tires and 8" front casters
- Nylon, padded armrests provide added user comfort
- Flip-back arms allow for easy transfers
- The front mounted release lever can be flipped up or down Precision sealed wheel bearings in front and rear
- Dual axle provides easy adjustment of seat height
- Composite foot plates with heel loops
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks



[Shop Now](#)

Pediatric Viper Plus Reclining Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
PL412RBDDA	12"	Detachable Desk Arm
PL414RBDDA	14"	Detachable Desk Arm

- State-of-the art hydraulic reclining mechanism
- Backrest reclines from 90° to 180°
- Double-coated chrome carbon steel frame
- Durable flame retardant nylon upholstery
- Offset rear wheel position provides additional stability and evenly distributes user's weight
- Precision sealed front and rear wheel bearings
- Front caster forks are adjustable in three positions
- Easy transition of seat to hemi-level
- Push-to-lock wheel locks
- Spreader bar functions as a back stabilizer and a push handle for caregiver
- Flip-up anti-tippers with wheels



[Shop Now](#)

Wheelchairs

Lynx Ultra Lightweight Wheelchair with Swing-Away Footrests



Item #	Seat Width	Description
K516FBADDA-SF	16"	Detachable Desk Arm
K518FBADDA-SF	18"	Detachable Desk Arm
K520FBADDA-SF	20"	Detachable Desk Arm

- Ultra Lightweight Aluminum Frame
- Flip-back, padded, height adjustable desk arms provides added patient comfort
- Comes with push to lock wheel locks
- Treaded tires
- Quick release axles
- Dual axle positions
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Adjustable seat depth from 16" - 18"
- Rear anti-tippers come standard



[Shop Now](#)

Lynx Ultra Lightweight Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
K516FBADDA-ELR	16"	Detachable Desk Arm
K518FBADDA-ELR	18"	Detachable Desk Arm
K520FBADDA-ELR	20"	Detachable Desk Arm

- Ultra Lightweight Aluminum Frame
- Flip-back, padded, height adjustable desk arms provides added patient comfort
- Comes with push to lock wheel locks
- Treaded tires
- Quick release axles
- Dual axle positions
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability
- Adjustable seat depth from 16" - 18"
- Rear anti-tippers come standard



[Shop Now](#)

Wheelchairs

Deluxe Sentra Full Reclining Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
STD22RBDDA	22"	Detachable Desk Arm

- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Durable, heavy-gauge naugahyde upholstery
- Carbon Steel frame with triple-coated chrome
- Composite Mag-style wheels with chrome hand rims
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels set back on frame prevent tipping
- Precision sealed wheel bearings in front and rear
- Rear anti-tippers standard
- Front caster forks are adjustable in two positions
- Carry pocket standard
- Headrest extension with cushioned head immobilizer
- Comes with push-to-lock wheel locks



[Shop Now](#)

Bariatric Sentra EC Heavy-Duty Wheelchair with Swing-Away Footrests



Item #	Seat Width	Description
STD20ECDDAHD-SF	20"	Detachable Desk Arm
STD20ECDFAHHD-SF	20"	Detachable Full Arm
STD22ECDDAHD-SF	22"	Detachable Desk Arm
STD22ECDFAHHD-SF	22"	Detachable Full Arm
STD24ECDDAHD-SF	24"	Detachable Desk Arm
STD24ECDFAHHD-SF	24"	Detachable Full Arm

- Carbon Steel frame with powder-coated silver vein
- Double embossed vinyl upholstery
- Easy transition of seat height to hemi-level
- Polyurethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in three positions
- Padded armrests provide additional comfort
- Tool-free adjustable length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks



[Shop Now](#)

Wheelchairs

Bariatric Sentra EC Heavy-Duty Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
STD20ECDDAHD-ELR	20"	Detachable Desk Arm
STD20ECDFAHd-ELR	20"	Detachable Full Arm
STD22ECDDAHD-ELR	22"	Detachable Desk Arm
STD22ECDFAHd-ELR	22"	Detachable Full Arm
STD24ECDDAHD-ELR	24"	Detachable Desk Arm
STD24ECDFAHd-ELR	24"	Detachable Full Arm

- Carbon Steel frame with powder-coated silver vein
- Double embossed vinyl upholstery
- Easy transition of seat height to hemi-level
- Polyurethane tires mounted on composite wheels provide durability and low maintenance
- 8" front casters, adjustable in three positions
- Padded armrests provide additional comfort
- Tool-free adjustable length riggings
- Aluminum foot plates
- Comes with push-to-lock wheel locks



[Shop Now](#)

Bariatric Sentra Extra Heavy-Duty Wheelchair with Swing-Away Footrests



Item #	Seat Width	Description
STD20DDA-SF	20"	Detachable Desk Arm
STD22DDA-SF	22"	Detachable Desk Arm
STD24DDA-SF	24"	Detachable Desk Arm

- Carbon Steel frame with silver vein finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Dual-position front fork
- Dual-position front frame
- Precision sealed wheel bearings in front and rear Padded armrests provide added user comfort
- Tool-free adjustable length riggings
- Comes with push-to-lock wheel locks



[Shop Now](#)

Wheelchairs

Bariatric Sentra Extra Heavy-Duty Wheelchair with Elevating Leg Rests



Item #	Seat Width	Description
STD20DDA-ELR	20"	Detachable Desk Arm
STD22DDA-ELR	22"	Detachable Desk Arm
STD24DDA-ELR	24"	Detachable Desk Arm

- Carbon Steel frame with silver vein finish
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable, heavy-gauge reinforced naugahyde upholstery
- Comes standard with carry pocket on backrest
- Dual-position front fork
- Dual-position front frame
- Precision sealed wheel bearings in front and rear padded armrests provide added user comfort
- Tool-free adjustable length riggings
- Comes with push-to-lock wheel locks



[Shop Now](#)

Bariatric Sentra EC Heavy-Duty, Extra-Extra-Wide Wheelchair



Item #	Seat Width	Description
STD26ECDDA	26"	Detachable Desk Arm
STD26ECDFA	26"	Detachable Full Arm
STD28ECDDA	28"	Detachable Desk Arm
STD28ECDFA	28"	Detachable Full Arm
STD30ECDDA	30"	Detachable Desk Arm
STD30ECDFA	30"	Detachable Full Arm

- Carbon steel frame with triple-coated chrome
- Reinforced steel gussets at all weight bearing points provide additional strength
- Durable reinforced nylon upholstery
- Steel, spoke wheels are strong and maintenance free
- Heavy-duty, flat-free 8" front casters are adjustable
- Extra-heavy-duty front forks
- Reinforced side frame, caster journals and dual cross bars
- Precision sealed wheel bearings in front and rear
- Comes with push-to-lock wheel locks



[Shop Now](#)

Transport Chairs

AeroFly™ Ultra-Lite Transport Chair



Item #	Seat Width	Color
AFT18GY-SF	18"	Grey Frame & Blue Stitching
AFT18BK-SF	18"	Black Frame & Red Stitching

- Ultralight magnesium alloy frame weighs just 13.3 lbs. (excluding rear wheels, 17 lbs. with rear wheels)
- Includes attendant brakes
- Fixed, desk length arms with molded armrest pads for comfort and durability
- Lightweight spoke wheels with treaded flat-free tires
- Large 13" treaded rear wheels for a smooth ride over many surfaces
- Quick release rear axles
- Fold down backrest
- Seat integrated carry handles
- Removable seat and back padded overlay cushions
- Swing away footrests with ergonomic release levers



Coming Soon

Poly-Fly High Strength, Lightweight Wheelchair/ Flyweight Transport Chair Combo



Item #	Seat Width	Color
TR18	18"	White Frame & Black Upholstery
TR20	20"	

- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair
- All-aluminum frame
- Comes with two sets of installed wheel locks
- Back folds down with deluxe back release
- Comes standard with swing-away footrests with composite foot plates and heel loops
- 8" casters and anti-tippers included
- Compatible with optional Swing-away elevated leg rests (Material# STDELR-TF)



Shop Now

Transport Chairs

Steel Transport Chair



Item #	Seat Width	Color
TR37E-SV	17"	Silver Vein Frame & Black Upholstery
TR39E-SV	19"	Silver Vein Frame & Black Upholstery

- Weighs 26 lb
- Attractive silver vein finish
- Available in 17" or 19" seat widths
- Back folds down for storage and transport
- Swing-away footrests come standard
- Padded armrests provide additional comfort
- 8" casters in front and rear
- Comes with attractive nylon upholstery that is durable and easy-to-clean
- Standard rear wheel locks
- Comes with upholstery handles to fold chair



[Shop Now](#)

Aluminum Transport Chair



Item #	Seat Width	Color
ATC17-BK	17"	Black Frame & Black Upholstery
ATC17-BL	17"	Blue Frame & Black Upholstery
ATC17-RD	17"	Red Frame & Black Upholstery
ATC19-BK	19"	Black Frame & Black Upholstery
ATC19-BL	19"	Blue Frame & Black Upholstery
ATC19-RD	19"	Red Frame & Black Upholstery

- Weighs only 19 lb
- Aluminum frame is lightweight and durable
- Composite 8" wheels are lightweight and maintenance free
- Come with swing-away footrests
- Seat belt included for added safety
- Padded armrests provide additional comfort
- Nylon upholstery is comfortable, lightweight and easy-to-clean
- Aluminum casted folding hinge allows back to fold down for easy storage and transport



[Shop Now](#)

Transport Chairs

Lightweight Expedition Aluminum Transport Chair



Item #	Seat Width	Color
EXP19LTBL	19"	Blue Frame & Black Upholstery
EXP19LTRD	19"	Red Frame & Black Upholstery

- Weighs only 19 lbs
- Available in blue or red frames
- Folds flat for easy transport
- Handy lever on chair back makes it easy and convenient to fold back down
- Comes with seat belt
- Companion activated wheel locks
- 12" rear "flat-free" wheels



[Shop Now](#)

Fly-Lite Aluminum Transport Chair



Item #	Seat Width	Color
DFL19-BL	19"	Blue Frame & Blue Plaid Upholstery
DFL19-BK	19"	Black Frame & Red Plaid Upholstery
DFL19-RD	19"	Burgandy Frame & Black Upholstery

- Weighs only 16.2 lb without footrests
- Aluminum frame is lightweight and durable
- Comes with carry pocket on backrest
- Includes seat belt for added safety
- Attractive plaid style nylon upholstery
- Comes with handy carry pocket on back
- Tool-free, height-adjustable swing-away footrests 8" rear wheels and 6" front casters
- Folding back canes allows backrest to fold down for storage and transport
- Comes with aluminum rear wheel locks
- Optional elevating leg rests available (Model# STDELR-TF)



[Shop Now](#)

Transport Chairs

Deluxe Fly-Weight Aluminum Transport Chair with Removable Casters



Item #	Seat Width	Color
RTLFW19RW-RD	19"	Red Frame & Black Upholstery

- Weighs only 19 lbs
- Aluminum frame for lightweight strength and durability
- Easy folding and storage in closet or trunk
- Composite wheels are lightweight and maintenance free
- 8" wheels in front and rear easily remove without tools
- Deluxe back release
- Deluxe, all-aluminum rear wheel locks
- Comes with carry pocket on backrest
- Tool-free, swing-away footrests included
- Seat belt included for added safety
- Flared arms allow the back to be folded virtually flat for easy storage and transport
- Soft, plush nylon upholstery

[Shop Now](#)



TranSport Aluminum Transport Chair



Item #	Seat Width	Color
TS19	19"	Blue Frame & Blue Upholstery

- Weighs 17.5 lbs without footrests
- Aluminum frame is lightweight and strong
- Frame magnets keep footrests secure against frame to prevent swinging open when the chair is carried
- Padded flip-back arms are removable
- Upholstery is comfortable, lightweight and easy-to-clean
- Two-way back release allows to push or pull lever
- The transport chair conveniently folds for travel and storage
- Back folds forward
- Comes with tool-free height adjustable swing-away footrests
- Aluminum rear wheel locks

[Shop Now](#)



Transport Chairs

Super Light, Folding Transport Chair

Item #	Seat Width	Color
SL18	18"	Silver Frame & Blue Upholstery

- Padded flip-back armrests
- Fold-down footrest
- Weighs only 19 lb
- Easy to push or transport
- Aluminum frame with lightweight nylon upholstery
- Composite wheels are lightweight and maintenance free
- Folds like a conventional folding chair
- Carry bag (item #835N) with shoulder strap included
- Standard rear wheel locks and seat belt

[Shop Now](#)



Travelite Transport Chair

Item #	Seat Width	Color
TC005GY	18"	Black Frame & Black Upholstery

- The Travelite transport chair is lightweight and sturdy aluminum transport chair
- It is easily transportable and comes with a compact carry bag with wheels making it ideal for traveling
- Folds for storage and transportation
- It is ideal for use indoors and outdoors
- 18" square nylon seat comes with lift-up straps
- Comes with wheel locks for 8" solid tires
- Comes with padded, flip-up desk arms and height adjustable footrests

[Shop Now](#)



Transport Chairs

Bariatric Steel Transport Chair



Item #	Seat Width	Color
BTR20-B	20"	Blue Frame & Black Upholstery
BTR20-R	20"	Red Frame & Black Upholstery
BTR22-B	22"	Blue Frame & Black Upholstery
BTR22-R	22"	Red Frame & Black Upholstery

- Reinforced steel frame provides added support
- Supports individuals up to 450 lb
- Dual, reinforced steel cross braces
- 12" rear flat free wheels
- Heavy-duty, nylon reinforced upholstery with a back carry pouch
- Removable, reversible desk-length arms
- Padded armrests provide additional comfort
- Companion activated wheel locks
- Available in blue or red frames
- Swing-away footrests included

[Shop Now](#)



22" Bariatric Aluminum Transport Chair



Item #	Seat Width	Color
ATC22-R	22"	Red Frame & Black Upholstery

- Weighs only 33 lbs
- Supports individuals up to 450 lbs
- Attachable push bar with padded foam grip
- 12" rear flat-free wheels
- Back folds down for convenient storage
- Padded armrests provide additional comfort
- Heavy-duty, nylon reinforced upholstery
- Comes with hand wheel lock levers and companion activated wheel locks
- Comes with seat belt
- Comes with swing-away footrests and anti-tippers

[Shop Now](#)



Accessories

Drive® DualEase® Gel Cushion

Item #	Product Description
14991	16" W x 16" L x 2" H
14992	18" W x 16" L x 2" H
14993	20" W x 16" L x 2" H

- Expertly engineered for pressure relief, skin protection, and personalized comfort. Perfect for individuals with limited mobility or anyone who values superior seated comfort.
- Dual-layer unique polymer gel structure features a different firmness level on each side - simply flip the cushion to customize your comfort and relieve pressure.
- Flexible gel structure reduces shear and friction, supporting those with limited seated tolerance or at risk for pressure ulcers.
- Open-cell design and cooling hyper-elastic gel promote airflow and prevent moisture, consistently providing a cool seating experience.
- Built-in carry handle makes transport easy, while the non-skid bottom ensures the cushion stays put during use.

[Shop Now](#)



Gel-U-Seat™ Skin Protection 3" Gel/Foam Wheelchair Cushion

Item #	Product Description
8046-1	16" W x 16" L x 3" H
8046-2	18" W x 16" L x 3" H
8046-3	20" W x 16" L x 3" H
8046-4	20" W x 18" L x 3" H
8046-5	22" W x 18" L x 3" H
8046-7	24" W x 18" L x 3" H

- The Gel-U-Seat™ skin protection cushion has a high density foam shell surrounding a uniquely designed gel bladder
- Dual chambered gel bladder provides proper immersion and envelopment for pressure redistribution
- Cover is water resistant and low shear stretch to reduce heat and moisture build-up
- Convenient 12" tie straps secure the cushion to the wheelchair
- Meets flammability standard TB117

[Shop Now](#)



Gel "E" Skin Protection 3" Gel/ Foam Wheelchair Cushion

Item #	Product Description
14891	16" W x 16" L x 2" H
14886	18" W x 16" L x 3" H
14878	18" W x 18" L x 3" H
14894	20" W x 16" L x 3" H
14902	22" W x 18" L x 3" H

- Pressure redistributing cushion designed to assist in the prevention, treatment and management of pressure injuries
- Viscous gel bladder is encased in polyurethane foam shell and provides optimal pressure redistribution, support and comfort
- Cover top is a urethane-coated nylon that is low-shear, vapor permeable and water-resistant
- Cover base is waterproof vinyl for durability
- Carry handle included for convenient transport
- Meets flammability standard TB117

[Shop Now](#)




Titanium Skin Protection and Positioning Gel/Foam Wheelchair Cushion

Item #	Product Description
FPT-1	16" W x 16" L x 3.5" H
FPT-2	18" W x 16" L x 3.5" H
FPT-1818	18" W x 18" L x 3.5" H
FPT-2018	20" W x 18" L x 3.5" H
FPT-4	22" W x 18" L x 3.5" H

- The Titanium cushion is specifically designed to offer superior pressure redistribution and positioning
- High-resilient, high-density foam combined with a dual chamber gel bladder encourages proper patient posture through immersion
- Middle contour acts as a built-in leg abductor
- Low-shear, fluid-resistant stretch cover helps to maintain sensitive skin tissue
- Non-skid bottom offers additional stability
- Meets flammability standard TB117

[Shop Now](#)



Molded Foam General Use Wheelchair Cushion

Item #	Product Description
14880	16" W x 16" L x 2" H
14887	18" W x 16" L x 2" H
14881	20" W x 16" L x 2" H
14908	18" W x 18" L x 2" H
14909	20" W x 18" L x 2" H

- A unique construction of high-resilient molded foam provides optimal support, stabilization and pressure redistribution
- Designed for low risk patients seeking comfort and support
- Removable, zippered, low-shear stretch cover is vapor permeable and water-resistant
- Cover has a "Non-Slip" vinyl bottom for comfort and safety
- Convenient 12" tie straps secure the cushion to wheelchair
- Meets flammability standard TB117

[Shop Now](#)

Adjustable Tension General Use Wheelchair Back Cushion

Item #	Product Description
14300	16"-21" Adjustable
14301	22"-26" Adjustable

- Versatile back support system that improves patient comfort and posture
- 1" thick high density foam pad provides optimal pressure redistribution
- Secures directly over wheelchair canes for added strength, stability and safety
- 5 Adjustment straps featured for multi-purpose back support and customization to individual seating needs
- Handy back pocket is convenient for carrying keys, water bottle, medications, tissues and more
- Breathable cover is fluid resistant and easy to clean
- Cushion folds with the wheelchair for easy transport and storage
- Meets flammability standard TB117

[Shop Now](#)



Seat Rail Extension Kit, 16" Wide

Item # STDSCS1618V

- Kit comes complete with seat upholstery and attached rail extension to add 2" to depth of the wheelchair seat Kit installation results in 18" seat depth For replacement of 16" wide, 16" deep upholstery

[Shop Now](#)



Seat Rail Extension Kit, 18" Wide

Item # STDSCS1818V

- Kit comes complete with seat upholstery and attached rail extension to add 2" to depth of wheelchair seat For replacement of 18" wide, 16" deep upholstery

[Shop Now](#)



Seat Rail Extension Kit, 20" Wide

Item # STDSCS2018V

- Kit comes complete with seat upholstery and attached rail extension to add 2" to depth of the wheelchair seat. Kit installation results in 20" seat depth. For replacement of 16" wide, 16" deep upholstery.

[Shop Now](#)



Seat Belt

Item

STDS850
STDS851
STDS855
STDS856

Description

Standard, Auto-Clasp
Standard, Hook and Loop
Bariatric, Auto-Clasp
Bariatric, Hook and Loop

- Provides extra security and safety while in wheelchair
- Easy to install

[Shop Now](#)



Universal Forward Stabilizers

Item # STDS860

- Can accommodate a secure fit with Sentra EC 20", Sentra HD 20", Sentra Reclining, Cruiser III, Cirrus IV, Viper, Viper Plus GT, Viper Plus Reclining and most other leading manufacturers wheelchairs.
- Easy to install

[Shop Now](#)

Amputee Adapters

Item # STDS828GT

- Allows the rear wheels of a wheelchair to be moved further back, adjusting the center of gravity
- Easy to install

[Shop Now](#)



Anti-Tippers

Item #

STDS802
STDS807
STDS814
STDS819
STDS829
STDS833N
STDS818

For Use With:

Sentra EC (20"-24"), Sentra HD
Cruiser III, Cruiser X4, Viper Plus GT
Viper Jr
Cruiser III, Cruiser X4
Sentra EC Extra-Wide (26"-30")
Blue Streak, Silver Sport 1 & 2, Polyfly
Silver Sport 2, Silver Sport Reclining,
Sentra Reclining

- Provides additional stability and help prevent a wheelchair from tipping backwards
- Easy to install

[Shop Now](#)



Swing-away Elevating Legrests

Item #

ALELR
STDELR-TF
HDELR
PHELRL
LK3JELR

For Use With:

Viper Plus GT
Sentra EC, Sentra HD, Sentra Reclining
Sentra HD
Sentra EC Extra-Wide (26"-30")
Blue Streak, Silver Sport 1, Silver Sport 2,
Cruiser III, Cruiser X4, Polyfly

- Allows user to elevate legs while seated in wheelchair
- Able to be swung away from chair when getting in or out of wheelchair
- Easy to install

[Shop Now](#)



Swing-Away Footrests

Item #

STDSF-TF

PH-SF

AL3-4SF

STDS3J24SF

For Use With:

Sentra EC (20"-24"), Sentra HD

Sentra EC Extra-Wide (26"-30")

Viper Plus GT

Blue Streak, Silver Sport 1, Silver Sport 2, Cruiser III, Cruiser X4, Polyfly

- Provides foot support while seated in wheelchair
- Able to be swung away from chair when getting in or out of wheelchair
- Easy to install

[Shop Now](#)

Articulating Elevating Legrests

Item # STDELR-AL

- Extend as they are elevated
- Compatible with most Drive Wheelchairs
- Easy to install

[Shop Now](#)



Wheelchair Cane/Crutch Holder

Item # STDS1034

- Securely holds cane or crutch while using wheelchair
- Easy to install

[Shop Now](#)



Anti-Fold Device

Item #

STDS806

STDS806-HD

Seat Size:

16" - 22"

24" - 30"

- Installs across bottom of chair to prevent folding

[Shop Now](#)



Universal Telescoping I.V. Pole

Item # STDS820

- Universal spacing sleeve can be added or removed to accommodate most Drive (except Sentra Recliner, Viper GT Recliner, Chrome Sport Recliner, Silver Sport Recliner, Viper Plus GT, Cougar, Cruiser X4, Sentra HD Extra Wide) and other leading manufacturers' wheelchairs

[Shop Now](#)

Overhead Adjustable Anti-Theft Device

Item # STDS834

- Prevents wheelchair from being loaded into a small space, such as a car trunk or back seat
- Easy to install

[Shop Now](#)



Oxygen "E" Cylinder Holders

Item # STDS803

- Holds "E" Cylinder Tank
- Protruding Bar holds Oxygen Cylinder in place
- Universal spacing sleeve can be added or removed to accommodate a secure fit with additional models

[Shop Now](#)

Flat-Free Wheel with Handrim

Item # STDS840

- For use with Sentra Reclining, Viper Reclining, Viper and Cirrus IV Wheelchairs

[Shop Now](#)



One-Arm Drive

Item # STDS4056

- Allows wheelchair to be propelled with one arm

[Shop Now](#)

Super-Hemi Kit

Item #

STDSSHK

STDSSHK-K3

For Use With:

Viper Plus GT

Cruiser III

- Lowers the seat to floor height by 2"

[Shop Now](#)



Limb Support

Item #

WASL

WASR

Production Description:

Limb Support, Left

Limb Support, Right

- Provides limb support for amputees
- Easy to install

[Shop Now](#)

Wheel Lock Extensions

Item #

STDS601

STDS801

For Use With:

Viper Plus GT, Lynx, Cirrus IV, Cougar

Blue Streak, Entra EC, Sentra HD,
Sentra Reclining, Sentra EC Extra-Wide,
Cougar, Cirrus IV

STDS801GT

Cruiser III, Revel, Silver Sport 1 & 2
(Excluding serial #3J)

- Extends wheel lock handles for easier access and use
- Easy to install

[Shop Now](#)



Wheelchair/Scooter Ramps

Item #	Ramp Length:
STDS1093	2"
STDS1094	3"
STDS1095	4"
STDS1096	5"
STDS1097	6"

- Designed to transition wheelchair or scooter from one surface to another
- Durable, lightweight construction
- Comes with carry bag
- Perforated slots keep water from collecting on ramp
- Variety of sizes to meet many elevation changes
- Weight capacity: 600 lbs

[Shop Now](#)



Oxygen Cylinder Carry Bag

Item # STDS6008-1

- The Oxygen Cylinder Carry Bag adjusts for use with "D" and "E" cylinders
- Top and bottom clasps easily attach for a secure fit to your folding power wheelchair
- Constructed of durable, easy-to-clean nylon, the carry bag also includes two carry handles

[Shop Now](#)



Universal Cup Holder

Item # STDS1040S

- New clamp allows cup holder to swivel to keep contents level
- Easily attaches to wheeled walkers, rollators, transport chairs/wheelchairs with tubing diameter from $\frac{5}{8}$ " to 1"
- Special clamp insert can be added or removed as necessary to accommodate different size tubing
- Made with durable, lightweight plastic

[Shop Now](#)





CORPORATE OFFICE  99 Seaview Blvd., Port Washington, NY 11050  877-374-8763  www.drivemedical.com

CANADA OFFICE  256 Aviva Park Drive, Vaughan, Ontario Canada L4L 9C7  866-752-1363

©2026 Medical Depot, Inc. dba Drive DeVilbiss Healthcare, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive DeVilbiss are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.