



2026

Over 100 Years
of Excellence

Our veterans deserve
the best. We are
honored to deliver it.



drive 
MEDICAL

DeVilbiss
by drive

©2026 Medical Depot, Inc. dba Drive Medical,
99 Seaview Boulevard, Port Washington, NY 11050.
All trademarks used in association with the sale of
products of Drive Medical are trademarks owned
by Medical Depot, Inc. All other trademarks, trade
names, service marks and logos referenced herein
belong to their respective companies.

Over 100 Years of Excellence

For more than 135 years, DeVilbiss has been a trusted name in respiratory innovation and patient care. Founded in Ohio in 1888, the DeVilbiss legacy began when a country doctor recognized the need for better medical equipment—and set out to create it.

Today, that legacy lives on as DeVilbiss by Drive, combining DeVilbiss's long-standing respiratory expertise with Drive Medical's scale, manufacturing strength, and commitment to quality. We proudly manufacture in Somerset, Pennsylvania, where skilled American workers assemble respiratory equipment our Veterans depend on.

From oxygen concentrators to wheelchairs, we are honored to serve our Veterans with TAA- and BAA-compliant products—backed by a heritage that spans generations.



Stationary Oxygen Concentrators

USA-assembled 5L & 10L concentrators backed by decades of engineering excellence



525 5 Liter Oxygen Concentrator

Item #	Description
525DS	5 Liter Oxygen Concentrator
525DS-Q	5 Liter Oxygen Concentrator, Quiet

- Built upon a legacy of reliable, field-proven concentrator models, the 5 Liter Oxygen Concentrator was designed with enhanced quality, durability and simplicity in mind
- Turn-Down Technology minimizes wear on internal components, reduces power consumption by 15% or more and extends the life of the concentrator



1025 10 Liter Oxygen Concentrator

Item #	Description
1025DS	10 Liter Oxygen Concentrator

- Assembled in the USA and built in the same shell as the 5 Liter, the 10 Liter Concentrator is also compatible with cylinder transfill systems
- With a weight of 42 lbs. and an identical footprint to the 525 unit, the 1025 is substantially smaller and lighter than competitive 10 Liter concentrators



iFill® Personal Oxygen Station

Item #	Description
535D	iFill® Personal Oxygen Station

- Offers a clinically focused and operationally efficient non-delivery oxygen solution a model that is great for patients and your business
- Makes it easy and convenient for your customers to fill portable oxygen cylinders in their homes and gives them greater freedom to come and go as they please

iGo²
SmartDose
TECHNOLOGY



Portable Oxygen Concentrator (POC)

Industry-leading technology that adjusts with every breath



iGo2 Portable Oxygen Concentrator

Item #	Description
125D	iGO2 POC
125D-XB	iGO2 POC, Extra Battery Pack

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Conserver Technology included
- Weighs less than 5 lb.

Oxygen Conservers

Optimize oxygen supply costs and patient outcomes with a range of conservers, including SmartDose® models



SmartDose® Evolution™ Oxygen Conserver

Item #	Description
OM-975	SmartDose® Evolution™ Oxygen Conserver

- Auto-Adjusting Dosing - Monitors breathing patterns for optimal oxygen delivery
- 7 Bolus Settings - Precise dosage adjustments for individual needs
- Extended Battery Life - Uses common AA batteries for up to 2 years of battery life; Auto shut-off



Bonsai® Velocity Pneumatic Oxygen Conserver

Item #	Description
OM-812	Bonsai® Velocity Pneumatic Oxygen Conserver

- Delivers a uniform oxygen pulse at increments: 1, 2, 3, 4, 5 and 6
- Two continuous flow (CF) settings preset @ 2 LPM and 4 LPM
- Delivers up to a 6:1 savings ratio from 14 –40 breaths per minute
- Light weight and quiet operation

Suction Units

Hospital-grade power meets whisper-quiet operation and truly portable systems



Vacu-Aide® Suction Machine



Item #	Description
7325D-D-EXF	Without battery, with external filter
7325D-D	Without battery, with internal filter
7325P-D-EXF	With battery, external filter & carry bag
7325P-D	With battery, internal filter & carry bag

- Features reduced weight and noise level with improved battery function.
- The suction unit utilizes the widely available standard canisters for added convenience and ease of use.

Heavy-duty Suction Machine



Item #	Description
18600	Heavy-Duty Suction Machine

- Combines the performance expected for an institution with the ease of operation needed for homecare use
- Includes 800 cc suction canister, 6' suction tube, 10" suction canister tubing, hydrophobic filter, plastic elbow connector and manual
- Vacuum levels up to 560 mmHg



Interested in more information on our TAA and BAA compliant product line? Contact your sales representative today!