

**Cirrus Plus EC/HD Folding
Power Wheelchair
Order Form/Price List**



Company Name:		Phone:	
Contact Person:		P.O. #	Date:
Bill to:		Ship to:	

Chair with Seating

Sling Seat

- CPN16FBA EC, 16" (W) x 16" (D).....\$6,222.00
- CPN18FBA EC, 18" (W) x 16" (D).....\$6,222.00
- CPN20FBA EC, 20" (W) x 16" (D).....\$6,784.50
- CPN22FBA HD, 22" (W) x 16" (D).....\$8,439.00

Armrests

- Adjustable Height, Removable, Full Length, Flip Up Armrests
(HCPCS E0973).....Included

Controller Arm

- Adjustable Length.....Included
- AA4000 Swingaway Controller Arm (HCPCS E1028)\$247.50

Front Riggings

- Adjustable FootplateIncluded
- CPELR Elevating Legrest Kit\$328.00
Includes left & right swing-away ELRs (HCPCS E0990).

Batteries & Charger

- 12V x 36AH (U1) Batteries (pair) (HCPCS E2365)Included
- 4A Battery ChargerIncluded

Accessories

- Positioning/Safety Belt (HCPCS E0978)Included
- AB1010 Armrest Bag\$25.00

Specifications

HCPCS

Cirrus Plus EC	K0822
Cirrus Plus HD	K0824

Maximum Weight Capacity

Cirrus Plus EC	300 lbs.
Cirrus Plus HD	400 lbs.

Performance

Maximum Speed	5 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°

Turning Radius

CPN16FBA	30"
CPN18FBA	32"
CPN20FBA	33"
CPN20FBA	35"

Dimensions

Overall Length	42.5"
Overall Width (CPN16FBA)	24"
Overall Width (CPN18FBA)	26"
Overall Width (CPN20FBA)	28"
Overall Width (CPN20FBA)	29"
Floor to Seat Height	18.5"
Ground Clearance	4"

Power

Motor	Two 24V x 320W x 4600rpm
Controller	PG VSi 50A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires

Casters	8" x 2" Solid
Drive Wheels	12.5" x 2.25" Flat Free

Weights

Base Weight	90 lbs.
Battery Weight	56 lbs.
Total Weight	146 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.