

\_\_\_ Order Form

February 2014

\_\_\_ Request Quote

Name: \_\_\_\_\_ Phone Number: \_\_\_\_\_

Business Name \_\_\_\_\_ Drop-Ship/Ship To Name: \_\_\_\_\_

Account #: \_\_\_\_\_ PO #: \_\_\_\_\_ Date: \_\_\_\_\_

Bill To: \_\_\_\_\_ Marked For: \_\_\_\_\_

Address: \_\_\_\_\_ City: \_\_\_\_\_ Zip: \_\_\_\_\_

Ship to: \_\_\_\_\_

Address: \_\_\_\_\_ City: \_\_\_\_\_ Zip: \_\_\_\_\_

**MSS TILT & RECLINE**

		<b>HCPCS CODE</b>	<b>MSRP</b>	<b>TOTAL PRICE</b>
<input type="checkbox"/>	MS 3000N MMS Tilt & Recline, Medium	E1399	\$1,275.00	\$ _____

Includes: Lateral thoracic supports, Adductors (2 sizes), Abductor, Headrest with Multi Axis Hardware, Zippered Chest Harness, Pelvic Belt, Angle Adjustable Armrests, Angle Adjustable Footrests, Toe Straps and a Tray.

**MSS TILT & RECLINE ACCESSORIES**

		<b>HCPCS CODE</b>	<b>MSRP</b>	<b>TOTAL PRICE</b>
<input type="checkbox"/>	MS 0014N High Base	E1399	\$140.00	\$ _____
<input type="checkbox"/>	MS 0015N Low Base	E1399	\$140.00	\$ _____
<input type="checkbox"/>	MS 3019 Removable Vinyl Covers	E1399	\$270.00	\$ _____

**SPECIFICATIONS**

Item #	Seat Width	Seat Depth	Back Height	Overall Width	Seat to Footrest Height	Seat to Floor Height	Distance Between Laterals	Headrest Height Adds
MS 3000N	9"-10.5"	3.5"-12"	14"-19"	18"	4"- 10.5"	14"	9"-15"	4"-10"
Item #	Product Weight		Carton Shipping Weight		Weight Limit			
MS 3000N	35 lbs.		40 lbs.		100 lbs.			
Item #	Seat to Floor Height (chair mounted on base)		Product Weight		Carton Shipping Weight		Weight Limit	
MS 001 4N	24"		6 lbs.		9 lbs.		100 lbs.	
MS 001 5N	17"		6 lbs.		9 lbs.		100 lbs.	